

# PROFESSIONAL **PG**x™
















# PERSONAL. PROVEN. PRECISE.

Patient: Addy Hadad Sample ID: 0000096272 Patient DOB: 12/25/1985 10101010 Accession ID: Joe Clinician Ordering Clinician: Sample Collection Date: 9/10/2018 Sample Type: Buccal Sample Received Date: 9/14/2018

Assay Ordered: Genomind Professional PGx 3.0 Report Date: 9/17/2018 10:17 AM

#### **Electronically Signed By**

David Robbins, PhD, DABCC, MT (AAB), Lab Director for Genomind, Inc.

#### **Literature Information Reviewed By**

David Krause, M.D., Chief Medical Officer for Genomind, Inc.

Genomind Professional PGx is intended to assist health care professionals in the selection of safe and appropriate pharmaceuticals and other treatment modalities for patients with mental illness and other brain disorders. This report is designed to be adjunctive to a complete patient assessment, including, but not limited to, proper diagnosis, clinical history, assessment of concomitant comorbidities and medications, family history, and other factors.

# **Personalized Consultation Available for Clinicians**

A complimentary consultation, performed by our expert psychopharmacologists, is included with all Genomind® Professional  $PGx^{TM}$  tests. Consultations can be scheduled directly from the <u>Genomind Clinician Portal</u>.

#### **CONTACT INFORMATION**

Phone: 877-895-8658

E-mail: customerservice@genomind.com

www.genomind.com

#### **TABLE OF CONTENTS**

- I. PHARMACODYNAMIC GENE VARIATIONS
- II. PHARMACOKINETIC GENE VARIATIONS
- III. GENE DRUG INTERACTION SUMMARY
- IV. DIAGNOSIS SPECIFIC SUMMARY\*
- V. TEST METHODOLOGY/LITERATURE REFERENCE
- VI. PATIENT'S GENOMIND RX METATYPE™ CARD

**Disclaimer:** The following report provides a summary of the pharmacokinetic and pharmacodynamic impact certain genes can have on particular drugs. This report is intended to serve as a guide for health care professionals to compare different medication options based on an individual patient's genetics. This report is not intended to recommend a particular course of treatment or medication for a patient. Prescribing health care professionals must use their independent medical judgment and are solely responsible for determining the most appropriate medication for their patients. The clinician must consider other relevant clinical factors in determining which is the most appropriate medication. The test results in this report are intended to be prognostic and not diagnostic. The understanding of the relationship between genetics and pharmacokinetics and pharmacodynamics changes periodically; this report will not be updated to reflect new information. A White Paper summarizing individual gene-drug associations, strength of evidence and effect size is available upon request from Genomind Customer Service.



#### GENOMIND DRUG INTERACTION GUIDE

You have access to G-DIG, the Genomind Drug Interaction Guide. You may utilize this CYP450 Gene-Drug-Environmental checker at any time by selecting the icon next to your patient's name on the Genomind Results Manager screen after signing onto portal genomind.com from your computer or smart phone. The G-DIG technology is included and available at no additional cost for any patient tested with Genomind Professional PGx.

0000096272 Page 2 of 16 © 2019 Genomind, Inc. Confidential

<sup>\*</sup>Diagnosis specific summaries are available for the diagnoses of depression, anxiety & related disorders, bipolar disorder, pain management and ADHD. The provided pages in this report are the closest fit for this individual's diagnosis, as provided to us. All 5 summaries, however, are available to you on the <u>Genomind Clinician Portal</u>.

# I. PHARMACODYNAMIC GENE VARIATIONS

GENE RESULT	THERAPEUTIC IMPLICATIONS	GUIDE	CLINICAL IMPACT
SLC6A4  L(A)/L(A) [High Activity]	Serotonin Transporter (SLC6A4) is a synaptic transporter protein responsible for serotonin reuptake  • Patients with the L(A)/L(A) genotype may have improved likelihood of remission and/or reduced side effects with SSRIs	Q	Therapeutic options: SSRIs if clinically indicated
ADRA2A  C/G [Improved response]	<ul> <li>Alpha-2A Adrenergic Receptor (ADRA2A) is a receptor which plays an important role in norepinephrine signaling</li> <li>Improved response to stimulants (mostly methylphenidate studies) for symptoms of attention deficit/hyperactivity disorder in children and adolescents as compared to those with the C/C genotype</li> </ul>	0	Therapeutic options: methylphenidate may be considered for attention deficit/hyperactivity disorder if clinically indicated
HTR2A  G/G [Normal response]	Serotonin Receptor 2A (HTR2A) is a serotonin receptor which is a target for several serotonergic drugs  • This genotype confers normal activity		No known significant clinical impact
BDNF Val/Val [Normal activity]	Brain-derived Neurotrophic Factor (BDNF) is a protein involved in neuronal development and neural plasticity  • This genotype confers normal activity		No known significant clinical impact
MTHFR  C677T: C/C A1298C: A/C [Normal activity]	<ul> <li>Methylenetetrahydrofolate Reductase (MTHFR) is an enzyme responsible for the conversion of folic acid to methylfolate, which is a cofactor needed for serotonin, norepinephrine, and dopamine synthesis</li> <li>This genotype confers normal activity</li> </ul>		No known significant clinical impact
COMT  Val/Met [Normal activity]	Catechol-O-Methyltransferase (COMT) is an enzyme responsible for breakdown of dopamine in the frontal cortex of the brain  COMT is involved in response to stimulants  This genotype confers normal activity		No known significant clinical impact
HLA-A *31:01  Not Detected [Normal]	<ul> <li>Major histocompatibility complex, class I, A (HLA-A) is part of a cluster of genes known as the Human Leukocyte Antigen complex</li> <li>Certain variants greatly increase risk of drug induced skin reactions</li> <li>This genotype is associated with normal risk of skin reactions with carbamazepine</li> </ul>		Normal risk of skin reactions with carbamazepine
Not Detected [Normal]	<ul> <li>Major histocompatibility complex, class I, B (HLA-B) is part of a cluster of genes known as the Human Leukocyte Antigen complex</li> <li>Certain variants greatly increase risk of drug induced skin reactions</li> <li>This genotype is associated with normal risk of skin reactions with carbamazepine, oxcarbazepine, phenytoin and fosphenytoin</li> </ul>		Normal risk of skin reactions with carbamazepine, oxcarbazepine, phenytoin/fosphenytoin
DRD2  C/C [Normal activity]	<ul> <li>Dopamine Receptor D2 (DRD2) is a receptor activated by dopamine in the brain</li> <li>DRD2 is involved in response to antipsychotics</li> <li>This genotype confers normal activity</li> </ul>		No known significant clinical impact



Alert/Caution



PGx Guided Options

0000096272 Page 3 of 16 © 2019 Genomind, Inc. Confidential

# I. PHARMACODYNAMIC GENE VARIATIONS

GENE RESULT	THERAPEUTIC IMPLICATIONS	GUIDE	CLINICAL IMPACT
MC4R  C/C [Normal activity]	<ul> <li>Melanocortin 4 Receptor (MC4R) is a receptor that plays a central role in the control of food intake</li> <li>MC4R is involved in antipsychotic-induced weight gain</li> <li>This genotype confers normal activity and average risk of weight gain</li> </ul>		No known significant clinical impact
5HT2C  C/C  [Standard weight gain risk]	<ul> <li>Serotonin Receptor 2C (5HT2C) is a receptor involved in the regulation of satiety</li> <li>Some 2nd generation antipsychotics act by blocking this receptor</li> <li>Patients with the C/C genotype have standard risk of weight gain with 2nd generation antipsychotics. C/C is the most common genotype</li> <li>Higher risk: clozapine; olanzapine</li> <li>Medium risk: aripiprazole; brexpiprazole; iloperidone; paliperidone; quetiapine; risperidone</li> <li>Lower risk: asenapine; cariprazine; lurasidone; ziprasidone</li> </ul>		Assess weight gain risk with various second generation antipsychotics
ANK3  C/C  [Normal activity]	<ul> <li>Sodium Channel (ANK3) is a protein that plays a role in sodium ion channel function and is involved in excitatory signaling in the brain</li> <li>This genotype confers normal activity</li> </ul>		No known significant clinical impact
CACNA1C  G/G [Normal activity]	Calcium Channel (CACNA1C) is a subunit of L-type voltage gated calcium channels which are involved in excitatory signaling in the brain  • This genotype confers normal activity		No known significant clinical impact
OPRM1  A/A  [Normal activity]	<ul> <li>μ-Opioid Receptor (OPRM1) is an opioid receptor which is affected by endogenous and exogenous opioids</li> <li>OPRM1 is involved in response to opioids</li> <li>This genotype confers normal activity</li> </ul>		No known significant clinical impact
GRIK1  A/A [Normal activity]	<ul> <li>Glutamate Receptor Kainate 1 (GRIK1) is an excitatory neurotransmitter receptor</li> <li>GRIK1 is involved in response to topiramate for alcohol abuse</li> <li>Patients of European descent with the A allele may be less likely to respond to topiramate for alcohol use disorder; future studies, however, are needed to confirm these findings</li> </ul>		No known significant clinical impact



Alert/Caution



PGx Guided Options

0000096272 Page 4 of 16 © 2019 Genomind, Inc. Confidential

# **II. PHARMACOKINETIC GENE VARIATIONS**

GENE RESULT	THERAPEUTIC IMPLICATIONS	GUIDE	CLINICAL IMPACT
CYP2D6 IM  *4/*10  Duplication [Intermediate activity]	<ul> <li>Intermediate metabolizer: Risk of elevated serum levels &amp; drug interactions, or decreased production of active metabolites</li> <li>A dose adjustment or alternate therapy may be considered</li> </ul>	1	Be advised that there may be altered exposure to medications metabolized by CYP2D6  Use G-DIG for a more complete drug-gene-environment
			interaction assessment
CYP1A2 EM *1Dc/*1Vc	Variations in the CYP1A2 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)
[Normal activity]	• Each of the CYP1A2 variants detected in this patient sample is well characterized, although this specific combination of alleles has not been formally named. We have adopted a modified (*)star allele naming system that identifies all the variants detected for this gene. (Adapted from Soyama et al 2005. PMID: 15770072; Gunes et al 2009. PMID: 19450128)		Use G-DIG for a more complete drug-gene-environment interaction assessment
CYP2B6 EM *1/*1	Variations in the CYP2B6 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)
[Normal activity]			Use G-DIG for a more complete drug-gene-environment interaction assessment
CYP2C9 EM	Variations in the CYP2C9 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)
*1/*1 [Normal activity]			Use G-DIG for a more complete drug-gene-environment interaction assessment
CYP2C19 EM	Variations in the CYP2C19 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)
*1/*1 [Normal activity]			Use G-DIG for a more complete drug-gene-environment interaction assessment
CYP3A4 *1/*1 CYP3A5 *7/*2	Variations in the CYP3A4/5 liver enzymes can result in altered drug metabolism and unexpected drug serum levels  • 3A5 non-expresser		Normal metabolism is expected (other factors may influence metabolism)
*7/*3 [Normal activity]	<ul> <li>CYP3A activity is determined by the sum activity of the CYP3A family of genes; in adults the most influential are 3A4 and 3A5</li> <li>This genotype confers normal activity</li> </ul>		Use G-DIG for a more complete drug-gene-environment interaction assessment



Alert/Caution



PGx Guided Options

0000096272 Page 5 of 16 © 2019 Genomind, Inc. Confidential

# **II. PHARMACOKINETIC GENE VARIATIONS**

GENE RESULT	THERAPEUTIC IMPLICATIONS	GUIDE	CLINICAL IMPACT
UGT1A4 EM *1a/*1a [Normal activity]	Variations in the UGT1A4 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)  Use G-DIG for a more complete drug-gene-environment interaction assessment
UGT2B15 EM *1/*1 [Normal activity]	Variations in the UGT2B15 liver enzyme can result in altered drug metabolism and unexpected drug serum levels  • This genotype confers normal activity		Normal metabolism is expected (other factors may influence metabolism)  Use G-DIG for a more complete drug-gene-environment interaction assessment
ABCB1 (rs2032583) A/A [Normal activity]	ATP Binding Cassette B1 (ABCB1) encodes for P-glycoprotein (P-gp). P-gp is a drug efflux pump that reduces the intestinal absorption and blood-brain barrier penetration of certain drugs  • This genotype is associated with normal activity of P-gp		Normal exposure is expected (other factors may influence drug exposure)  Use G-DIG for a more complete drug-gene-environment interaction assessment
ABCB1 (rs1045642) G/G [Normal activity]	ATP Binding Cassette B1 (ABCB1) encodes for P-glycoprotein (P-gp). P-gp is a drug efflux pump that reduces the intestinal absorption and blood-brain barrier penetration of certain drugs  • This genotype is associated with normal activity of P-gp		Normal exposure is expected (other factors may influence drug exposure)  Use G-DIG for a more complete drug-gene-environment interaction assessment



Alert/Caution



PGx Guided Options

0000096272 Page 6 of 16 © 2019 Genomind, Inc. Confidential

# **III. GENE DRUG INTERACTION SUMMARY**

CLASS	MEDI	CATION	PHARMACODYNAMIC ASSOCIATIO	ONS	PHARMACODYNAMIC GENE	DRUG EXPOSURE	PHARMACOKINETIC GENE
	ANTI	DEPRESSANTS					
	₽≡	Citalopram (Celexa®)	Higher odds of remission or re	esponse	SLC6A4		2C19, P-gp
		Escitalopram (Lexapro®)	Higher odds of remission or re	esponse	SLC6A4		2C19, P-gp
SSRIs		Fluoxetine (Prozac®)	Higher odds of remission or re	esponse	SLC6A4	$\uparrow$	2D6, 2C9
SS		Fluvoxamine (Luvox®)	Higher odds of remission or re	esponse	SLC6A4	$\uparrow$	2D6, 1A2, P-gp
		Paroxetine (Paxil®)	Higher odds of remission or re	esponse	SLC6A4	$\uparrow$	2D6, P-gp
		Sertraline (Zoloft®)	Higher odds of remission or re	esponse	SLC6A4		2C19, 2B6
		Desvenlafaxine (Pristiq®)					
SNRIs		<b>Duloxetine</b> (Cymbalta®)				$\uparrow$	1A2, 2D6
SN		Levomilnacipran (Fetzima®)					3A4/5
		Venlafaxine[1] (Effexor®)				$\uparrow$	2D6, 2C19, 3A4/5, P- gp
		Bupropion[1] (Wellbutrin®)					2B6
		Esketamine (Spravato®)					2B6, 3A4/5
		Mirtazapine (Remeron®)				$\uparrow$	2D6, 3A4/5, 1A2
Other		Nefazodone					3A4/5
		Trazodone (Desyrel®, Oleptro®)				$\uparrow$	3A4/5, 2D6
		Vilazodone (Viibryd®)					3A4/5
	R	Vortioxetine (Trintellix®)				$\uparrow$	2D6, 3A4/5
	1	Alert/Caution	PGx Guided Options		ced Drug Exposure 1A2 Inducers	$\uparrow\downarrow$	Drug Exposure

 $\ensuremath{\mathbb{R}}\xspace$  [1] See Gene Drug Interaction Summary footnotes for more information

0000096272 Page 7 of 16 © 2019 Genomind, Inc. Confidential

# **III. GENE DRUG INTERACTION SUMMARY**

LASS	MEDIC	CATION	PHARMACODYNAMIC ASSOCIATI	ONS	PHARMACODYNAMIC GENE	DRUG EXPOSURE	PHARMACOKINETIC GENE
	ANTI	DEPRESSANTS					
		Amitriptyline (Elavil®)				$\uparrow$	2D6, 2C19, P-gp
		Amoxapine (Asendin®)				$\uparrow$	2D6
		Clomipramine (Anafranil®)				$\uparrow$	2D6, 2C19, 1A2
		Desipramine (Norpramin®)				$\uparrow$	2D6
TCAs		<b>Doxepin</b> (Sinequan®)				$\uparrow$	2D6, 2C19
		Imipramine (Tofranil®)				$\uparrow$	2D6, 2C19
		Nortriptyline (Pamelor®)				$\uparrow$	2D6, P-gp
		Protriptyline (Vivactil®)				$\uparrow$	2D6
		Trimipramine (Surmontil®)				$\uparrow$	2D6, 2C19, P-gp
		Phenelzine (Nardil®)					
MAOIs		Selegiline (Eldepryl®, Emsam®)					2B6
		Tranylcypromine (Parnate®)					
	МОО	D STABILIZERS/ANTIC	ONVULSANTS				
	₽≣	Carbamazepine (Equetro®, Tegretol®)					3A4/5
		Gabapentin (Neurontin®)					
		Lamotrigine (Lamictal®)					UGT1A4
		<b>Lithium</b> (Lithobid®, Eskalith®)					
	₽₿	Oxcarbazepine (Trileptal®, Oxtellar®)					
		Pregabalin (Lyrica®)					
		Topiramate (Topamax®)					
		Valproate (Depakote®, Depakene®)					<b>2</b> C9
	A	Alert/Caution	PGx Guided Options		ed Drug Exposure .A2 Inducers	$\uparrow\downarrow$	Drug Exposure

0000096272 Page 8 of 16 © 2019 Genomind, Inc. Confidential

# **III. GENE DRUG INTERACTION SUMMARY**

CLASS	MEDI	CATION	PH/	ARMACODYNAMIC ASSOCIA	TIONS	PHARMACODYNAM GENE	IC DRUG EXPOSURE	PHARMACOKINETIC GENE
	ANTI	PSYCHOTICS						
	P	<b>Aripiprazole</b> (Abilify®)					$\uparrow$	2D6, 3A4/5, P-gp
		Asenapine (Saphris®)						1A2, UGT1A4
	B	Brexpiprazole (Rexulti®)					$\uparrow$	2D6, 3A4/5
		Cariprazine (Vraylar®)						3A4/5
2nd Generation Antipsychotics		Clozapine (Clozaril®)					$\uparrow$	1A2, 2D6, 3A4/5, P-g
ıtipsyc	ß	Iloperidone (Fanapt®)					$\uparrow$	2D6, 3A4/5
ion Ar		<b>Lurasidone</b> (Latuda®)						3A4/5
enerat		Olanzapine (Zyprexa®)						1A2, P-gp
2nd G		Paliperidone (Invega®)						
		Pimavanserin (Nuplazid®)						3A4/5
		Quetiapine (Seroquel®)						3A4/5
		Risperidone (Risperdal®)					$\uparrow$	2D6, 3A4/5, P-gp
		<b>Ziprasidone</b> (Geodon®)						
	1	Alert/Caution	Q	PGx Guided Options		Reduced Drug Exposure with 1A2 Inducers	$\uparrow\downarrow$	Drug Exposure

 $\ensuremath{\mathbb{R}}$  [1] See Gene Drug Interaction Summary footnotes for more information

0000096272 Page 9 of 16 © 2019 Genomind, Inc. Confidential

# **III. GENE DRUG INTERACTION SUMMARY**

ANTI	PSYCHOTICS					
	Chlorpromazine (Thorazine®)				$\uparrow$	2D6
	Fluphenazine (Prolixin®)				$\uparrow$	2D6
	Haloperidol (Haldol®)				$\uparrow$	2D6, 3A4/5
	Loxapine (Adasuve®, Loxitane®)					3A4/5, 1A2
	Perphenazine (Trilafon®)				$\uparrow$	2D6
R	Pimozide (Orap®)				$\uparrow$	2D6, 3A4/5
R	Thioridazine (Mellaril®)				$\uparrow$	2D6
	Thiothixene (Navane®)					1A2
	Trifluoperazine (Stelazine®)					1A2, UGT1A4
ANXI	OLYTICS					
	Alprazolam (Xanax®)					3A4/5
	Buspirone (Buspar®)					3A4/5
	Chlordiazepoxide (Librium®)					3A4/5, UGT2B1
	Clonazepam (Klonopin®)					3A4/5
	Clorazepate (Tranxene®)					UGT2B15
	<b>Diazepam</b> (Valium®)					2C19, 3A4/5, UGT2B15
	Hydroxyzine (Vistaril®)					
	Lorazepam (Ativan®)					UGT2B15
	Oxazepam (Serax®)					UGT2B15
	Temazepam (Restoril®)					UGT2B15
(	₽ ₽	(Haldol®)  Loxapine (Adasuve®, Loxitane®)  Perphenazine (Trilafon®)  Pimozide (Orap®)  Thioridazine (Mellaril®)  Thiothixene (Navane®)  Trifluoperazine (Stelazine®)  NXIOLYTICS  Alprazolam (Xanax®)  Buspirone (Buspar®)  Chlordiazepoxide (Librium®)  Clonazepam (Klonopin®)  Clorazepate (Tranxene®)  Diazepam (Valium®)  Hydroxyzine (Vistaril®)  Lorazepam (Ativan®)  Oxazepam (Serax®)  Temazepam	(Haldol®)  Loxapine (Adasuve®, Loxitane®)  Perphenazine (Trilafon®)  Pimozide (Orap®)  Thioridazine (Mellaril®)  Thiothixene (Navane®)  Trifluoperazine (Stelazine®)  INXIOLYTICS  Alprazolam (Xanax®)  Buspirone (Buspar®)  Chlordiazepoxide (Librium®)  Clonazepam (Klonopin®)  Clorazepate (Tranxene®)  Diazepam (Valium®)  Hydroxyzine (Vistaril®)  Lorazepam (Ativan®)  Oxazepam (Serax®)  Temazepam (Restoril®)	(Haldol®)  Loxapine (Adasuve®, Loxitane®)  Perphenazine (Trilafon®)  Thioridazine (Mellaril®)  Thioridazine (Navane®)  Trifluoperazine (Stelazine®)  INXIOLYTICS  Alprazolam (Xanax®)  Buspirone (Buspar®)  Chlordiazepoxide (Librium®)  Clonazepam (Klonopin®)  Clorazepate (Tranxene®)  Diazepam (Valium®)  Hydroxyzine (Vistaril®)  Lorazepam (Ativan®)  Oxazepam (Serax®)  Temazepam (Restoril®)	(Haldol*)  Loxapine (Adasuve*, Loxitane*)  Perphenazine (Trilafon*)  Pimozide (Orap*)  Thioridazine (Melaril*)  Thiotikene (Navane*)  Trifluoperazine (Stelazine*)  NXIOLYTICS  Alprazolam (Xanax*)  Buspirone (Buspar*)  Chlordiazepoxide (Librium*)  Clonazepam (Klonopin*)  Clorazepate (Tranxene*)  Diazepam (Valium*)  Hydroxyzine (Vistaril*)  Lorazepam (Ativan*)  Oxazepam (Serax*)  Temazepam (Serax*)  Temazepam (Restoril*)	(Haldol*)  Loxapine (Adasuve*, Loxitane*)  Perphenazine (Trilafon*)  Pimozide (Orap*)  Thiothixene (Navane*)  Trifluoperazine (Stelazine*)  NXIOLYTICS  Alprazolam (Xanax*)  Buspirone (Buspar*)  Chlordiazepoxide (Librium*)  Clonazepam (Konopin*)  Clorazepate (Tranxene*)  Diazepam (Valtium*)  Hydroxyzine (Vistari*)  Lorazepam (Ativan*)  Oxazepam (Ativan*)  Oxazepam (Serax*)  Temazepam (Restorii*)

[1] See Gene Drug Interaction Summary footnotes for more information

0000096272 Page 10 of 16 © 2019 Genomind, Inc. Confidential

# **III. GENE DRUG INTERACTION SUMMARY**

LASS	MEDIC	ATION	PHARMACODYNAMIC ASSOCIATION	ONS	PHARMACODYNAMIC GENE	DRUG EXPOSURE	PHARMACOKINETI GENE
	ADHE	MEDICATIONS					
		Amphetamine- Dextroamphetamine (Adderall®, Evekeo®)				$\uparrow$	2D6
3		<b>Dexmethylphenidate</b> (Focalin®)	Higher odds of response		ADRA2A		
		Dextroamphetamine (Dexedrine®, Procentra®, Zenzedi®)				$\uparrow$	2D6
5		<b>Lisdexamfetamine</b> (Vyvanse®)				$\uparrow$	2D6
2		Methamphetamine (Desoxyn®)				$\uparrow$	2D6
		Methylphenidate (Ritalin®, Concerta®, Daytrana®, Metadate®)	Higher odds of response		ADRA2A		
	₽Ħ	Atomoxetine (Strattera®)				$\uparrow$	2D6
Other		Clonidine (Kapvay®)				$\uparrow$	2D6
		Guanfacine (Intuniv®)					3A4/5
	SUPP	LEMENTS					
		L-methylfolate (Deplin®)					
	SLEEP	MODULATORS					
		Armodafinil (Nuvigil®)					3A4/5
		Eszopicione (Lunesta®)					3A4/5
		Modafinil (Provigil®)					3A4/5
		Ramelteon (Rozerem®)					1A2
		Suvorexant (Belsomra®)					3A4/5, 2C19
		Zaleplon (Sonata®)					3A4/5
		<b>Zolpidem</b> (Ambien®)					3A4/5
	1	Alert/Caution	PGx Guided Options		ed Drug Exposure .A2 Inducers	$\uparrow\downarrow$	Drug Exposure

 $\ensuremath{\mathbb{R}}$  [1] See Gene Drug Interaction Summary footnotes for more information

# **III. GENE DRUG INTERACTION SUMMARY**

LASS	MEDIC	CATION	PHARMACODYNAMIC ASSOCIATIONS	PHARMACODYNAMIC GENE	DRUG EXPOSURE	PHARMACOKINETIC GENE
	PAIN					
		Acetaminophen (Tylenol®)				
	R	Celecoxib (Celebrex®)				2C9
?		<b>Diclofenac</b> (Voltaren®, Cataflam®)				2C9
		Flurbiprofen (Ansaid®)				2C9
5		Ibuprofen (Advil®, Motrin®)				2C9
		Ketorolac (Toradol®)				
•		Meloxicam (Mobic®)				2C9
		Naproxen (Aleve®, Naprosyn®)				1A2, 2C9
		Piroxicam (Feldene®)				2C9
		Alfentanil (Alfenta®)				3A4/5
	₽ij	Codeine[1]			$\uparrow$	2D6, P-gp
		Fentanyl (Duragesic®)				3A4/5, P-gp
		Hydrocodone[1] (Vicodin®, Norco®, Lorcet®)			$\uparrow$	2D6, 3A4/5
		Hydromorphone (Dilaudid <sup>®</sup> )				
)		Meperidine (Demerol®)				2B6, 3A4/5
		Methadone (Dolophine®, Methadose®)				3A4/5, 2B6
7		Morphine (MS Contin®, Kadian®)				P-gp
		Oxycodone (Oxycontin®)			$\uparrow$	2D6, 3A4/5, P-g <sub>l</sub>
		Oxymorphone (Opana®)				
		<b>Tapentadol</b> (Nucynta®)				
	⊞	Tramadol[1] (Ultram®)			$\uparrow$	2D6, 3A4/5, P-g <sub>l</sub>

#### III. GENE DRUG INTERACTION SUMMARY

LASS	MEDIC	ATION	PHARMACODYNAMIC ASSOCIATIO	NS	PHARMACODYNAMIC GENE	DRUG EXPOSURE	PHARMACOKINETIC GENE
N	MISC	ELLANEOUS					
		<b>Dextromethorphan/Quinidine</b> (Nuedexta®)				$\uparrow$	2D6, 3A4/5
		Baclofen (Lioresal®)					
		<b>Buprenorphine/Naloxone</b> (Suboxone®)					3A4/5
		Buprenorphine (Butrans®)					3A4/5
		Cannibidiol (CBD) (Epidiolex®)					3A4/5, 2C19
		Carisoprodol (Soma®)					2C19
		Cyclobenzaprine (Flexeril®)					1A2
	B	<b>Deutetrabenazine</b> (Austedo®)				$\uparrow$	2D6
		Metaxalone (Skelaxin®)					
		Methocarbamol (Robaxin®)					
		Naltrexone (Revia®, Vivitrol®)					
•	R	Phenytoin/Fosphenytoin (Dilantin®, Cerebyx®)					2C19, 2C9
		<b>Tizanidine</b> (Zanaflex®)					1A2
	R	Valbenazine (Ingrezza®)				$\uparrow$	3A4/5, 2D6
	1	Alert/Caution	PGx Guided Options		ced Drug Exposure LA2 Inducers	$\uparrow\downarrow$	Drug Exposure

[1] See Gene Drug Interaction Summary footnotes for more information


# **GENE DRUG INTERACTION SUMMARY FOOTNOTES**

[1] Prodrug or highly active metabolite - requiring activation by the liver; CYP450 IMs/PMs may experience lower efficacy due to reduced conversion to the active metabolite and higher levels of the parent drug; CYP450 UMs may experience increased conversion of the parent drug, and higher levels of the active metabolite

- Medication has FDA biomarker guidance available
  - https://www.fda.gov/downloads/Drugs/ScienceResearch/UCM578588.pdf
- - <a href="https://cpicpgx.org/guidelines/">https://cpicpgx.org/guidelines/</a>
  - https://www.pharmgkb.org/page/dpwg
- \*References for the drug interaction summary are available upon request

0000096272 Page 13 of 16 © 2019 Genomind, Inc. Confidential

# **IV. ADHD SUMMARY**



#### V. TEST METHODOLOGY/LITERATURE REFERENCE

#### **TEST METHODOLOGY**

This test was developed and performance characteristics were validated in the Genomind clinical laboratory. This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA). FDA does not require review of this test because it is a Laboratory Developed Test (LDT). This test is used for clinical purposes and should not be regarded as investigational or for research use. This laboratory is certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical laboratory testing. Genomind performed the testing using standard and custom TaqMan reagents for all variants. The test results are intended to be used as prognostic and not diagnostic and are not intended as the sole means for patient management decisions.

Test Methodology Limitations: Factors influencing the amount and quality of DNA extracted include but are not limited to the amount of buccal cells extracted, patient oral hygiene, collection technique, and the presence of dietary or microbial sources of nucleic acids and nucleases. DNA quality and quantity are subject to matrix dependent influences. PCR inhibitors, extraneous DNA and nucleic acid degrading enzymes are all factors which may affect the evaluation of assay results. Some single nucleotide polymorphism (SNP) assays are problematic due to multiple base repeats and other sequence aberrations which may hinder proper amplification and analysis. DNA purity can influence the assay. SLC6A4 contains many polymorphisms and the assay was developed and validated according to the current available scientific information. For pharmacogenetics tests like Genomind Professional PGx, undetected genetic and/or non-genetic factors such as drug-drug interactions may impact the phenotype. The Genomind Professional PGx report is based on a current understanding of the clinical relevance of the variant identified, penetrance, phenotype predictions, and recurrence risks.

Variants tested include 5HT2C rs3813929; ABCB1 C3435T rs1045642; ABCB1 rs2032583; ADRA2A rs1800544; ANK3 rs10994336; BDNF rs6265; CACNA1C rs1006737; COMT rs4680; CYP1A2 \*1B, \*1C, \*1D, \*1E, \*1F, \*1K and \*11; CYP2B6 \*4, \*5, and \*6; CYP2C19 \*2, \*3, \*4, \*5, \*6, \*7, \*8, \*9, \*10, \*17, and \*35; CYP2C9 \*2, \*3, \*4, \*5, \*6, \*8, \*11, \*13, and \*27; CYP2D6 \*2, \*3, \*4, gene deletion (\*5), gene duplication, \*6, \*7, \*8, \*9, \*10, \*11, \*12, \*14, \*15, \*17, \*29 and \*41; CYP3A4 \*22; CYP3A5 \*3, \*6, \*7; DRD2 rs1799732; GRIK1 rs2832407; HLA-B\*15:02 presence with reflex testing for presence of HLA-B\*15:13 for all positive samples and Sanger sequencing for all double positive samples; HLA-A\*31:01 rs1061235; HTR2A rs7997012; MC4R rs489693; MTHFR rs1801131 and rs1801133; OPRM1 rs1799971; SLC6A4 rs25531 and rs63749047; UGT2B15 rs1902023; and UGT1A4 rs2011425. Other known variants that are not listed are not detected and will not be included in the test report.

Version 3.0 [05/20/2019]

#### LITERATURE REFERENCES

THE LITERATURE INFORMATION UPON WHICH THIS REPORT RELIES WAS AGGREGATED AND REVIEWED BY GENOMIND, INC. SUMMARIES OF THESE NUMBERED REFERENCES BELOW ARE AVAILABLE UPON REQUEST OF GENOMIND'S COMPREHENSIVE LITERATURE SUMMARY [V2019-05].

Gene	References
5HT2C	10-24
ADRA2A	25-33
ANK3	34-57
BDNF	58-78
CACNA1C	34, 37-39, 42, 47-55, 79-97
COMT	28, 98-131
DRD2	132-138
GRIK1	139-143
HLA-A *31:01	145-149
HLA-B *15:02	146-147, 150-157
HTR2A	71, 158-164
MC4R	21, 23, 165-172

Gene	References		
MTHFR	173-185		
OPRM1	186-197		
SLC6A4	163, 198-214		
ABCB1	215-233		
UGT1A4	235-238		
UGT2B15	238-241		
CYP1A2	20, 164, 247-250, 252, 258-281, 286-288		
CYP2B6	247-252, 273, 284, 288-305		
CYP2C9	247-257, 261, 288, 306-312		
CYP2C19	15, 242, 244, 247-252, 254-255, 273, 284, 288, 306, 311-323		
CYP2D6	15, 20, 244-245, 247-252, 254-255, 258, 261, 273-274, 288, 306, 311-314, 322, 324-344		
CYP3A4/5	15, 20, 247-252, 258, 261, 273-274, 282-285		

0000096272 Page 15 of 16 © 2019 Genomind, Inc. Confidential

#### VI. PATIENT'S GENOMIND RX METATYPE™ CARD

Your Genomind Rx MetaType™ wallet card includes information on six liver enzymes that are responsible for the metabolism of most drugs, and identifies your unique enzyme profile (your genotype). It is intended for use by your current, additional or future healthcare providers. This genetic information is mentioned in the FDA prescribing information of many drugs, and may provide useful prescribing recommendations. The websites on the back of the card provide more information.

PRO	DFESSIONAL PGx	Bob Sample #0000096272	
	MetaType™ C		
Gene	Genotype	Phenotype	Clinical Meaning*
CYP1A2	*1Dc/*1Vc	Extensive	Normal Metabolism
CYP2B6	*1/*1	Extensive	Normal Metabolism
CYP2C19	*1/*1	Extensive	Normal Metabolism
CYP2C9	*1/*1	Extensive	Normal Metabolism
CYP2D6	*4/*10 DUP	Intermediate	↓ Metabolism of some drugs
CYP3A4	*1/*1, *7/*3	Normal	Normal Metabolism



Issued Date: 09/17/2018

Issued Date: 09/17/2018

#### FOR USE BY HEALTHCARE PROFESSIONALS ONLY

Most medicines are metabolized by liver enzymes. Like blood types, you have a specific genetic profile which can affect the rate of metabolism, and may influence the dose of medicines prescribed for you. You may wish to inform your healthcare provider(s) about your metabolism status, shown on the reverse. More information about specific gene/drug interactions can be found at:

https://drug-interactions.medicine.iu.edu/Clinical-Table.aspx https://www.pharmgkb.org/guidelines https://www.fda.gov/downloads/Drugs/ScienceResearch/UCM578588.pdf

\*Do not discontinue or change the dose of any medicine without the advice of your healthcare provider. In addition to genetics, other factors may influence your metabolizer status.



	NOMIND® DFESSIONAL PGx*	Bob Sample #0000096272	
Rx	MetaType™ C		
Gene	Genotype	Phenotype	Clinical Meaning*
CYP1A2	*1Dc/*1Vc	Extensive	Normal Metabolism
CYP2B6	*1/*1	Extensive	Normal Metabolism
CYP2C19	*1/*1	Extensive	Normal Metabolism
CYP2C9	*1/*1	Extensive	Normal Metabolism
CYP2D6	*4/*10 DUP	Intermediate	$\downarrow$ Metabolism of some drugs
CYP3A4	*1/*1, *7/*3	Normal	Normal Metabolism



FOR USE BY HEALTHCARE PROFESSIONALS ONLY

Most medicines are metabolized by liver enzymes. Like blood types, you have a specific genetic profile which can affect the rate of metabolism, and may influence the dose of medicines prescribed for you. You may wish to inform your healthcare provider(s) about your metabolism status, shown on the reverse. More information about specific gene/drug interactions can be found at:

https://drug-interactions.medicine.iu.edu/Clinical-Table.aspx https://www.pharmgkb.org/guidelines https://www.fda.gov/downloads/Drugs/ScienceResearch/UCM578588.pdf

\*Do not discontinue or change the dose of any medicine without the advice of your healthcare provider. In addition to genetics, other factors may influence your metabolizer status.