



**DX-800 PROCESS ANALYZER
USER'S MANUAL**

©2001 Dionex Corporation

Document No. 031325
Revision 03
November 2001

©2001 by Dionex Corporation
All rights reserved worldwide.
Printed in the United States of America.

This publication is protected by federal copyright law. No part of this publication may be copied or distributed, transmitted, transcribed, stored in a retrieval system, or transmitted into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, manual, or otherwise, or disclosed to third parties without the express written permission of Dionex Corporation, 1228 Titan Way, Sunnyvale, California 94088-3603 U.S.A.

DISCLAIMER OF WARRANTY AND LIMITED WARRANTY

THIS PUBLICATION IS PROVIDED “AS IS” WITHOUT WARRANTY OF ANY KIND. DIONEX CORPORATION DOES NOT WARRANT, GUARANTEE, OR MAKE ANY EXPRESS OR IMPLIED REPRESENTATIONS REGARDING THE USE, OR THE RESULTS OF THE USE, OF THIS PUBLICATION IN TERMS OF CORRECTNESS, ACCURACY, RELIABILITY, CURRENTNESS, OR OTHERWISE. FURTHER, DIONEX CORPORATION RESERVES THE RIGHT TO REVISE THIS PUBLICATION AND TO MAKE CHANGES FROM TIME TO TIME IN THE CONTENT HEREINOF WITHOUT OBLIGATION OF DIONEX CORPORATION TO NOTIFY ANY PERSON OR ORGANIZATION OF SUCH REVISION OR CHANGES.

TRADEMARKS

AutoSuppression™ and DX-LAN are trademarks, and EluGen®, PeakNet®, and SRS® are registered trademarks of Dionex Corp.

Microsoft® and Windows NT® are registered trademarks of Microsoft Corp.

NOWPAK® is a registered trademark of NOW Technologies, Inc.

Teflon® and Tefzel® are registered trademarks of E.I. du Pont de Nemours & Co.

Wonderware® is a registered trademark and InTouch™ is a trademark of Wonderware Corp.

PRINTING HISTORY

Revision 01, May 1998

Revision 02, July 1999

Revision 03, November 2001

Contents

1 • Introduction

1.1	About This Manual	1-2
1.1.1	Safety Messages and Notes	1-3
1.1.2	Safety Labels	1-5
1.1.3	Related Manuals	1-5
1.2	DX-800 Safety Practices	1-6
1.2.1	General Precautions	1-6
1.2.2	Compressed Gas or Liquid Cylinder Precautions	1-6
1.2.3	Mechanical Precautions	1-7
1.2.4	Electrical Precautions	1-7
1.3	Analyzer Hardware Configuration Overview . . .	1-8
1.3.1	Channel Modules	1-9
1.3.2	SS80 Sample Selector (Optional)	1-14
1.3.3	Accessories	1-15
1.4	PeakNet-PA Software	1-16

2 • AE80 Enclosure

2.1	Front Door	2-2
2.2	Electrical System	2-4
2.3	Fluid and Pneumatic Systems	2-6
2.3.1	External Features	2-6

2.3.2	Internal Features	2-6
2.4	Environmental Controls	2-7

3 • CC80 Component Controller

3.1	Overview	3-1
3.2	Operating Features	3-2
3.2.1	Front Panel	3-2
3.2.2	Rear and Side Panels	3-5
3.2.3	DIP Switches	3-6
3.3	Operating Modes	3-9
3.3.1	Local Mode	3-9
3.3.2	Remote Mode	3-9
3.4	Power-Up Diagnostics	3-10

4 • SP80 Sample Preparation Panel

4.1	Overview	4-1
4.2	Configuration: SP81 for Concentration or Direct Injection	4-2
4.3	Configuration: SP82 for Dilution	4-4
4.4	Configuration: SP83 for Dilution with Reagent Addition	4-6
4.5	Configuration: SP84 for Concentration with Reagent Addition	4-8
4.6	Precision Displacement Pumps	4-10
4.7	Dilution Vessel	4-11

4.8 Post-Column Reagent Pump (Optional) 4-11

5 • LC80 Liquid Chromatography Panel

5.1 Overview 5-1

5.2 Chromatography Components 5-2

5.2.1 Load/Inject (LI) Valve 5-3

5.2.2 Detector Cell 5-3

5.2.3 Consumables 5-3

5.2.4 CH-4 Column Heater (Optional) 5-4

5.2.5 Column Switching Valve (Optional) 5-4

6 • EG40-PA Eluent Generator

6.1 Overview 6-1

6.2 Main Components 6-2

6.3 Eluent Generator Control 6-4

7 • LM80 Liquids Manager

7.1 LM80 Control Panel 7-2

7.2 Pneumatic Requirements 7-3

8 • SS80 Sample Selector

8.1 Overview 8-1

8.2 Main Components 8-2

9 • Operation and Maintenance

9.1	Installation Checklist	9-1
9.1.1	Prepare the Site and Facilities	9-1
9.1.2	Connect the Facilities	9-1
9.1.3	Connect the Communications Cables	9-2
9.1.4	Connect the Sample Inlet to the Analyzer Channels	9-2
9.2	Initial Startup	9-3
9.2.1	Turn on the Power	9-3
9.2.2	Start PeakNet-PA	9-3
9.2.3	Flush the Flow Path	9-3
9.2.4	Install Consumable Components	9-4
9.2.5	Channel Calibration	9-4
9.2.6	Initial PeakNet-PA Setup	9-8
9.3	Routine Startup and Operation	9-9
9.3.1	Routine Startup	9-9
9.3.2	Routine Operation	9-10
9.4	Short-Term Shutdown	9-11
9.5	Long-Term Shutdown	9-12
9.6	Maintenance	9-13
9.6.1	Daily Maintenance	9-13
9.6.2	Weekly Maintenance	9-14
9.6.3	Biweekly Maintenance	9-15
9.6.4	Monthly Maintenance	9-15
9.6.5	Quarterly Maintenance	9-15

10 • Troubleshooting

10.1	Troubleshooting Strategies	10-1
10.2	Liquid Leaks	10-2
10.2.1	Precision Displacement Pump Leaks . . .	10-3
10.2.2	Dilution Vessel Leaks	10-4
10.2.3	SS80 Leaks	10-5
10.2.4	EluGen Cartridge Leaks	10-5
10.2.5	Degas Assembly Leaks	10-6
10.3	Air and Gas Leaks	10-6
10.4	Excessive System Backpressure	10-7
10.5	Channel Stops Running	10-8
10.6	Module(s) Does Not Power Up	10-9
10.7	CC80 Sample LED Displays Spinning Segments	10-11
10.8	CC80 Sample LED Displays <i>EE</i>	10-11
10.9	CC80 Analyzer Leak LED Is Flashing	10-12
10.10	Precision Displacement Pump Does Not Prime	10-12
10.11	No Sample Delivered to Loading Pump	10-13
10.12	Loading Pump Delivers Inconsistent Volume	10-14
10.13	Irregularity in Loading Pump	10-15
10.14	Dilution Pump Does Not Pump	10-16
10.15	Dilution Pump Delivers Inconsistent Volume	10-17
10.16	Dilution Vessel Does Not Empty	10-18
10.17	Column Heater Does Not Heat	10-20

DX-800 Process Analyzer

10.18 Inoperative Sample Select Valve	10-20
10.19 Lack of Flow at Selected Sample Outlet	10-20
10.20 No Peaks Detected	10-21
10.21 Spurious Peaks	10-23
10.22 Poor Peak Resolution	10-24
10.23 Small Peaks Detected	10-26
10.24 Peak Height Greater Than Expected	10-26
10.25 Poor Peak Area (or Height) Precision	10-26
10.26 Nonreproducible Peak Area and/or Retention Time	10-29
10.27 Poor Retention Time Precision	10-29
10.28 Abnormal Shift in Retention Time	10-30
10.29 Poor Linear Curve	10-32
10.30 Baseline Drift	10-33
10.31 Baseline Noise—Conductivity Detection System	10-34
10.32 Baseline Noise—Absorbance Detection System	10-34
10.33 High Background—Conductivity Detection System	10-34
10.34 High Background—Absorbance Detection System	10-35
10.35 Loss of Sensitivity	10-36

11 • Service

11.1 Eliminating a Fluid System Restriction 11-2

11.2 Cleaning/Replacing Precision Displacement
Pump Check Valves 11-2

11.3 Changing a CC80 Fuse 11-5

12 • TTL and Relay Control

12.1 TTL Input Control 12-1

 12.1.1 Configuring TTL Input Actions 12-2

 12.1.2 TTL Input Signal Mode 12-3

12.2 Relay Output Control 12-4

12.3 TTL and Relay Connections 12-8

12.4 Alternate TTL Input Signal and Relay Output
Modes 12-10

A • Specifications

A.1 Electrical A-1

A.2 Environmental A-1

A.3 Physical A-2

A.4 Pumps A-2

 A.4.1 Dilution Pump A-2

 A.4.2 Loading Pump A-3

 A.4.3 PC80 Post-Column Reagent Pump
 (Optional) A-3

DX-800 Process Analyzer

A.5	Valves	A-3
A.5.1	Check Standard (CS) Valve	A-3
A.5.2	Column Switching Valve (Optional)	A-3
A.5.3	Diluent (DI) Valve	A-3
A.5.4	Diluent Select (DS) Valve	A-4
A.5.5	Dilution Vessel (DV) Valve	A-4
A.5.6	Gas Valve	A-4
A.5.7	Load/Inject (LI) Valve	A-4
A.5.8	Metering (ME) Valve	A-4
A.5.9	Sample Select (SM) Valve	A-4
A.5.10	Sample/Standard (SS) Valve	A-5
A.5.11	Standard (ST) Valve	A-5
A.6	Dilution Vessel	A-5
A.7	SS80 Sample Selector (Optional)	A-5
A.8	Column Heater (Optional)	A-6

Index

1 • Introduction

The DX-800 Process Analyzer is designed for reliability, accuracy, low maintenance, and regulatory compliance. Its IC and HPLC capabilities allow the determination of species that are not possible with other process analytical techniques. Multicomponent characterization of a sample can be performed in a single analysis, and multiple samples can be scheduled for automatic analysis.

The DX-800 is typically plumbed for one of four sample preparation configurations: concentration, concentration with reagent addition, dilution, or dilution with reagent addition. Operation is controlled by Dionex PeakNet-PA software running under Microsoft® Windows NT® 4.0 (or above).

PeakNet-PA incorporates a user interface, based on Wonderware® InTouch™, that provides an industry-standard approach to process monitoring, industrial I/O, and process control. PeakNet-PA instrument control capabilities include:

- Scheduled sampling
- Synchronized operation of multiple channels (see below)
- Automatic alarm handling, using preprogrammed conditional responses

A *channel* is a system configured to perform a specific analysis. Each channel contains an analytical pump, a detector, and a CC80 Component Controller.

An *analyzer* is defined as a sample or sampling system configured with from one to four analysis channels. The DX-800 can be multiplexed to different sample sources, using the SS80 Sample Selector.

1.1 About This Manual

Chapter 1 Introduction	Introduces the DX-800 subsystems. Explains conventions used throughout the manual. Presents safety-related information.
Chapter 2 AE80	Describes the key features of the AE80 enclosure.
Chapter 3 CC80	Describes the function and key components of the CC80 Component Controller.
Chapter 4 SP80	Describes the function and key components of the SP80 Sample Preparation panel, including the optional PC80 Post-Column Reagent Pump Kit.
Chapter 5 LC80	Describes the function and key components of the LC80 Liquid Chromatography panel.
Chapter 6 EG40-PA Eluent Generator	Describes the function and key components of the optional EG40-PA Eluent Generator.
Chapter 7 LM80	Describes key features of the LM80 Liquids Manager.
Chapter 8 SS80	Describes the function and key components of the optional SS80 Sample Selector.
Chapter 9 Operation and Maintenance	Provides start-up, operation, and shutdown instructions. Includes routine preventive maintenance procedures.
Chapter 10 Troubleshooting	Lists problems and presents step-by-step procedures for how to isolate and eliminate them.
Chapter 11 Service	Provides step-by-step instructions for routine service and parts replacement procedures.
Chapter 12 TTL/Relay Control	Describes TTL and relay control functions. Provides connection instructions.
Appendix A Specifications	Lists the DX-800 specifications and installation site requirements.

1.1.1 Safety Messages and Notes

This manual contains warnings and precautionary statements that can prevent personal injury and/or damage to the DX-800 when properly followed. Safety messages appear in bold type and are accompanied by icons, as shown below.



Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.



Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.



Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. Also used to identify a situation or practice that may seriously damage the instrument, but will not cause injury.



Indicates that the function or process of the instrument may be impaired. Operation does not constitute a hazard.

Messages d'avertissement en français



Signale une situation de danger immédiat qui, si elle n'est pas évitée, entraînera des blessures graves à mortelles.



Signale une situation de danger potentiel qui, si elle n'est pas évitée, pourrait entraîner des blessures graves à mortelles.



Signale une situation de danger potentiel qui, si elle n'est pas évitée, pourrait entraîner des blessures mineures à modérées. Également utilisé pour signaler une situation ou une pratique qui pourrait gravement endommager l'instrument mais qui n'entraînera pas de blessures.

Warnhinweise in Deutsch



Bedeutet unmittelbare Gefahr. Mißachtung kann zum Tod oder schwerwiegenden Verletzungen führen.



Bedeutet eine mögliche Gefährdung. Mißachtung kann zum Tod oder schwerwiegenden Verletzungen führen.



Bedeutet eine mögliche Gefährdung. Mißachtung kann zu kleineren oder mittelschweren Verletzungen führen. Wird auch verwendet, wenn eine Situation zu schweren Schäden am Gerät führen kann, jedoch keine Verletzungsgefahr besteht.

Notes

Informational messages also appear throughout this manual. These are labeled NOTE and are in bold type:

NOTE

NOTES call attention to certain information. They alert the user to an unexpected result of an action, suggest how to optimize instrument performance, etc.

1.1.2 Safety Labels

The DX-800 meets EN 61010-1:1993 (safety), CAN/CSA-C22.2 No. 1010.1-92 (safety), UL 3101-1/10.93 (safety), EN 50082-1:1992 (susceptibility) and EN 55011:1991 (emissions). The TUV GS, C, US Mark safety labels and the CE Mark label on the DX-800 attest to compliance with these standards.

The symbols below appear on the DX-800 or on DX-800 labels.

 Alternating current

 Protective ground conductor terminal

 Power supply is on

 Power supply is off

1.1.3 Related Manuals

The following related manuals are on the Dionex Manuals and Literature CD-ROM (P/N 053891).

- *Installation of Dionex Ferrule Fittings*
(Document No. 034213)
- *PeakNet-PA Software User's Guide*
(Document No. 031322)
- Operator's manuals for the analytical pumps, detectors and suppressors configured with the DX-800
- *EG40 Eluent Generator System Manual*
(Document No. 031373)

The following is shipped with the DX-800: NOWPAK® II documentation.

1.2 DX-800 Safety Practices

Do not operate the DX-800 for any purpose other than that for which it is designed and described in the user documentation. If there is a question regarding appropriate usage, contact Dionex before proceeding.

1.2.1 General Precautions

- Periodically check all liquid lines for leaks. Clean up spills and use deionized (DI) water to rinse dried reagents off system components.
- Make sure that gas and liquid lines cannot become kinked, punctured, or otherwise damaged.
- Do not allow liquid wastes to accumulate. Follow a regulated, approved waste disposal program. Never dispose of wastes containing organic solvents through the municipal sewage system. Neutralize all acidic and caustic wastes before disposal.

1.2.2 Compressed Gas or Liquid Cylinder Precautions

- Periodically check all pressure regulators and verify that pressure settings are within the recommended limits.
- Compressed gas cylinders are initially pressurized at 14 to 15 MPa (2200 to 2500 psi). Use a regulator to reduce the delivered air pressure to 0.3 to 0.5 MPa (50 to 75 psi).
- Fasten all cylinders securely to an immovable structure.
- Do not store or move a cylinder unless the safety cap is in place.
- Store or move cylinders in a vertical position only. Do not move the cylinders with regulators attached.
- Store cylinders in a well-ventilated area, away from heat or ignition sources.
- Clearly label each cylinder with its contents.

- Use only approved regulators and tubing connections of the appropriate material and purity.

1.2.3 Mechanical Precautions

- The piston-drive mechanism of the analytical pump contains hazardous moving parts. Before servicing, turn off the main power switch and unplug the pump.
- The precision displacement pump(s) on the SP80 Sample Preparation panel contains a piston-drive mechanism with moving parts. Before servicing, turn off the air supply to the pump and unplug the pump from the SP80 distribution board.

1.2.4 Electrical Precautions

- Replace blown fuses with the size and rating stipulated for each component.
- Verify that the selected operating voltage for the analyzer(s) is the same as the actual power line voltage.



The power supply cord is used as the main disconnect device. Make sure the outlet is located near the DX-800 and is easily accessible.



Le cordon d'alimentation principal est utilisé comme dispositif principal de débranchement. Veillez à ce que la prise de base soit située/installée près du DX-800 et facilement accessible.



Das Netzkabel ist das wichtigste Mittel zur Stromunterbrechung. Stellen Sie sicher, daß sich die Steckdose nahe am Gerät befindet und leicht zugänglich ist.

1.3 Analyzer Hardware Configuration Overview

Every DX-800 analyzer includes at least one analysis channel. Each channel is composed of several components configured inside an enclosure. Figure 1-1 depicts a typical DX-800 channel.




Figure 1-1. Typical DX-800 Channel

1.3.1 Channel Modules

This section summarizes the major functional modules of a DX-800 channel.

AE80 Enclosure

The AE80 provides the liquid and gas connections, electrical connections, electrical emissions shielding, and environmental protection necessary in a process environment.

CC80 Component Controller

The CC80 controls the valves and pumps that perform sample selection, standard preparation, sample preparation, and chromatography. The CC80 also controls the CH-4 Column Heater and PC80 Post-Column Reagent Pump, if these options are installed.

The CC80 front panel features an interactive plumbing schematic for status and local control of the SP80 and LC80 (described below), controls for selecting the sample source and column heater temperature, and eight alarm indicators.

When the CC80 power is turned on, the Moduleware and BIOS version numbers display briefly on the front panel.

The CC80 is installed at the top of the enclosure, above the analytical pump and the detector. For more information, see Chapter 3.

SP80 Sample Preparation Panel

The SP80 is purchased in one of four configurations to perform sample concentration, concentration with reagent addition, dilution, or dilution with reagent addition. All four configurations automatically prepare calibration standards. Components installed on the front of the SP80 panel include pumps, valves, and a dilution vessel.

Electronics, air valve manifolds, and gas regulators are installed on the rear of the SP80. The panel swings forward to permit access to these components.

The SP80 is located inside the enclosure, directly below the analytical pump and the detector. For more information, see Chapter 4.

LC80 Liquid Chromatography Panel

The LC80 panel provides mounting for the load/inject valve, columns, conductivity detector cell, suppressor (for conductivity detection only), and the optional column heater, column switching valve, and eluent generator cartridge.

If absorbance detection is used, the absorbance detector cell is installed inside the detector.

The LC80 distribution board is located on the back of the swing-out panel. The LC80 is secured to the inside of the enclosure door. For more information, see Chapter 5.

EG40-PA Eluent Generator (Optional)

An eluent generator can be added for generating high-purity acid or base eluents online from deionized water. The EluGen® cartridge is installed on the LC80 Liquid Chromatography panel. The controller, DX-LAN™ network interface, and other components are installed inside the AE80 enclosure. For more information, see Chapter 6.

LM80 Liquids Manager

The LM80 consists of a rack and a control panel, attached to the bottom of the channel enclosure (see Figure 1-2). A polypropylene holder in the rack accommodates two 1-liter and/or 2-liter plastic bottles for standards and reagents. The control panel includes a pressure regulator and gauge, plus valves to control the gas supply that pressurizes the reservoirs and eluent containers.

Eluents and solvents can be delivered using one or more NOWPAK II containers placed below the channel enclosure. From one to three NOWPAKs can be placed on the floor or in the polypropylene drip tray provided. For more information, see Chapter 7.




Figure 1-2. LM80 Liquids Manager

Analytical Pumps

The analytical pump is either a gradient or an isocratic pump. The pump is installed below the CC80 Component Controller.

The pump front panel display is visible through a window in the door of the DX-800 enclosure. The display backlight automatically turns off after 2 hours if no front panel buttons have been pressed. To turn it on again, press the blue **Display Refresh** button on the enclosure door or any button on the pump front panel. (The **Display Refresh** button is functional only after the PeakNet-PA Method has been loaded into the pump.)

Gradient Pump

The gradient pump is a dual processor-controlled gradient delivery system designed to blend and pump mixtures of up to four mobile phase (or eluent) components at precisely controlled flow rates. The selected mobile phase composition may be delivered isocratically, or as a linear or exponential gradient.

For operating, troubleshooting, and service information, see the pump operator's manual. Online operation of the pump is similar to the laboratory operation described in the manual.

Isocratic Pump

The isocratic pump is a dual processor-controlled isocratic delivery system designed to pump mobile phase (or eluent) at a precisely controlled flow rate.

For operating, troubleshooting, and service information, see the pump operator's manual. Online operation of the pump is similar to the laboratory operation described in the manual.

Detectors

Each channel includes a conductivity detector, absorbance detector, or electrochemical detector. The detector is installed below the CC80 Component Controller.

The detector's front panel display is visible through a window in the enclosure door. The display backlight automatically turns off after 2 hours if no front panel buttons have been pressed. To turn it on again, press the blue **Display Refresh** button on the enclosure door or any button on the detector front panel. (The **Display Refresh** button is functional only after the PeakNet-PA Method has been loaded into the detector.)

Conductivity Detector

The conductivity detector is a microprocessor-driven precision conductivity detector. The detector incorporates digital autoranging, auto offset, and temperature compensation. Optional thermal control of the conductivity cell, combined with AutoSuppression™, provides unsurpassed signal-to-noise ratios, baseline stability, and dynamic range.

For operating, troubleshooting, and service information, see the detector operator's manual. Online operation of the detector is similar to the laboratory operation described in the manual.

Absorbance Detector

The absorbance detector is a dual-beam, variable-wavelength photometer. Advanced fiber-optics technology provides simple and reliable operation. Full spectral capability is provided by two light sources: a deuterium lamp for ultraviolet detection and a tungsten lamp for visible wavelength operation.

For operating, troubleshooting, and service information, see the detector operator's manual. Online operation of

the detector is similar to the laboratory operation described in the manual.

Electrochemical Detector

The electrochemical detector provides the three major forms of electrochemical detection: conductivity, DC amperometry, and integrated amperometry. An additional mode is voltammetry, which can be used to determine the potentials used for amperometry operation.

For operating, troubleshooting, and service information, see the detector operator's manual. Online operation of the detector is similar to the laboratory operation described in the manual.

1.3.2 SS80 Sample Selector (Optional)

The SS80 Sample Selector can select one of 7, 14, or 21 sample sources, depending on the number of valves inside the module. For a single sample, no SS80 is required.

Sample lines can be continuously flushing or static, and can be returned to the process or to waste. The SS80 is in a separate enclosure to isolate the bulk of the liquid flow from the analytical instrumentation. For more information, see Chapter 8.

1.3.3 Accessories

- Column heater. Installing a CH-4 Column Heater (P/N 051890) on the LC80 Liquid Chromatography panel ensures stable column temperatures for applications that are temperature-sensitive or require elevated temperatures. The heater holds one column with these dimensions: 6 or 8 mm OD x 100, 150, or 250 mm ID. The column heater temperature range is from ambient + 5 °C to 80 °C ± 1°C.
- PC80 Post-Column Reagent Pump. The PC80 delivers reagent for post-column detection methods. The PC80 kit, available in 115 V (P/N 050305) and 230 V (P/N 050307) versions, includes a single-piston pump, knitted reaction coil, and pulse damper. The PC80 is installed at the base of the enclosure behind the SP80 Sample Preparation panel.
- Column switching valve. A 10-port, 2-position, electrically-actuated PEEK valve (P/N 051824) can be added to the LC80 Liquid Chromatography panel for control of the liquid flow in column switching methods.

Contact Dionex for more information about DX-800 options.

1.4 PeakNet-PA Software

PeakNet-PA software provides automated control of DX-800 instrumentation, data acquisition and analysis, and results reporting. PeakNet-PA can control up to four analyzers, provided there are no more than eight total channels. A maximum of four channels of instrumentation can be configured to a single analyzer.

PeakNet-PA is shipped from the factory already installed on a workstation. Communication between the workstation and the DX-800 is via the Dionex DX-LAN interface.

For operating instructions, refer to the *PeakNet-PA Software User's Guide* (Document No. 031322) or the PeakNet-PA online Help.

2 • AE80 Enclosure

The standard AE80 enclosure is configured as a NEMA 12 enclosure. NEMA 4 and 4X enclosures are available as options. NEMA enclosures are constructed of 14-gauge steel, with all seams continuously welded and ground smooth. These enclosures are designed for use in nonhazardous locations, as indicated in the table below.

NOTE

The information in this table is not intended to be a complete representation of NEMA standards for enclosures nor those of the Electrical and Electronic Manufacturers Association of Canada (EEMAC).

Enclosure Type	Location
NEMA 12	Intended for use primarily to provide a degree of protection against dust, falling dirt, and dripping noncorrosive liquids.
NEMA 4	Intended for use primarily to provide a degree of protection against windblown dust and rain, splashing water, and hose-directed water.
NEMA 4X	Intended for use primarily to provide a degree of protection against corrosion, windblown dust and rain, splashing water, and hose-directed water.

The external dimensions of the enclosure (without an air conditioner) are 93 x 58 x 58 cm (37 x 23 x 23 in). Allow at least 1 m (40 in) of clearance in front of the enclosure for opening the door and performing maintenance tasks. The 5-cm (2-in) flanges on the top and bottom are required to mount the enclosure on the wall.

IMPORTANT

The user is responsible for ensuring the structural integrity of the installation site and for supplying mounting hardware for the enclosure.

DX-800 Process Analyzer

2.1 Front Door

- The front panel displays on the analytical pump and the detector are visible through two small windows in the enclosure door. The display brightness on the pump or detector modules can be adjusted from the **MODULE SETUP** screen (a **MAIN MENU** option).
- The red **Emergency Off** switch controls power to all components except for the intrinsically-safe power circuit routed to the front door.

IMPORTANT

If an emergency occurs, turn off the power by pushing the Emergency Off switch completely in. After resolving the situation, push the Power Reset button to restore power.




Figure 2-1. AE80 Enclosure Door

- Push the **Power Reset** button to start up the DX-800 system when the power is turned on initially or to restore power after an emergency shutdown. As long as the power is on and the system is ready to operate, the **Power Reset** lamp remains illuminated continuously.
- The blue **Display Refresh** button controls the LCD backlights on the analytical pump and the detector. These displays automatically turn off after 2 hours if no front panel buttons have been pushed; to turn on the backlights again, press **Display Refresh**. (The **Display Refresh** button is functional only after a PeakNet-PA Method has been downloaded to the modules.)
- To open the enclosure door, lift up on the door handle and turn it one-quarter turn to the left. To close the door, turn the handle to the right and push in.



SHOCK HAZARD—A shock hazard exists inside the enclosure when the door is opened.



Various types of chemicals are used in the DX-800, depending on the application that is being performed. Follow all appropriate hazardous materials and safety guidelines for chemicals when operating the DX-800.



DANGER D'ÉLECTROCUTION—Un danger d'électrocution existe dans l'enceinte lorsque la porte est ouverte.'



Différents types de produits chimiques sont utilisés dans le DX-800, selon l'application à effectuer. Respectez toutes les directives de sécurité sur les matières dangereuses pour les produits chimiques lors de l'utilisation du DX-800.



STROMSCHLAGGEFAHR—Bei geöffneter Tür besteht im Gehäuseinnern Gefahr durch elektrischen Schlag.


DX-800 Process Analyzer



Je nach Anwendung, die gerade läuft, werden im DX-800 verschiedenartige Chemikalien verwendet. Beachten Sie beim Betrieb des DX-800 alle entsprechenden Sicherheitsrichtlinien bezüglich gefährlicher Stoffe für die verwendeten Chemikalien.

2.2 Electrical System

Electrical connections are made to an electrical I/O panel on top of the AE80 enclosure. There are two versions of electrical panels: one with BNC DX-LAN connectors (see Figure 2-2) and one with an RJ-45 (10BASE-T) DX-LAN connector (see Figure 2-3).



*Figure 2-2. Top View of AE80 Enclosure
Electrical I/O Panel with BNC Connectors*




Figure 2-3. Top View of AE80 Enclosure Electrical I/O Panel with RJ-45 Connector

- The power in connector provides a connection to the main power. Connect the modular power cord (IEC 320 C13) from here to a grounded, single-phase power source.



SHOCK HAZARD—If a grounded receptacle is not used, a shock hazard may result. Do not operate or connect to AC power mains without earthed ground connections.



The power cord is used as the main disconnect device. Make sure the outlet is located near the enclosure and is easily accessible.



DANGER D'ÉLECTROCUTION—Pour éviter toute électrocution, il faut utiliser une prise de courant avec prise de terre. Ne l'utilisez pas et ne le branchez pas au secteur C.A. sans utiliser de branchement mis à la terre.



Le cordon d'alimentation principal est utilisé comme dispositif principal de débranchement. Veillez à ce que la prise de base soit située/installée près du DX-800 et facilement accessible.



STROMSCHLAGGEFAHR—Zur Vermeidung von elektrischen Schlägen ist eine geerdete Steckdose zu verwenden. Das Gerät darf nicht ohne Erdung betrieben bzw. an Wechselstrom angeschlossen werden.



Das Netzkabel ist das wichtigste Mittel zur Stromunterbrechung. Stellen Sie sicher, daß sich die Steckdose nahe am Gerät befindet und leicht zugänglich ist.

- The TTL/relay connector provides a passage for TTL and relay cables connected to the CC80 rear panel.
- The SS80 connector provides a connection from the analyzer to an SS80 Sample Selector. Connect a 9-pin shielded cable (P/N 050363) from here to the connector on the side of the SS80.
- The DX-LAN connector(s) provides connections to the host computer of the PeakNet-PA workstation.

2.3 Fluid and Pneumatic Systems

2.3.1 External Features

All liquid and gas lines are connected to the plumbing interface panel on the bottom of the AE80 enclosure, behind the LM80 control panel (see Figure 1-2).

2.3.2 Internal Features

A drip tray in the bottom of the AE80 enclosure traps liquid spills and leaks. Trapped liquid normally exits the enclosure via a drain tube in the drip tray.

A sensor in the drip tray reports the presence of liquid to PeakNet-PA and the CC80 electronics. When a leak is reported, the **Analyzer Leak** LED on the CC80 front panel begins flashing. Respond to leak alarms promptly in order to prevent overflows.

2.4 Environmental Controls

Each channel is equipped with a ventilation blower, an air conditioner, or a purge unit, depending on the ambient conditions of the installation site.

- A blower is installed when the ambient temperature range is 4 to 40 °C (40 to 105 °F).
- An optional air conditioner will allow operation, as a sealed enclosure, in ambient temperatures up to 55 °C (130 °F).

Ventilation Blower

The blower dissipates heat by drawing in air through a grill and filter, routing it to the top of the AE80 enclosure, and expelling it through a vent in the top of the enclosure. The blower operates continuously while the enclosure power is on.

Two air filters (one on each side, near the bottom of the enclosure) prevent the blower from pulling dirt and dust into the enclosure. Clean the filters with warm water whenever a fine layer of dust or lint is visible on the surface. Establish a cleaning schedule to meet this requirement, taking local air quality into account.

Air Conditioner (Optional)

The air conditioner is installed on the side of the enclosure. The air conditioner cools to below ambient and keeps the enclosure isolated to prevent internal components from coming in contact with damaging airborne materials. Periodically remove the filter from the bottom of the air conditioner and clean it.

Purge-and-Pressurization Unit (Optional)

When the installation site is a Class 1, Division 2 location requiring intrinsically safe equipment, a purge-and-pressurization unit is installed on top of the AE80 enclosure.

The purge-and-pressurization unit maintains a positive pressure inside the enclosure. The over-pressure relief valve and orifice are located on the side of the unit.

DX-800 Process Analyzer

3 • CC80 Component Controller

3.1 Overview

The CC80 Component Controller (CC80) is installed at the top of the AE80 enclosure. The CC80 provides control and status for these components:

- Valves and pumps on the SP80 Sample Preparation panel and LC80 Liquid Chromatography panel
- SS80 Sample Selector valves (optional)
- CH-4 Column Heater (optional)
- PC80 Post-Column Reagent Pump (optional)

The center of the CC80 front panel is an interactive flow chart that represents the system components. Besides providing a convenient overview of the plumbing configuration, the flow chart buttons and indicator lights serve these purposes:

- Display the real-time state, or operating position, of valves and pumps. An indicator light beside each button indicates the status.
- Enable direct control of any component on the SP80 and LC80 panels, or the SS80 valves.
- Initiate a drain sequence for the dilution vessel if a method is ended or aborted before the dilution vessel is emptied.

Finally, certain alarm conditions are reported to the CC80. Four hardware alarm conditions are preassigned, while the user can define four additional alarms in PeakNet-PA software. When an alarm condition occurs, the corresponding annunciator on the CC80 front panel begins flashing.

3.2 Operating Features

3.2.1 Front Panel

There are four versions of the CC80 front panel, corresponding to the four SP80 Sample Preparation configurations:

- CC81/SP81 for concentration or direct injection
- CC82/SP82 for dilution or direct injection
- CC83/SP83 for dilution with reagent addition
- CC84/SP84 for concentration with reagent addition

Figure 3-1 illustrates the CC81/SP81 front panel for concentration. For illustrations of the other front panels, see Chapter 4.




Figure 3-1. CC81/SP81 Front Panel (for Concentration)

3 • CC80 Component Controller


Button	Function
	<p>Toggles between Local and Remote control modes. In Local mode, components displayed in the flow chart are controlled via the CC80 front panel. In Remote mode, components displayed in the flow chart are controlled from PeakNet-PA. The light indicates the selected mode. The CC80 is in Local mode when powered up.</p>
	<p>Provides on/off control of the CC80 and components controlled by the CC80.</p>
<p data-bbox="402 888 513 919">Sample</p> 	<p>These buttons are functional only when the SS80 Sample Selector is installed. The buttons select SS80 valve positions that correspond to specific samples. If no SS80 is installed, only one sample is being supplied to the analyzer.</p> <p>Press a Select button to increase (or decrease) the sample number. (Numbers blink while an entry is in progress.) The selected number is displayed above the buttons. To confirm the selected sample number, press Enter.</p> <p>When a sample select valve is switching positions, two of the LED segments rotate in a clockwise direction. If no button in this section is pressed, the display reverts to the currently selected sample number.</p>
<p>FLOW CHART</p>	<p>The buttons in the flow chart vary, depending on the SP80 configuration. For descriptions of all four versions, see Chapter 4.</p>

Table 3-1. CC80 Front Panel Button Functions

DX-800 Process Analyzer


Button	Function
	<p><i>If no alarm condition exists, annunciator lights in this section are not illuminated.</i></p> <p>When a problem is detected, an annunciator that identifies the source of the problem begins flashing. Press Alarm Reset to turn off the annunciator. If the problem has not been fixed, the annunciator will start flashing again.</p> <p>The preassigned Analyzer Leak and Sampling Leak annunciators indicate leaks in the AE80 enclosure and the SS80 Sample Selector, respectively. The preassigned Loading Pump and Dilution Pump annunciators indicate pump failure due to cavitation or starvation of the pump. Activation of Alarms 1-4 can be assigned in PeakNet-PA by the user.</p>
	<p>Sets the temperature of the column heater (optional) from ambient + 5 °C to 80 °C ± 1 °C. Press a Select button to increase or decrease the temperature by one degree. (Numbers blink while an entry is in progress.) The selected temperature is displayed in °C above the buttons.</p> <p>To confirm the selected temperature, press Enter. If the column heater has not reached the selected temperature, Temp Not Ready flashes.</p>

Table 3-1. CC80 Front Panel Button Functions (continued)

3.2.2 Rear and Side Panels

- Five AC outlets provide connections to the analytical pump, detector, column heater, and fan, cooler, or air conditioner.
- The main power receptacle provides a connection to the main power inlet on top of the AE80 enclosure (see Figures 2-2 and 2-3).
- Electrical cables, including the DX-LAN cable and the TTL/Relay cable, exit the CC80 via openings in the rear panel. The cables are connected to the top of the enclosure. See Chapter 12 for details about using TTLs and relays.
- A 2-pole, 20 amp circuit breaker is installed on the left side panel of the CC80. If the breaker is tripped, reset it by flipping up the switch.

3.2.3 DIP Switches

The CC80 contains two sets of DIP switches. The CC80 reads the state of the switches each time the power is turned on. The switches are typically set at the factory. If an option is installed in the field, the appropriate DIP switch must be reset.

DIP switch #1 records the following information:

- The number of valves in the SS80
- The configuration of the SP80
- The options that are installed in the DX-800

DIP switch #2 controls the following functions:

- Enabling or disabling an audible tone when an alarm condition is reported to the CC80
- Setting the TTL input and relay output modes (see Section 12.4)
- Enabling or disabling the metering (ME) valve functions

Follow this procedure to reset the DIP switches:

1. Press the **Emergency Off** button on the front door of the AE80 enclosure to turn off the main power. Disconnect the power cord.



SHOCK HAZARD—Disconnect the power cord before raising the top cover of the CC80.



DANGER D'ÉLECTROCUTION—Débranchez le cordon d'alimentation électrique avant de soulever le couvercle du CC80.



STROMSCHLAGGEFAHR—Ziehen Sie den Netzstecker, bevor Sie den Deckel des CC80 öffnen.

3 • CC80 Component Controller

2. Open the front door of the AE80 enclosure. Loosen the two retaining screws in the module chassis. Pull out the module chassis just until the first stop on the sliders is engaged.



Do not pull the chassis beyond the slider stops. The chassis can become disengaged from the enclosure.



Ne tirez pas le châssis au-delà des butées du mécanisme coulissant. Le châssis peut se détacher de l'enceinte.



Ziehen Sie die Montageplatte nicht über die Arretierungen hinaus. Sie könnte sich sonst vom Gehäuse lösen.

3. Loosen, but do not remove, the four screws securing the CC80 top cover in place.
4. Raise the top cover and slide it back to expose the DIP switches on the back of the front panel display board.
5. Check Figure 3-2, and then set the switches as required. The switches are numbered 1 through 8. The *off* position (logic-false) is down and the *on* position (logic-true) is up.
6. Replace the CC80 top cover. Push the module chassis back into the enclosure and tighten the retaining screws.
7. Close the enclosure door.
8. Plug in the power cord and press the **Power Reset** button on the front door of the AE80 to turn on the power.

DX-800 Process Analyzer




Figure 3-2. DIP Switch Settings

3.3 Operating Modes

The CC80 has two operating modes:

- In Local mode, operation is controlled from the CC80 front panel.
- In Remote mode, operation is controlled from PeakNet-PA software.

3.3.1 Local Mode

When the CC80 power is turned on, it is in Local mode. In Local mode, the CC80 accepts operating commands from the front panel buttons. This allows direct control of individual valves and pumps, which is useful for troubleshooting and service procedures.

When the CC80 is in Remote mode, press the **Local/Remote** button on the CC80 front panel to return to Local mode. If a PeakNet-PA Method is running at the time, the CC80 automatically aborts the Method. The light indicates which mode is selected.

3.3.2 Remote Mode

In Remote mode, PeakNet-PA sends operating commands from the PeakNet-PA workstation to the CC80 via the DX-LAN interface. The CC80 is automatically switched to Remote mode when a Schedule is downloaded from PeakNet-PA or when a Method is downloaded as a result of running a manual sample.

In Remote mode, sample preparation and chromatography Methods are executed. The CC80 performs the sample preparation Method, and sends a run command to the analytical pump and the detector. PeakNet-PA is responsible for downloading Methods and controlling the start of Method execution. To run a single Method, run a manual sample from the Analyze program.

3.4 Power-Up Diagnostics

Each time the CC80 power is turned on or reset (which occurs when new Moduleware is downloaded), this sequence of events occurs:

- All CC80 front panel LEDs are illuminated for 2 seconds.
- The BIOS version number is displayed for 4 seconds, with the integer in the **Sample** section and the decimal fraction in the **Column Temperature** section.
- The Moduleware version number is displayed for 4 seconds, with the integer in the **Sample** section and the decimal fraction in the **Column Temperature** section.
- The CC80 automatically begins running a series of diagnostic tests (see below).

If the CC80 fails a diagnostic test, contact Dionex. In the U.S., call Dionex Technical Support at 1-800-346-6390. Outside the U.S., call the nearest Dionex office.

Diagnostic Test	Purpose
CPU Test	Test the CPU's internal configuration and checksum.
Load/Inject Valve Test	Verify that the valve can be set to the load and inject positions.
Metering Valve Test	Verify that the ME valve can be set to both operating positions. If the ME valve or SP80 panel is not installed, this test must be disabled with DIP switch 2, position 5 (see Figure 3-2).
LED Test	Turn on all LEDs and 7-segment display segments for 2 seconds.

4 • SP80 Sample Preparation Panel

4.1 Overview

The SP80 Sample Preparation panel is equipped with the pump(s), valves, dilution vessel, and interconnecting tubing required for one of these four configurations:

- Concentration or direct injection (SP81)
- Dilution or direct injection (SP82)
- Dilution with reagent addition (SP83)
- Concentration with reagent addition (SP84)

While the DX-800 can automatically prepare and analyze calibration standards in all four configurations, the distinguishing feature of each is the type of sample preparation performed. For a summary of the components required for each configuration, see Sections 4-2 through 4-5.

The SP80 panel is located in the AE80 enclosure, below the analytical pump and the detector (see Figure 1-1). Components are located on the front and rear of the swing-out panel. To access the rear, loosen the thumbscrew on the right side of the panel and swing it open.

4.2 Configuration: SP81 for Concentration or Direct Injection

This configuration is used for trace ion analysis or direct injection, in which the sample is pre-concentrated on a concentrator column, or pumped to the sample loop, before the chromatographic analysis.

All materials in the flow paths are of the highest purity and permit trace ion analysis even at the low part-per-trillion (ppt) level.

Figure 4-1 shows the CC80 front panel flow chart for the concentration configuration.




Figure 4-1. CC81 Front Panel: SP81 for Concentration or Direct Injection

4 • SP80 Sample Preparation Panel

Components for Concentration		
Component	Description	Valve Position and/or Function*
Check Standard (CS) Valve	3-way solenoid valve	Selects sample (0) or check standard (1).
Sample (SM) Valve	3-way solenoid valve	Directs sample or check standard to waste (0) or to SS valve and loading pump for analysis (1).
Sample/Standard (SS) Valve	3-way solenoid valve	Selects sample (0) or calibration standard (1).
Dilution Vessel (DV) Valve	3-way solenoid valve	Purges dilution vessel to waste (0) or directs calibration standard to SS valve (1).
Diluent (DI) Valve	2-way solenoid valve	Provides on/off control of diluent (typically DI water) to dilution pump.
Metering (ME) Valve	Rheodyne 10-port valve	Measures stock standard for delivery to dilution vessel (for preparation of calibration standard).
Gas Valve	3-way solenoid valve	Vents (0) or pressurizes (1) dilution vessel.
Dilution Pump	Precision displacement pump	Delivers diluent to dilution vessel.
Loading Pump	Precision displacement pump	Loads samples or standards on concentrator column.
Regulator/Valve Manifold (on rear of SP80)	4-way air valve manifold	Regulates air flow to dilution and loading pumps.
Regulator/Valve Manifold (on rear of SP80)	3-way gas valve manifold	Regulates flow of high purity gas (typically helium) to dilution vessel.
* 0 (the default) is normally open, or non-energized. 1 is normally closed, or energized.		

4.3 Configuration: SP82 for Dilution

This configuration is used for applications in which the samples are injected directly into the chromatographic system or for applications in which the samples must be diluted prior to chromatographic analysis. Dilution factors up to 1/25,000 can be achieved. A channel configured for dilution is often used when assaying process samples for major constituents.

Figure 4-2 shows the CC80 front panel flow chart for the dilution configuration.




Figure 4-2. CC82 Front Panel: SP82 for Dilution or Direct Injection

4 • SP80 Sample Preparation Panel

Components for Dilution		
Component	Description	Valve Position and/or Function*
Check Standard (CS) Valve	3-way solenoid valve	Selects sample (0) or check standard (1).
Sample (SM) Valve	3-way solenoid valve	Directs sample to SS valve (0) or to ST valve (1).
Standard (ST) Valve	3-way solenoid valve	Selects stock standard for calibration standard preparation (0) or sample for dilution (1).
Sample/Standard (SS) Valve	3-way solenoid valve	Selects undiluted sample (0) or diluted sample (1).
Dilution Vessel (DV) Valve	3-way solenoid valve	Purges dilution vessel to waste (0) or directs diluted sample/standard to SS valve (1).
Diluent (DI) Valve	2-way solenoid valve	Provides on/off control of diluent (typically DI water) to dilution pump.
Metering (ME) Valve	Rheodyne 10-port valve	Measures sample or standard for delivery to dilution vessel.
Gas Valve	3-way solenoid valve	Vents (0) or pressurizes (1) dilution vessel.
Dilution Pump	Precision displacement pump	Delivers diluent to dilution vessel.
Regulator/Valve Manifold (on rear of SP80)	4-way air valve manifold	Regulates air flow to dilution and loading pumps.
Regulator/Valve Manifold (on rear of SP80)	3-way gas valve manifold	Regulates flow of high purity gas (typically helium) to dilution vessel.
* 0 (the default) is normally open, or non-energized. 1 is normally closed, or energized.		

4.4 Configuration: SP83 for Dilution with Reagent Addition

This configuration is used when the samples require dilution and reagent addition prior to chromatographic analysis.

Figure 4-3 shows the CC80 front panel flow chart for the dilution with reagent addition configuration.




Figure 4-3. CC83 Front Panel: SP83 for Dilution with Reagent Addition

4 • SP80 Sample Preparation Panel

Components for Dilution with Reagent Addition		
Component	Description	Valve Position and/or Function*
Sample (SM) Valve	3-way solenoid valve	Directs sample to ST valve (0) or to SS valve (1).
Standard (ST) Valve	3-way solenoid valve	Selects sample (0) or stock standard (1) for dilution.
Sample/Standard (SS) Valve	3-way solenoid valve	Selects undiluted sample (0) or diluted sample (1).
Dilution Vessel (DV) Valve	3-way solenoid valve	Purges dilution vessel to waste (0) or directs diluted sample/standard to SS valve (1).
Diluent Select (DS) Valve	3-way solenoid valve	Selects diluent (0) or reagent (1).
Diluent (DI) Valve	2-way solenoid valve	Provides on/off control of diluent (typically DI water) to dilution pump.
Metering (ME) Valve	Rheodyne 10-port valve	Measures sample or standard for delivery to dilution vessel.
Gas Valve	3-way gas/liquid solenoid valve	Vents (0) or pressurizes (1) dilution vessel.
Dilution Pump	Precision displacement pump	Delivers diluent to dilution vessel.
Regulator/Valve Manifold (on rear of SP80)	4-way air valve manifold	Regulates air flow to dilution and loading pumps.
Regulator/Valve Manifold (on rear of SP80)	3-way gas valve manifold	Regulates flow of high purity gas (typically helium) to dilution vessel.
* 0 (the default) is normally open, or non-energized. 1 is normally closed, or energized.		

4.5 Configuration: SP84 for Concentration with Reagent Addition

In addition to performing concentration (see Section 4.2), this configuration is used to add reagent to samples prior to analysis. A typical application is the acidification of samples prior to chromatographic analysis for trace transition metals.

Figure 4-4 shows the CC80 front panel flow chart for the concentration with reagent addition configuration.




Figure 4-4. CC84 Front Panel: SP84 for Concentration with Reagent Addition

4 • SP80 Sample Preparation Panel

Components for Concentration with Reagent Addition		
Component	Description	Valve Position and/or Function*
Sample (SM) Valve	3-way solenoid valve	Directs sample (0) or check standard (1) to DS valve.
Diluent Select (DS) Valve	3-way solenoid valve	Selects diluent (0) or sample/check standard (1).
Diluent (DI) Valve	2-way solenoid valve	Provides on/off control of diluent (typically DI water) to dilution pump.
Standard (ST) Valve (on rear of SP80)	3-way solenoid valve	Selects reagent (0) or stock standard (1).
Metering (ME) Valve	Rheodyne 10-port valve	Measures stock standard or reagent for delivery to dilution vessel.
Gas Valve (on rear of SP80)	3-way solenoid valve	Vents (0) or pressurizes (1) dilution vessel.
Sample/Standard (SS) Valve	3-way solenoid valve	Selects no flow (0) or selects sample or standard from dilution vessel (1).
Dilution Vessel (DV) Valve	3-way solenoid valve	Purges dilution vessel to waste (0) or directs prepared sample or standard to SS valve (1).
Dilution Pump	Precision displacement pump	Delivers diluent to dilution vessel.
Loading Pump	Precision displacement pump	Loads samples or calibration standards on concentrator column.
Regulator/Valve Manifold (on rear of SP80)	4-way air valve manifold	Regulates air flow to dilution and loading pumps.
Regulator/Valve Manifold (on rear of SP80)	3-way gas valve manifold	Regulates flow of high purity gas (typically helium) to dilution vessel.
* 0 (the default) is normally open, or non-energized. 1 is normally closed, or energized.		

4.6 Precision Displacement Pumps

The SP80 contains one or two precision displacement pumps, depending on the configuration.

- The loading pump is required for loading sample onto a concentrator column. The pump is present in the SP80 configurations for concentration and concentration with reagent. The loading pump is not used when a sample loop is connected to the load/inject (LI) valve.
- The dilution pump is required for pumping diluent (typically deionized water) through the metering (ME) valve to the dilution vessel.

A precision displacement pump is an air-driven, volume displacement pump. Each pump is controlled by two three-way air valves on the rear of the SP80 panel. Air pressure at 275 kPa (40 psi) is applied to two chambers on either end of a cylindrical block. When the air pressure is sufficient to overcome the liquid pressure on the other side of the piston, the piston pulls or pushes liquid through the pump head.

Two check valves control the direction of liquid flow. With each pump stroke, the piston draws liquid into the pump. When the piston reaches maximum travel, the valve states are reversed, so that the piston moves back and pushes out the liquid. When the piston reaches maximum travel in the other direction, the pump is ready to perform the next stroke. A small rod on the piston trips optical sensors on the pump controller board when the piston has reached maximum travel in either direction.

With a fixed driving air pressure, the liquid flow rate will decrease with increasing backpressure. This is advantageous, given that the desired flow rate is 10 to 15 mL/min for the sample dilution pump and less than 3 mL/min for the sample loading pump.

The PeakNet-PA Configuration Editor includes calibration Methods for the precision displacement pumps. To begin, select the Calibration Wizard command from the Configure menu and follow the step-by-step calibration instructions.

4.7 Dilution Vessel

The dilution vessel is used to prepare calibration standards and to dilute samples (in the dilution configuration). The 250 mL dilution vessel is made of high purity polypropylene (PP), with an option for perfluoroalkoxy (PFA) Teflon®.

Pressurize the dilution vessel with high purity nitrogen or helium (filtered, dry, and oil-free) regulated to 170 to 240 kPa (25 to 35 psi). A pressure relief valve on the rear of the SP80 is designed to open if this reaches 340 kPa (50 psi).

4.8 Post-Column Reagent Pump (Optional)

The PC80 Post-Column Reagent Pump Kit is available in both a 115 V version (P/N 050305) and a 230 V version (P/N 050307). The kit includes the following items:

- Dionex Reagent Pump (RP-1)
- Pulse damper
- Knitted reaction coil

The post-column components are installed in the area behind the SP80 Sample Preparation panel. The PC80 pump is turned on with either a CC80 Method setup parameter (Remote mode) or from the CC80 front panel (Local mode).

The pump power is also automatically turned off if the analytical pump encounters a pressure limit alarm and shuts off.

Post-Column Reagents

The LM80 Liquids Manager holds two plastic bottles for standards or reagents. Two-liter plastic bottles (P/N 044129) are shipped with the LM80; 1-liter bottles (P/N 044128) are also available.

When preparing post-column reagents, use only ASTM Type II (18.0 MΩ/cm resistance or 1 μS) deionized water and reagents of adequate purity. Improperly prepared reagent is a common cause of baseline drift, high background, and column contamination.

DX-800 Process Analyzer

5 • LC80 Liquid Chromatography Panel

5.1 Overview

The LC80 Liquid Chromatography panel is equipped with the channel chromatography components—the load/inject (LI) valve, the columns, and (for conductivity detection) the suppressor and conductivity cell. For channels using UV/Vis absorbance detection, the absorbance detector cell is inside the detector. Optional chromatography components, such as the CH-4 column heater or EluGen Cartridge, are also mounted on the LC80.

The LC80 panel is located on the inside of the enclosure door. For servicing, the panel is hinged to permit easy access to the distribution board on the rear of the panel. Electrical connections for the LC80 components are made to the distribution board.

5.2 Chromatography Components

Chromatography hardware components, such as the load/inject (LI) valve and conductivity cell (if required), are installed on the LC80 panel at the factory. Chromatography consumables, including columns and suppressors, are installed during the analyzer installation and setup.

The numerous mounting holes and slots on the LC80 panel maximize the installation options. Figure 5-1 illustrates an example component layout. Many other component configurations are possible.




Figure 5-1. Example LC80 Component Layout

5.2.1 Load/Inject (LI) Valve

The load/inject valve is a six-port, two-position, electrically-actuated Rheodyne valve. The valve can be configured with a sample loop (10 to 1000 μ L) or a concentrator column.

5.2.2 Detector Cell

For applications requiring conductivity detection, the conductivity cell for the conductivity detector or electrochemical detector is housed in a DS3 Detection Stabilizer, which is mounted on the LC80 panel.

The DS3 improves baseline stability by heating the cell to a selectable set point above ambient and preventing temperature fluctuations.

5.2.3 Consumables

Guard and Analytical Columns

The LC80 has several sets of column clips for mounting guard and analytical columns so as to minimize the total extra-column volume.

Concentrator Columns

Concentrator columns are used for trace analysis to enrich analyte concentration.

Suppressor

The suppressor is used with conductivity detection to neutralize the eluent and enhance analyte conductivity.

Use the manual regenerant shutoff valve at the bottom of the panel to turn off the flow of regenerant (typically deionized water) to the suppressor whenever the suppressor is being changed or whenever eluent flow from the analytical pump has been turned off. When the regenerant flow is stopped, the

DX-800 Process Analyzer

power to the suppressor must also be turned off. For detailed operating instructions, refer to the suppressor manual.

EluGen Cartridge

If the DX-800 is configured with the optional EG40-PA Eluent Generator, an EluGen cartridge is mounted on the LC80 panel. The mounting bracket is supplied with the EG40-PA.

5.2.4 CH-4 Column Heater (Optional)

The CH-4 column heater (P/N 051890) accommodates one 6 or 8 mm OD x 100, 150, or 250 mm ID column. The heater operates between ambient + 5 °C to 80 °C and is used for temperature-sensitive chromatographic methods. Select the heater temperature from the CC80 front panel or from PeakNet-PA software.

5.2.5 Column Switching Valve (Optional)

A 10-port, two-position, electrically-actuated rotary valve (P/N 051824) is available for use with column switching chromatography methods. The high pressure valve is controlled by an output from the analytical pump.

6 • EG40-PA Eluent Generator

6.1 Overview

The EG40-PA Eluent Generator is an optional device that generates high-purity acid or base eluents online, using only deionized water as the carrier. The device consists of the EG40-PA electronics components, a degas tubing assembly, and either an EGC-KOH EluGen Cartridge or an EGC-MSA EluGen Cartridge.

The DX-800 analytical pump delivers deionized water to the EluGen cartridge, which generates the eluent. The eluent exits the cartridge and flows through a degas tubing assembly that removes electrolysis gases created during eluent generation. After degassing, the eluent flows to the inject valve. Figure 6-1 diagrams the basic system flow.

For detailed system flow information for different applications, refer to the eluent generator manual.




Figure 6-1. EG40-PA System Flow Diagram

6.2 Main Components

The eluent generator consists of the following main components:

- The disposable EluGen cartridge is mounted on the LC80 Liquid Chromatography Panel (see Figure 6-2).
- A vent tube exits the side of the EluGen cartridge and is routed out the AE80 through the plumbing I/O panel. See the caution on the next page for venting instructions.
- A high-pressure degas tubing assembly is plumbed between the EluGen cartridge and the inject valve. The degas assembly is on the rear of the LC80 panel. A system waste line from the assembly is routed out the AE80 through the plumbing I/O panel. See the caution on the next page for waste disposal and venting instructions.



Figure 6-2. DX-800 with EG40-PA Installed



The EG40-PA generates eluent by means of electrolysis, which results in the production of small amounts of oxygen or hydrogen gas. To ensure that the gas is not trapped in a closed container and allowed to concentrate, install a 1.3-cm (0.52-in) ID black gas separator waste tube (P/N 045460) to an uncapped waste reservoir. See the suppressor manual for installation instructions. Direct the waste line from the EG40-PA degas assembly to the waste tube. In addition, direct the clear vent tubing from the gas vent port on the EluGen cartridge to the gas separator waste tube or to any open, well-ventilated location.



Le EG40-PA produit des éluants par électrolyse, résultant en la production de petites quantités de gaz d'oxygène et d'hydrogène. Pour veiller à ce que les gaz ne soient pas emprisonnés dans un contenant fermé où ils pourraient se concentrer, installez un tube d'évacuation du séparateur de gaz noir ID (réf. 045460) de 1,3 cm (0,52 po) dans un réservoir d'évacuation ouvert (non bouché). Consultez le manuel du dispositif de suppression pour obtenir des instructions d'installation. Orientez la conduite d'évacuation de l'ensemble de dégazage de l'EG40-PA vers le tube d'évacuation. De plus, dirigez le tuyau d'aération clair de l'évent des gaz sur la cartouche EluGen vers le tube d'évacuation du séparateur de gaz ou vers n'importe quel lieu découvert et bien aéré.



Der EG40-PA erzeugt Eluenten durch Elektrolyse. Dabei entstehen kleine Mengen von Sauerstoff- und Wasserstoffgas. Verbinden Sie einen Gasabscheiderschlauch (ID = 1,3 cm) Bestell-Nr. 045460) mit einem offenen (unverschlossenen) Abgasbehälter, damit sich kein Gas in einem geschlossenen Behälter sammelt und aufkonzentriert. Hinweise zur Installation finden Sie im Suppressor-Handbuch. Führen Sie die Abgasleitung von der Entgasungseinheit des EG40-PA zum Abgasschlauch. Führen Sie außerdem den Entlüftungsschlauch vom Entlüftungspport der EluGen-Cartridge zum Abgasschlauch des Gasabscheiders oder zu einer beliebigen anderen, gut belüfteten Stelle.

DX-800 Process Analyzer

- An EG40-PA controller is installed inside the AE80 at the lower rear of the enclosure (see Figure 6-2). The main power switch for the eluent generator is on the controller. The power LED is illuminated when the power is on. See the eluent generator manual for details about the other LEDs on the controller.
- TTL input and output connectors are present, but are not typically used with the DX-800.
- DX-LAN interface network connections are below the controller. The DX-LAN network allows communication between the eluent generator and the PeakNet-PA workstation.

6.3 Eluent Generator Control

To control the eluent concentration, add the eluent generator to the PeakNet-PA system (channel) configuration. Then, add commands for controlling the eluent concentration to the PeakNet-PA Analyzer Method.

PeakNet-PA also monitors the EluGen cartridge use and remaining lifetime and displays a warning when the cartridge should be replaced.

Refer to the eluent generator manual and the PeakNet-PA manual and online Help for details.

7 • LM80 Liquids Manager

The LM80 Liquids Manager is located below the AE80 enclosure. The external location permits servicing of eluents, standards, and reagents without opening the enclosure door and subjecting the analytical instrumentation to the environment.

Key features of the LM80 include:

- A control panel, described in the following section.
- One removable polypropylene holder, which accommodates two plastic bottles for standards or reagents. Two-liter bottles (P/N 044129) are shipped with the channel; 1-liter bottles (P/N 044128) can be ordered.




Figure 7-1. LM80 Liquids Manager

DX-800 Process Analyzer

- One 20-liter aluminum NOWPAK II container (P/N 052882) with Teflon liner for eluent, mobile phase, or solvent. A spring-loaded check valve prevents overpressurization of the container. Additional NOWPAKs may be purchased.

NOTE

The Teflon liner (P/N 052885) can be reused if it is being refilled with the same solution. However, if contamination is suspected or if the NOWPAK cannot be pressurized, replace the liner.

The NOWPAK sits directly below the AE80 enclosure. To contain spills, place the container in the provided drip tray. The tray material varies depending on the liquids being used (for example, a polypropylene tray is used for ion chromatography eluents).

- All liquid connections are made with 1/8-in OD PFA Teflon tubing and 10-32 ferrule fittings.
- The gas connection uses a 1/4-in pressfit fitting.

7.1 LM80 Control Panel

- The pressure gauge indicates the pressure applied to the reservoirs and NOWPAK eluent container(s). The recommended operating pressure is 70 to 80 kPa (10 to 12 psi).
- Pressurizing gas is directed to the pressure regulator. A pressure relief valve behind the control panel will open at 100 kPa (15 psi). If the valve opens during operation, turn off the pressure momentarily to allow the valve to reset itself.



The pressure relief valve prevents overpressurization of the LM80 containers, which might damage the containers and injure the user. Never operate the LM80 without the relief valve.



La soupape de détente empêche la surpression des conteneurs du LM80, surpression qui pourrait endommager les conteneurs et blesser l'utilisateur. N'utilisez jamais le LM80 sans la soupape de détente.



Das Überdruckventil verhindert einen Überdruck in den Behältern des LM80. Überdruck kann die Behälter beschädigen und zu Verletzungen des Anwenders führen. Betreiben Sie den LM80 daher niemals ohne Überdruckventil.

- The **A** and **B** controls provide on/off control of gas to the reservoirs and NOWPAK(s). To apply pressure, turn the knob to the **GAS** position. To turn off the gas, turn the knob to the **VENT** position.

NOTE

A second (identical) control panel is installed for applications, such as transition metals, that require two gas supplies (one for reagents and one for standards and eluents). This prevents cross-contamination between reagents such as nitric acid and ammonium acetate buffered PAR.

7.2 Pneumatic Requirements

The reagent and standard reservoirs and NOWPAK eluent containers require a pressurized supply of nitrogen or helium, regulated to between 70 and 80 kPa (10 to 12 psi). The gas purity should be appropriate for the application.

After pressurizing the reservoirs and eluent containers, check the LM80 pressure gauge and the supply tanks (if used) after 15 to 30 minutes. If the pressure is not between 70 and 80 kPa (10 and 12 psi), reset it. It may take several hours for the pressure to stabilize, depending on how much eluent the NOWPAK contains.

To maintain the desired pressure, the reservoirs should be installed within 3 meters (10 feet) of the AE80 enclosure and no more than 0.5 to 1 meter (2 to 3 feet) below the bottom of the enclosure.

DX-800 Process Analyzer

8 • SS80 Sample Selector

8.1 Overview

Multiple sample selection is achieved with the SS80 Sample Selector (P/N 050332). The SS80 is a stand-alone module, so the bulk of the sample flow is isolated from the analytical instrumentation. The basic SS80 contains a multiport sample valve that selects one of seven sample sources. With the addition of a second and third valve, the SS80 can select one of 14 and 21 samples, respectively. Contact Dionex for information about installing additional valves.




Figure 8-1. SS80 with Two Sample Valves

DX-800 Process Analyzer

Sample lines enter the SS80 through the seal bushing. One sample at a time is selected for analysis. One port is selected for diversion to the analysis channel(s). If more than one sampling valve is configured, sample streams are cascaded to the first valve before sample exits the SS80 through the common line to the channel(s).

Samples that are not selected flow continuously out of the valve and exit the SS80. Continuously flowing samples can either be directed to waste or returned to the process. The continuous flow ensures that all samples are fresh when selected and that a representative sample is delivered to the analyzer.

NOTE

If continuous flow is not required, the exit ports on the sample select valve(s) can be plugged with 1/4-28 fittings, provided that the incoming sample pressure does not exceed 0.34 MPa (50 psi).

8.2 Main Components

The SS80 is a wall-mounted unit. To accommodate the control cable length and to minimize flush times between samples, locate the SS80 within 8 meters (25 feet) of the Process Analyzer. The SS80 enclosure has the following main components:

- A mechanical compartment that contains the sample select valves and the tubing.
- A leak sensor in the bottom of the compartment reports liquid leaks to the CC80. When a leak is detected, the **Sampling Leak** annunciator on the CC80 front panel begins flashing. Corrective action can be programmed in PeakNet-PA software.
- An electronics compartment at the rear of the module that contains the valve motors and the distribution board.
- The SS80 main power receptacle is on the rear panel.
- The 9-pin shielded wire (P/N 050363) from the SS80 is connected to the **SS80** connector on top of the AE80 enclosure.

9 • Operation and Maintenance

NOTE

Do not operate the DX-800 for any purpose other than that for which it is designed and described in the user documentation. If there is a question regarding appropriate usage, contact Dionex before proceeding.

9.1 Installation Checklist

Complete the following installation steps before initial operation of the DX-800.

9.1.1 Prepare the Site and Facilities

1. Select the location for the analyzer channel(s), the PeakNet-PA workstation, and the SS80 Sample Selector (optional) at the installation site. Refer to *Installation Requirements and Customer Responsibilities* (Document No. 031176) for appropriate locations and distances.
2. Provide the facilities specified in the installation requirements document. Refer to the document for appropriate voltages, currents, pressures, and flow rates.

9.1.2 Connect the Facilities

1. Connect electrical power to each channel, the PeakNet-PA workstation, and the printer. A power cord is provided with each channel. If the wall receptacle is not a NEMA L7-20R receptacle, arrangements must be made with the Dionex representative prior to installation.

DX-800 Process Analyzer

2. Confirm that LC80 regulators are turned down and valves are turned to the vent position. Connect gas sources to each channel. Tubing and fittings are supplied in the DX-800 installation kit. Shutoff valves for each of the gas inlets should be closed.
3. Connect water sources to each channel. Tubing and fittings are supplied in the DX-800 installation kit. Shutoff valves for each of the sample inlets should be closed.
4. Direct the waste lines into the drain. Refer to the installation requirements document for location and flow rates.

9.1.3 Connect the Communications Cables

1. Attach and dress the DX-LAN cable to the PeakNet-PA workstation and each of the analyzer channels.
2. Connect the SS80 control and power cable from the SS80 to the channel configured with the second power supply in the CC80.
3. Connect the TTL in and/or Relay out cable(s) to each of the channels (if applicable).

9.1.4 Connect the Sample Inlet to the Analyzer Channels

1. Locate the Sample Inlet line for each channel. If multiple channels are configured with the analyzer, route a 3-mm (1/8-in) OD line from the channel farthest from the sample source and install three-way manifolds from this line at each of the remaining channels.
2. Connect the sample line from the analyzer channel(s) to the SS80 Sample Selector or sample source. Refer to the installation requirements document for correct pressures and flow rates.

3. If an SS80 is installed, connect it to the sample panel. Refer to the installation requirements document for correct pressures and flow rates.

9.2 Initial Startup

9.2.1 Turn on the Power

1. Confirm that the circuit breaker on the left side of the CC80 is in the closed/on (up) position.
2. Press the **Power Reset** button on the front door of the channel enclosure(s).
3. Confirm that the power switches for the CC80, analytical pump, analytical detector, and (if applicable) EG40-PA or post-column reagent pump are on.

9.2.2 Start PeakNet-PA

1. Turn on the PeakNet-PA workstation and start PeakNet-PA. Refer to the PeakNet-PA manual for details.
2. Run the Configuration Editor to confirm that all modules are recognized; define the channels (systems), and save the configuration. This will be required for the pump calibration (see Section 9.2.5).

9.2.3 Flush the Flow Path

1. Open the gas inlet valves. Adjust the pressures of the regulators on the rear of the SP80 panel(s) to the following:
 - 275 kPa (40 psi) for the precision displacement pump control
 - 175 kPa (25 psi) for the dilution vessel evacuation
2. Set the LM80 controls to vent. Adjust the pressure of the LM80 to 70 kPa (10 psi).

DX-800 Process Analyzer

3. Prepare eluents, reagents (if used), and standards. Fill containers and pressurize to the appropriate pressures. Refer to the NOWPAK documentation provided with the installation kit for details about filling the eluent containers. Open water supplies to each channel.
4. Manually actuate the SP80 valves to flush the lines with water, diluent, sample, and standards.
5. Fill the dilution vessel with water or diluent and then drain the liquid. Depending on the application, it may be necessary to repeat this process several times or soak the vessel overnight in order to remove trace contaminants.
6. Before installing consumable components, set the analytical pump flow rate to 0.25 mL/min (for a microbore pump) or 1 mL/min (for a standard bore pump), prime and start the pump to flush the chromatography flow path. If a suppressor is being used with external water for regenerant, flush these lines, also.
7. Set the analytical pump to flow in “flow mode.” While pumping through the chromatography flow path (without consumables) confirm that the total backpressure remains below 690 kPa (100 psi).

9.2.4 Install Consumable Components

Install the column(s), suppressor (if used), and any other consumable components according to the product manuals.

9.2.5 Channel Calibration

Precision Displacement Pump Calibration

The Calibration Wizard in PeakNet-PA (accessible from the Configuration Editor) details a step-by-step procedure for calibrating the dilution pump and the loading pump (if present) on the SP80 panel. The purpose of the calibration procedure is to accurately determine the pump stroke volume,

using the same hardware components (tubing, fittings, columns, etc.) used for routine analysis. After replacing any component in the pump flow path that alters the backpressure in the system (concentrator columns, sample loops, or tubing), the pump(s) must be recalibrated.

Standard Loop Calibration

Before running the initial analysis, calibrate the fixed-volume standard loop on the metering (ME) valve. Calibration of the standard loop (and the dilution pump) will determine the correct calibration standard concentration and ensure accurate analytical results.

1. Standard loops are usually made from 0.5-mm (0.020-in) ID PEEK tubing with Dionex 10-32 ferrule fittings. For analyses in the ppm to ppb concentration range, use a loop with a volume of 20 to 100 μL . For analyses of 10 to 100 ppb, the loop volume should be 100 to 250 μL .

The table below indicates the tubing length for several standard loop sizes. *These values are approximations* because tubing IDs vary. After checking the table, cut a piece of tubing to the required length. Be very careful to cut the end square to the axis of the tubing, with no angle. Poorly cut tubing will cause fittings to leak.

Loop Size (μL)	Length of Tubing (cm)
10	4.93
25	12.33
50	24.67
100	49.34
150	74.01
200	98.68
250	123.35
500	246.70
1000	493.40

DX-800 Process Analyzer

2. Install a 10-32 PEEK bolt (P/N 043275) and PEEK ferrule fitting (P/N 043276) on both ends of the tubing cut in Step 1.
3. Install a plug (P/N 042772) in one end of a black coupler (P/N 042627), and then install the coupler on one end of the standard loop. Repeat on the other end of the loop.
4. Weigh the loop on an analytical balance, to the nearest 0.001 gram. Record the weight. Remove the loop from the balance and remove the plugs from the couplers.
5. Using a syringe (P/N 016640) and luer adapter (P/N 024305), fill the standard loop with deionized water. Do not introduce any air into the loop.
6. Reinstall the plug on the end of the coupler from which the water exited. Remove the syringe and luer adapter from the other end of the loop and install the plug in it.
7. Examine the outside of the loop for water droplets. **Carefully** dry any water and then weigh the loop to the nearest 0.0001 gram, if possible.
8. Subtract the weight of the empty loop (Step 4) from the weight of the filled loop (Step 7); the difference is the weight of the water in the standard loop. Repeat Steps 7-8 until four to five consecutive weighings ± 0.009 are achieved.

9 • Operation and Maintenance

9. Multiply the weight of the water by 1000 to obtain the standard loop volume in microliters (μL). The table below lists examples of dilution factors and final concentrations of the diluted standard.

Standard Loop Size (μL)	Dilution Volume (mL)	Dilution Factor	Calibration Standard Concentration (mg/L)	Diluted Standard Concentration ($\mu\text{g/L}$)
10	50	5000	10	2
20	50	2500	10	4
25	50	2000	10	5
50	50	1000	10	10
100	50	500	10	20
150	50	333	10	30
200	50	250	10	40
250	50	200	10	50
250	25	100	10	100

Calculate the diluted standard concentration as follows:

$$V_1 C_1 = V_2 C_2$$

$$C_2 = \frac{V_1 C_1}{V_2}$$

$$D_f = \frac{V_1}{V_2} \frac{25 \mu\text{L}}{50 \text{ mL}} = 2000$$

where:

Standard loop = V_1

Dilution volume = V_2

Dilution factor = D_f

Calibration standard concentration = C_1

Diluted standard concentration = C_2

10. Install the loop between ports 1 and 4 of the metering (ME) valve.

DX-800 Process Analyzer

Analytical Pump Calibration

Calibrate the analytical pump as instructed in the pump manual.

Detector Calibration

Calibrate the detector as instructed in the detector manual.

9.2.6 Initial PeakNet-PA Setup

From PeakNet-PA:

- Run Analyze and configure the analyzer and channels.
- Write and save Methods.
- Download and check Methods.
- Run a standard and define components.
- Calibrate the Method.
- Define and run a Schedule.

See the PeakNet-PA user's guide and online Help for details.

9.3 Routine Startup and Operation

Routine operation consists of first confirming that all hardware is operating properly and then running PeakNet-PA Methods and Schedules to control sample analysis. Refer to the PeakNet-PA manual and the online Help for complete operating instructions.

9.3.1 Routine Startup

1. Confirm that power is on for all modules and accessories. Press the **Power Reset** button on the front door, if necessary.
2. Confirm that the PeakNet-PA workstation is on and that the PeakNet-PA MainMenu is displayed.
3. Confirm that all water and gas utilities are on and adjusted to their proper pressures.
4. Confirm that all water, eluents, standards, and reagents are supplied.
5. Start the Analyze program, load an appropriate Schedule or Method, and verify the following:
 - The eluent flow rates are correct.
 - The detector cells are on and the suppressor is powered (if used).



Always turn on the flow to the suppressor (from the analytical pump) before turning on the detector. Operating the suppressor with no flow going to it will damage the suppressor.

- The post-column flow rate is correct. Adjust the flow rate, if required.
- The dilution vessel is empty. If the vessel contains liquid, press the CC80 **Drain** button.

DX-800 Process Analyzer

- Let each channel stabilize for 20 to 30 minutes. Verify that the detector background has stabilized before beginning the analysis.

9.3.2 Routine Operation

Once each of the channels has stabilized and the Schedule has been started, use the following checklist to monitor operation.

1. Check PeakNet-PA and the CC80 front panel(s) for alarms. If an alarm LED is lighted, eliminate the cause of the problem.
2. Check for any liquid or gas leaks. Isolate and eliminate any leaks.
3. Check levels and flows for all eluents, standards, and reagents. Replenish them as needed.
4. Use a Check Standard in the Schedule to monitor and trend system performance.

9.4 Short-Term Shutdown

Follow these procedures to prepare the DX-800 for a shutdown of a few days to less than two weeks. For long-term shutdowns, see Section 9.5.

1. On the Analyzer Status screen, press the Standby button for the analyzer being shut down.
2. Push the **Emergency Off** button on the front door of each channel.
3. Turn the LM80 gas controls to Vent and turn off the LM80.
4. Shut off the water and gas facilities to each of the channels.
5. Leave the columns installed and filled with eluent.
6. Empty standards bottles if stability is questionable.
7. Exit PeakNet-PA and shut down the PeakNet-PA workstation (optional).

9.5 Long-Term Shutdown

Follow these procedures to prepare the DX-800 for a shutdown of two weeks or more. For short-term shutdowns, see Section 9.4.

1. On the Analyzer Status screen, press the Standby button for the analyzer being shut down.
2. Empty and rinse eluent, standard, and reagent bottles.
3. Prepare the columns, suppressors, and EluGen Cartridges for long-term storage as instructed in the product manuals.
4. Flush the pumps, valves, post-column system, and interconnecting tubing with deionized water. Blow out the lines with high purity nitrogen or helium.
5. Turn the LM80 gas controls to Vent and turn off the LM80.
6. Shut off the water and gas facilities to each of the channels.
7. Push the **Emergency Off** button on the front door of each channel.
8. Exit PeakNet-PA and shut down the PeakNet-PA workstation.

9.6 Maintenance

This section describes routine maintenance procedures that users can perform. Any other maintenance procedures must be performed by Dionex personnel.

Establish a routine maintenance program based on the guidelines here, as well as information in user manuals for other elements of the system: the analytical pump, detector, columns, etc. Following a strict maintenance schedule ensures proper operation of the DX-800.

NOTE

Dionex recommends recording the date on which each routine maintenance procedure is performed. Besides ensuring that the procedures are accomplished, a maintenance log is very helpful when troubleshooting the system.

9.6.1 Daily Maintenance

Completion Time: 10-15 min for a dual-channel analyzer

Component or Feature	Action
Gas pressures	Check house pressure; the cylinder must have enough pressure to supply gas for the day. Air/N ₂ regulator for SP80 pumps=275 kPa (40 psi) LM80 regulator=70 kPa (10 psi) Helium pressure for dilution vessel=175 kPa (25 psi)
Reagent supplies	Check all liquid levels; replenish if necessary. Eluent for the day=1 L minimum (2 L recommended) Stock standard solution for the day=1 L minimum Regenerant water pressure=100 kPa (15 psi) Deionized water pressure to enclosure=100 to 140 kPa (15 to 20 psi) Sample line pressures=70 to 140 kPa (10 to 20 psi minimum); 100 kPa (15 psi recommended)
Sample lines	Make sure all waste lines flow freely.

DX-800 Process Analyzer

Air and liquid lines	Check for leaks or spills. Isolate and repair leaks; clean up spills. Rinse dried chemicals from components with DI water.
All pumps	Check for crimping; replace damaged lines. Check for piston seal leaks; replace defective seals.
Conductivity detector	Record the total conductivity readings at the beginning of a run.
SRS, DS3, EGC	Check for leaks. If Autosuppression is being used, check that bubbles are flowing from the SRS regenerant outlet line.
Chromatography	Check Trend plots and chromatograms for trending problems (missed peaks, etc.).
Printer	Make sure there is paper. Check the ink or toner cartridge.

9.6.2 Weekly Maintenance

Completion Time: 30-40 min for a dual-channel analyzer

Component or Feature	Action
Standard solutions	Prepare new solution for check standard and calibration standard.
Analytical pump	Rinse pistons. Record pump pressure when load/inject valve (LI) is in the load position.
SS80 valve	Check for leaks.
Gas and drain connections	Check all connections, including the drain manifold and the fluid connection panels. Check for accumulated liquid on the inside bottom cover of the enclosure and underneath the enclosure. Fix leaks promptly.
Power and signal connections	Visually inspect all connections and cables. Secure loose connections; move pinched or strained cables.

9.6.3 Biweekly Maintenance

Completion Time: 1-2 hrs for a dual-channel analyzer

Component or Feature	Action
Reagent reservoirs	Thoroughly rinse all reagent reservoirs with deionized water to remove precipitates.
Eluents, reagents	Prepare new eluents and reagents.
Eluent trap columns (if used)	Replace columns (may be required weekly).

9.6.4 Monthly Maintenance

Completion Time: 1-2 hrs for a dual-channel analyzer

Component or Feature	Action
In-line filters	Replace all in-line filters.
Guard columns	If eluent pressure increases by 1.4 MPa (200 psi), replace the bed support in the guard column inlet. If the pressure does not return to near the original for this column, replace the guard column.
Air filter	Clean with warm water whenever a fine layer of dust or lint is visible. Establish a cleaning schedule, taking local air quality into account.
Enclosure	Clean with a mild soap solution and then rinse with water. Chips and scratches are invitations to corrosion.

9.6.5 Quarterly Maintenance

Completion Time: 1-2 hrs for a dual-channel analyzer

Component or Feature	Action
Pump seals	Replace pump seals.
Rotary valves	Replace rotors and stators in the load/inject (LI) and metering (ME) valves.
Analytical pump	Calibrate flow and pressure.
EluGen Cartridge	Replace if necessary. Refer to the lifetime remaining value in PeakNet-PA.

DX-800 Process Analyzer

10 • Troubleshooting

This chapter is a guide to troubleshooting problems that may occur while operating the DX-800. Turn to the section of this chapter that best describes the operating problem; there, the possible causes of the problem are listed in order of probability.

If you are unable to eliminate a problem, contact Dionex for help. In the U.S., call Dionex Technical Support at 1-800-346-6390. Outside the U.S., call the nearest Dionex office.

NOTE

For troubleshooting help specific to other aspects of the system, such as the analytical pump or chromatography consumable components, refer to the relevant user manual.

10.1 Troubleshooting Strategies

For any system, the initial troubleshooting objective is to isolate the source of the problem to a specific module or component, or to an aspect of the analysis. Once this has been done, corrective action can be taken. This manual cannot address every possible symptom and failure; however, the techniques discussed here can be applied to problem solving throughout the DX-800 system.

A solid understanding of system operation is necessary for troubleshooting. Rather than immediately assuming that a problem exists, first check the user manuals to verify that correct operating procedures are being followed.

Also, it is essential that users keep a log of all maintenance-related activities (when eluents are prepared, when columns are changed, etc.), since this can provide valuable insights. For example, if the chromatogram on an anion system seems to have undergone a radical and sudden change, check the log for the date that the latest batch of eluent was placed in service. If the problem was first observed with the new batch, verify that the eluent was properly

prepared. Maintaining a written record of problems and their resolution can help solve similar problems in the future.

In summary, an effective troubleshooting strategy requires that users:

1. Understand the operation of the entire system.
2. Maintain a maintenance log.
3. Isolate the problem to either the hardware or chemistry.
4. Refer to the troubleshooting and service sections of the appropriate user manual.

10.2 Liquid Leaks



If leaking liquid creates a hazard, stop the leak immediately by turning off the flow at the source.



Si une fuite de liquide crée un danger, arrêtez immédiatement la fuite en fermant l'écoulement à la source.



Wenn eine Gefährdung durch austretende Flüssigkeit besteht, stoppen Sie die Leckage unmittelbar, indem Sie den Fluß an der Quelle abstellen.

NOTE

When cutting tubing and preparing fittings, avoid crimping the tubing. Crimped tubing is a common cause of high backpressure.

- **Leaking fitting**

Make sure that all liquid line connections are tight. If a fitting continues to leak, replace it. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

- **Broken liquid line**

Cut the tubing at the break and install a new fitting. (If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.) Or, if cutting the tubing would make it too short, replace it. The new tubing must be the same type, and have the same internal diameter, as the tubing it replaces.

- **Blocked or improperly installed waste line**

1. Make sure the waste lines from the valves, detector cell, and SRS (if installed) are not crimped or otherwise blocked.
2. Make sure the waste lines from the channel are not elevated at any point after they exit the AE80 enclosure. Waste lines should be clear and open to the atmosphere.
3. If 10-32 ferrule fittings are installed on Teflon (PTFE or PFA) tubing, make sure the fittings are not overtightened, thus pinching off the tubing. If in doubt, cut off the end of the tubing and reconnect the fitting. See *Installation of Dionex Ferrule Fittings* for instructions.

10.2.1 Precision Displacement Pump Leaks

The loading pump and dilution pump are precision displacement pumps. Depending on the DX-800 configuration, one or both of these pumps is installed on the SP80 Sample Preparation panel.

- **Loose check valve fittings**

With the pump flowing, use a wrench to tighten the inlet and outlet check valve fittings just until the leak stops.

- **Defective piston seal**

A defective piston seal allows leaks between the pump head and end plate and/or from the piston rinse ports on the pump head. Check the pump head for leaks and replace the seal if necessary.

10.2.2 Dilution Vessel Leaks

- **Loose connections**

Make sure all fittings are tightened securely.

- **Dilution vessel contains excess liquid**

1. The dilution vessel may contain partially diluted sample or standard from a previous analysis. To purge the vessel, press the **Local/Remote** button on the CC80 front panel to select Local mode, and then press the **Drain** button. The drain operation is programmed to take 3 minutes.

If the vessel empties in less than 3 minutes, you may manually terminate the drain operation. Before doing so, check the dilution vessel waste line: If there is no liquid in the line, the vessel is empty. To terminate the drain operation, press the **Drain** button again, or press the **Gas Valve** button (this vents the dilution vessel).

2. The pressure applied to the dilution vessel may be insufficient to empty it. Make sure the high purity gas supply is regulated to 170 to 240 kPa (25 to 35 psi).
3. Make sure that no more than 250 mL of liquid is pumped into the dilution vessel. In the normally open (default) position, the dilution vessel (DV) valve purges the dilution vessel to waste. Check the CC80 Method to verify that the valve remains on (open) while the dilution pump is running.

10.2.3 SS80 Leaks

When a sample line leaks, the **Sampling Leak** LED on the CC80 front panel begins flashing.



If leaking liquid creates a hazard, turn off the **SS80** power switch immediately, and then stop the leak by turning off the flow at the source.



Si une fuite de liquide crée un danger, éteignez immédiatement le courant du **SS80** avec son interrupteur puis arrêtez la fuite en fermant l'écoulement à la source.



Schalten Sie den **SS80** sofort am Netzschalter aus, wenn eine Gefährdung durch austretende Flüssigkeit besteht. Stoppen Sie die Leckage unmittelbar, indem Sie den Fluß an der Quelle abstellen.

- **Leaking fitting(s)**

Tighten any leaking fittings. If a fitting continues to leak, replace it. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

- **Damaged tubing**

If a piece of tubing is crimped or otherwise damaged, replace it. If this is inconvenient, patch the tubing by cutting out the bad section and inserting a new piece of tubing with a coupler (P/N 040240) on each end.

NOTE

Routinely patching tubing increases the possibility of leaks.

10.2.4 EluGen Cartridge Leaks

The cartridge must be replaced. Refer to the eluent generator manual for instructions.

10.2.5 Degas Assembly Leaks

The degas assembly must be replaced. Refer to the eluent generator manual for instructions.

10.3 Air and Gas Leaks

Air leaks, which can cause excessive air consumption, are usually audible. Gas leaks cause sluggish liquid delivery, unreliable pump operation, and excessive gas consumption. Minor gas leaks can sometimes be felt, while major gas leaks are usually audible.

To detect a minor gas leak, shut off the gas at the source and then check the pressure gauge for a drop in pressure. Repeat as often as necessary until the leak is found.

NOTE

The use of Snoop or other dilute soap solutions for leak detection will contaminate the tubing. Use water if desired.

- **Leaking fitting**

If the fitting is stripped, cross-threaded, or otherwise damaged, replace it. If the fitting is not damaged, securely tighten it. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

If the leak continues, cut the tube off the fitting and replace the fitting.

- **Leaking standard or reagent reservoir**

Sluggish liquid delivery is usually due to a helium or nitrogen leak from a reservoir. Follow the strategy above to eliminate leaks at fittings and caps.

- **Damaged tubing**

Over time, gas tubing can become compromised by chemical fumes (for example, eluent vapor in the air lines to the NOWPAK containers). If this happens, replace the tubing.

- **Defective air piston O-ring/seal on precision displacement pump**

Although only one O-ring/seal set may be defective, both sets must be replaced at the same time. Contact Dionex for assistance.

10.4 Excessive System Backpressure

After referring to the column documentation for the recommended operating pressure, check the analytical pump **MAIN** screen for the current system backpressure. A system backpressure of 3.4 MPa (500 psi) above the recommended value for the column is excessive.

Excessive backpressure can cause leaks, or even irreparable damage to system components. Crimped tubing is a common cause of high backpressure.

- **Excessive flow rate through columns**

Make sure the pump is set to the flow rate specified by the analytical protocol. If the pump needs to be recalibrated, go to the **FLOW CALIBRATION** screen (from the **DIAGNOSTIC MENU**, select **CALIBRATION MENU**, and then select the **FLOW CALIBRATION** option) and follow the on-screen instructions.

- **Restriction in chromatography flow path**

1. Follow these steps to isolate the source of the high backpressure:
 - a. Disconnect the analytical pump eluent line from the load/inject (LI) valve. Turn on the pump and monitor the operating pressure; it should not exceed 0.3 MPa (50 psi) for a standard bore system or 1.4 MPa (200 psi) for a microbore system.

- b. Begin reinstalling system components, starting with the load/inject valve, while monitoring the system pressure. When the analytical column is connected, the pressure should increase to the recommended value. Components other than the analytical column should add less than 2.1 MPa (300 psi) to the system backpressure.
- If the load/inject valve is the cause of the high backpressure, refer to the manual included with the Rheodyne rebuild kit for corrective action.
- If the concentrator or guard column is the cause, replace the column.
- If the analytical column is the cause, refer to the column manual for corrective action.
- If the SRS is the cause, refer to the SRS manual for corrective action.
- If the EluGen cartridge is the cause, refer to the eluent generator manual for corrective action.

10.5 Channel Stops Running

- **Dilution vessel contains residual liquid**

After termination of a Method that involves dilution, some partially diluted sample or standard may remain in the dilution vessel. To purge the vessel, press the **Local/Remote** button on the CC80 front panel to select Local mode, and then press the **Drain** button. The drain operation is programmed to take 3 minutes.

If the vessel empties in less than 3 minutes, you may manually terminate the drain operation. Before doing so, check the dilution vessel waste line: If there is no liquid in the line, the vessel is empty. To terminate the drain operation, press the **Drain** button again, or press the **Gas Valve** button (this vents the dilution vessel).

10.6 Module(s) Does Not Power Up



Electrical system circuits carry dangerous voltages. Disconnect all power before working on them.



Les circuits du système électriques ont des tensions dangereuses. Débranchez toute l'alimentation électrique avant de travailler sur les circuits.



Elektrische Schaltkreise führen gefährliche Spannungen. Entfernen Sie alle Stromversorgungen, ehe Sie daran arbeiten.

- **Internal power cord(s) not connected**

Make sure the power cords are connected from the analytical pump and the detector to the appropriate AC outlets on the CC80 rear panel.

- **PC80 Post-Column Pump power cord not connected**

Make sure the PC80 power cord is connected to the appropriate AC outlet on the CC80 rear panel.

- **System power cord not connected**

1. Make sure the **Power Reset** lamp on the front door of the AE80 enclosure is illuminated.
2. Make sure the modular power cord is connected from the **POWER IN** connector on top of the AE80 enclosure to the main power.

- **Eluent generator power turned off**

Verify that the power switch on the EG40-PA electronics module is on. When the power switch is on, the power LED is illuminated. To access the EG40-PA electronics module, loosen the thumbscrew on the right side of the SP80 Sample Preparation panel and swing the panel out. The electronics module is at the

DX-800 Process Analyzer

lower rear of the AE80 enclosure. For additional troubleshooting information, refer to the eluent generator manual.

- **Main power turned off**

After pressing the **Emergency Off** button on the front door of the AE80 enclosure, remember to press the **Power Reset** button to restore power.

- **CC80 circuit breaker tripped**

1. Open the front door of the AE80 enclosure. Loosen the two retaining screws in the module chassis. Pull out the chassis just until the first stop on the sliders is engaged.
2. To reset the circuit breaker, flip up the switch on the left side of the CC80.
3. Push the module chassis back into the enclosure and tighten the retaining screws.
4. Close the enclosure door.



If the breaker continues to trip, the circuit may be shorted or overloaded. Disconnect all power and contact Dionex for assistance.



Si le disjoncteur continue de sauter, le circuit peut être court-circuité ou surchargé. Débranchez toute l'alimentation électrique et contactez Dionex pour obtenir de l'aide.



Wenn der Unterbrecherschalter kontinuierlich auslöst, kann es sein, daß die Schaltung kurzgeschlossen oder überlastet ist. Entfernen Sie alle Stromversorgungen und wenden Sie sich an Dionex.

- **Blown fuse**

Five IEC127 3.15 amp fast-blow fuses (P/N 954745) are installed inside the CC80. Replace the open fuse (see Section 11.3).



If the fuse continues to blow, the circuit may be shorted or overloaded. Disconnect all power and contact Dionex for assistance.



Si le fusible continue de sauter, le circuit peut être court-circuité ou surchargé. Débranchez toute l'alimentation électrique et contactez Dionex pour obtenir de l'aide.



Wenn die Sicherung weiterhin durchbrennt, kann es sein, daß die Schaltung kurzgeschlossen oder überlastet ist. Ziehen Sie den Netzstecker und wenden Sie sich an Dionex.

10.7 CC80 Sample LED Displays Spinning Segments

- **Sample select valve is between positions**

There is no problem; the rotation of the LED display segments indicates that the selected valve is switching positions.

10.8 CC80 Sample LED Displays *EE*

- **Malfunctioning sample select valve**

The sample select valve is frozen in position or the valve sensor is not working properly. Press the **Select** button to select a different valve, and then reselect the first valve. If *EE* is displayed again, the valve is broken and must be replaced. Contact Dionex for assistance.

10.9 CC80 Analyzer Leak LED Is Flashing

The **Analyzer Leak** LED flashes when a leak is detected on the LC80 Liquid Chromatography panel or SP80 Sample Preparation panel.

- **Leaking fitting**

Find and eliminate the source of the leak (see Section 10.2).

- **Loose electrical connection**

1. Verify that the leak sensor cable is plugged into the leak detector connector (J7) on the LC80 distribution board. The board is located behind the LC80 panel. To access this area, loosen the two screws securing the LC80 to the module chassis and swing the panel open.
2. If the **Analyzer Leak** LED immediately begins flashing again, the leak sensor is out of calibration.

To recalibrate the sensor:

1. Open the door of the AE80 enclosure. Press the CC80 **Power** switch to turn off the power. After a few seconds, turn on the power again.
2. Wait until all CC80 front panel LEDs are illuminated, and then press the **Display Refresh** button on the enclosure door to calibrate the leak sensor.
3. Close the enclosure door.

10.10 Precision Displacement Pump Does Not Prime

The loading pump and the dilution pump are precision displacement pumps. Depending on the DX-800 configuration, one or both of these pumps is installed on the SP80 Sample Preparation panel.

- **Pump not plugged in**

One at a time, check the cable connections from the pump to the distribution board on the rear of the SP80 panel. The 5-pin connector plugged into the pump controller board (under the white splash cover) should be connected to the following

location on the distribution board: LP Sensor (for the loading pump) or DP Sensor (for the dilution pump).

- **Pump contains trapped air**

When air becomes trapped in the pistons, the stroke speed increases significantly as the pump attempts to push out the trapped air.

Before priming the pump, make sure the valve immediately before it is open by pressing the appropriate button on the CC80 front panel: the **SS Valve** button (for the loading pump) or the **DI Valve** button (for the dilution pump). In the open position, the gentle flow of water through the valve to the pump is effective at removing trapped air.

If the stroke speed remains unusually fast, even though the valve is open, turn off the pump for 2 to 3 minutes (leaving the valve open). Turn on the pump again. The pump speed should return to normal within a few seconds, confirming that the pump is primed.

10.11 No Sample Delivered to Loading Pump

The CC80 monitors the number of pump strokes per second. If the rate exceeds three strokes per second after 20 seconds of operation, it indicates that no liquid is being pumped. If this occurs, the **Loading Pump** LED on the CC80 front panel will begin flashing continuously.

Press **Alarm Reset** to turn off the LED and fix the problem. If the LED starts flashing again, the problem was not eliminated. Contact Dionex for assistance.

- **No sample supplied to channel**

1. Make sure the dilution vessel contains sample or standard.
2. Make sure the air supply for the dilution vessel is regulated to between 34 and 69 kPa (5 to 10 psi).

- **No liquid passing through SS80 valve**

Positions 1 and 2 of DIP switch #1 record the number of sample select valves installed in the SS80. This information is reported to the CC80 Moduleware and PeakNet-PA software. Verify that the switch settings are correct (see Section 3.2.3).

- **Liquid leaks**

Check for leaks in the liquid lines or fittings in the sample flow path. Tighten or replace leaking fittings.

- **Malfunctioning valve**

Press the **Local/Remote** button on the CC80 front panel to select Local mode. One at a time, press the flow chart button that represents each valve in the flow path and listen for the click that occurs when the valve is actuated. Any valve that is not actuated should be replaced. Contact Dionex for assistance.

- **Damaged pump controller board**

The loading pump must be replaced. Contact Dionex for assistance.

10.12 Loading Pump Delivers Inconsistent Volume

- **Excessive fluctuation in sample inlet pressure**

Regulate the sample inlet pressure to between 70 and 280 kPa (10 to 40 psi).

- **Defective pump seal**

A defective seal allows liquid leaks between the pump head and the end plate to which it is attached. Check the pump head for leaks and replace the seal, if necessary.

- **Gas leaks from dilution vessel**

Make sure all fittings are tightened securely.

- **Load/inject (LI) valve leaks**

Tighten all fittings. If a fitting continues to leak, replace it.

10.13 Irregularity in Loading Pump

- **Dilution vessel gas leak**

An irregularity in the pump stroke speed or rhythm when pumping from the dilution vessel may indicate a gas leak.

Follow this procedure to check for a gas leak from the dilution vessel:

1. Press the **Local/Remote** button on the CC80 front panel to select Local mode.
2. Press the **DV Valve** button to seal off the bottom of the dilution vessel.
3. Press the **SS Valve** button.
4. Press the **Gas Valve** button to pressurize the dilution vessel.
5. Press the **ME Valve** button.
6. Turn on the gas supply for the dilution vessel.
7. Using deionized water, verify that there are no gas leaks from the fitting on top of the dilution vessel.

NOTE

Do not use Snoop or other dilute soap solutions for leak detection, as they will contaminate the tubing. Use water, if desired.

10.14 Dilution Pump Does Not Pump

The CC80 monitors the number of pump strokes per second. If the rate exceeds three strokes per second after 20 seconds of operation, it indicates that no liquid is being pumped. If this occurs, the **Dilution Pump** LED on the CC80 front panel will begin flashing.

Press **Alarm Reset** to turn off the LED and fix the problem. If the LED starts flashing again, the problem was not eliminated. Contact Dionex for assistance.

- **No diluent supply**

1. Make sure there is a supply of diluent (usually deionized water) to the channel.
2. Check for leaking fittings on the diluent lines. Tighten or replace fittings as necessary.

- **Inadequate air supply (stroke speed slower than normal)**

The dilution pump is air-driven. For optimal pump performance, the facility gas supply (air or nitrogen) must be regulated to a minimum of 280 kPa (40 psi).

- **Diluent (DI) valve is off (or closed)**

1. Check the CC80 Method to verify that the diluent (DI) valve is on (open) while the dilution pump is running.
2. Verify correct valve operation as follows:
 - a. Press the **Local/Remote** button on the CC80 front panel to select Local mode.
 - b. Press the **DI Valve** button on the CC80 flow chart to open the diluent valve.
 - c. Press the **Dilution Pump** button to turn on the pump. Check the waste line to verify that diluent is being directed to waste.

- **Blockage in air or liquid connections**

Make sure the lines to and from the dilution pump are not crimped or otherwise blocked.

If 10-32 ferrule fittings are installed on Teflon (PTFE or PFA) tubing, make sure the fittings are not overtightened, thus pinching off the tubing. If in doubt, cut off the end of the tubing and reconnect the fitting. See *Installation of Dionex Ferrule Fittings* for instructions.

- **Leaking pump seal**

A defective seal allows liquid leaks between the pump head and the end plate to which it is attached. Check the pump head for leaks and replace the seal if necessary.

- **Damaged pump controller board**

The dilution pump must be replaced. Contact Dionex for assistance.

10.15 Dilution Pump Delivers Inconsistent Volume

- **Variation in diluent supply pressure**

Make sure the diluent supply is regulated to between 170 and 240 kPa (25 to 35 psi).

- **Metering (ME) valve leaks**

Refer to the manual included with the Rheodyne rebuild kit for corrective action.

- **Blockage in air or liquid connections**

Make sure the lines to and from the dilution pump are not crimped or otherwise blocked.

If 10-32 ferrule fittings are installed on Teflon (PTFE or PFA) tubing, make sure the fittings are not overtightened, thus pinching off the tubing. If in doubt, cut off the end of the tubing and reconnect the fitting. See *Installation of Dionex Ferrule Fittings* for instructions.

- **Dirty or worn pump check valve**

Clean the check valves (see Section 11.2).

- **Defective pump seal**

A defective seal allows liquid leaks between the pump head and the end plate to which it is attached. Check the pump head for leaks and replace the seal if necessary.

10.16 Dilution Vessel Does Not Empty

- **Fitting allows gas (helium) leaks**

If the fitting is stripped, cross-threaded, or otherwise damaged, replace it. If the fitting is not damaged, tighten it securely. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

If the leak continues, cut the tube off the fitting and replace the fitting.

- **Dilution vessel inadequately pressurized**

1. Disconnect the fitting on the helium inlet line to the dilution vessel lid and check for flow from the line. If there is no flow, check the pressure for the helium supply.

To check the helium supply pressure:

- a. Loosen the thumbscrew on the right side of the SP80 Sample Preparation panel and swing the panel open.
 - b. The pressure gauge is located on the rear of the SP80 panel, directly behind the dilution vessel. The gauge should read between 170 and 240 kPa (25 and 35 psi); if it does not, contact Dionex for assistance.
2. Check the CC80 Method to verify that the gas valve is in the pressurize (1, or normally open) position. If it is not, edit the Method.
 3. Make sure the waste line is clear and open to the atmosphere.

- **Pressure relief valve compromised**

The pressure relief valve is on the rear of the SP80 panel, behind the dilution vessel lid. If the vessel is insufficiently pressurized, the relief valve may be activated prematurely.

To check the pressure relief valve operation:

1. Loosen the thumbscrew on the right side of the SP80 and swing the panel open. The regulator is in the upper left corner of the panel, below the dilution vessel.
2. Set the pressure regulator to 345 kPa (50 psi).
3. Using a screwdriver, turn the valve adjustment screw to the right just until a hissing sound is audible. When the hissing starts, reset the pressure to 276 kPa (40 psi). If the hissing starts again, the pressure relief valve is broken and must be replaced. Contact Dionex for assistance.

- **Dilution vessel not drained after Method is aborted**

1. When a Method that involves dilution is terminated while the dilution vessel is being filled or emptied, some diluted sample or standard may remain in the vessel.

To ensure complete drainage:

- a. Press the **Local/Remote** button on the CC80 front panel to select Local mode.
- b. Press the **Drain** button on the CC80 flow chart to pressurize the dilution vessel and direct residual liquid to waste. The drain operation is programmed to take 3 minutes.

If the vessel empties in less than 3 minutes, you may manually terminate the drain operation. Before doing so, check the dilution vessel waste line: If there is no liquid in the line, the vessel is empty. To terminate the drain operation, press the **Drain** button again or press the **Gas Valve** button (this vents the dilution vessel).

- Dilution vessel (DV) valve not operating correctly

In the normally open (default) position, the dilution vessel valve purges the dilution vessel to waste. Check the CC80 Method to verify that the valve remains on (open) while the dilution pump is running.

10.17 Column Heater Does Not Heat

Verify that the column heater cable is plugged into the column heater connector on the LC80 distribution board. The board is located behind the LC80 panel. To access this area, loosen the two screws securing the LC80 to the module chassis and swing the panel open.

10.18 Inoperative Sample Select Valve

- Incorrect CC80 DIP switch setting

Positions 1 and 2 of DIP switch #1 record the number of sample select valves installed in the SS80. This information is reported to the CC80 Moduleware and PeakNet-PA software. Verify that the switch settings are correct (see Section 3.2.3).

10.19 Lack of Flow at Selected Sample Outlet

- Inadequate sample pressure and/or flow to SS80

The minimum sample flow to the SS80 is 15 mL/min; the minimum sample inlet pressure is 69 kPa (10 psi).



Various types of chemicals are used in the DX-800, depending on the application that is being performed. Follow all appropriate hazardous materials and safety guidelines for chemicals when operating the DX-800.



Différents types de produits chimiques sont utilisés dans le DX-800, selon l'application à effectuer. Respectez toutes les directives de sécurité sur les matières dangereuses pour les produits chimiques lors de l'utilisation du DX-800.



Je nach Anwendung, die gerade läuft, werden im DX-800 verschiedenartige Chemikalien verwendet. Beachten Sie beim Betrieb des DX-800 alle entsprechenden Sicherheitsrichtlinien bezüglich gefährlicher Stoffe für die verwendeten Chemikalien.

- **Leaking fitting and/or sample select valve**

Make sure all liquid line connections are tight. If this does not stop the leak, replace the fitting. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

- **Blocked tubing**

1. Replace all kinked or pinched tubing.
2. If particulates are blocking the sample lines, replace the affected tubing and filter the incoming sample (5 µm particulate size).

- **Sample select valve does not respond when energized**

Using the CC80 front panel **Select** buttons, select each valve position in turn and check for flow at the SS80 sample outlet. If a valve is unresponsive, it should be replaced. Contact Dionex for assistance.

10.20 No Peaks Detected

- **Detector not ready**

Verify that the cell (or lamp) is on and that there is flow through the cell.

- **Eluent generator not working**

Verify that the power switch on the EG40-PA electronics module is on. When the power switch is on, the power LED is illuminated. To access the EG40-PA electronics module, loosen the thumbscrew on the right side of the SP80 Sample Preparation panel and swing the panel out. The electronics module is at the

DX-800 Process Analyzer

lower rear of the AE80 enclosure. For additional troubleshooting information, refer to the eluent generator manual.

- **No sample injected**

1. Make sure the CC80 Method includes an inject command at the appropriate step.
2. Verify correct operation of the load/inject (LI) valve as follows:
 - a. Press the **Local/Remote** button on the CC80 front panel to select Local mode. Press the **Load/Inject** button on the flow chart a few times to verify that the valve switches between the load and inject positions.
 - b. If the valve does not switch positions, make sure the valve cable is connected to the LIV connector on the LC80 distribution board. The distribution board is located on the rear of the LC80 panel. To access this area, loosen the two screws securing the LC80 to the module chassis and swing the panel open.
 - c. If the valve is correctly connected to the distribution board but fails to switch positions, contact Dionex.
3. If no peaks are detected when analyzing a sample, verify that there is sample flow between the SS80 and the channel.
4. If no peaks are detected when preparing and analyzing a calibration standard, check the following:
 - a. Verify that the reservoir contains stock standard. The gas supply for the reservoir should be regulated to between 70 and 80 kPa (10 to 12 psi).
 - b. Verify that the metering (ME) valve is programmed correctly in the CC80 Method and is operating correctly.
 - c. Verify that the prepared calibration standard is flowing from the dilution vessel to the load/inject (LI) valve. Make sure the vessel is pressurized to between 170 and 240 kPa (25 to 35 psi).

10.21 Spurious Peaks

Spurious peaks in a chromatogram may be late-eluting peaks from a previous injection or may result from a contaminated valve or poor sample loading technique. Spurious peaks sometimes co-elute with peaks of interest, resulting in nonreproducible peak area or height.

To verify the cause, first select a run time that ensures that the peak elutes with the injection. If this is a late-eluting peak, there are two ways to accomplish this:

- Adjust the run time to permit the peak to elute with the injection, *or*
- Increase the eluent (or mobile phase) concentration, causing the peak to elute earlier. (Be sure to equilibrate the system with the composition required for the isocratic run or used to start the gradient separation.)

- Insufficient time between sample injections

Wait until the previous sample has been completely eluted before making another injection.

- Insufficient rinse between samples

Edit the CC80 Method to increase the time allowed for rinsing the sample line from the analyzer to the SS80 (or the sample source for a single-sample analyzer). In general, rinse the lines with at least five tubing volumes. For ultra-trace analyses or certain samples, larger rinse volumes may be required. The sample tubing is 1.5 mm (0.060-in) ID, with a volume of about 0.55 mL/ft.

To test for sample carryover, run the highest calibration standard or most concentrated sample, followed by a blank. (For ultra-low samples, use the highest calibration standard **only**.) The carryover should be less than 1%.

- **Analytical column degraded**

Clean the column as instructed in the column manual. When the column is used with a weak eluent system and samples contain an appreciable level of polyvalent anions or cations, the polyvalent anions or cations may contaminate the column. If this occurs, retention times for analytes will decrease and spurious, normally inefficient peaks may appear at unexpected times.

- **Baseline upset**

1. First, run a gradient without making an injection. Study the baseline; if there are spurious peaks, the analytical column may be contaminated (see above).
2. Run a second gradient. This time, switch the load/inject (LI) valve, but not the injection sample or standard (the sample loop should contain deionized water or eluent).

A baseline upset, especially at the beginning of a chromatogram, is probably caused by actuation of the load/inject valve or column switching valve. Clean the valve.

NOTE

A minor baseline disturbance at the start or end of a chromatogram may be disregarded unless it interferes with the quantitation of peaks of interest.

10.22 Poor Peak Resolution

- **Retention times too short**

1. Verify that the selected flow rate is the one specified by the analytical protocol. When analytes elute too fast, their resolution is compromised.
2. If an eluent generator is installed, verify that the correct concentration and flow rate combination is programmed in the PeakNet-PA Method. Refer to the eluent generator manual for details.

- **Incorrect eluent composition or concentration**

1. Prepare fresh eluent. An excessively strong eluent causes peaks to elute more quickly (and vice versa).
2. The gradient pump proportioning valve may be malfunctioning (applicable when the gradient pump is proportioning the eluent by combining eluents from more than one container).

To check the proportioning valve operation:

Run an isocratic Method, using a container with the correct eluent composition. If retention times and resolution are recovered, the proportioning valve is defective. Replace the valve as instructed in the gradient pump manual.

3. Validate the gradient pump performance. For instructions, refer to the PeakNet-PA manual or online Help. Pump validation is included in the PeakNet-PA qualification programs. If this does not isolate the cause of the problem, contact Dionex for assistance.

- **Column contamination**

Clean the column as instructed in the column manual. Column contamination sometimes results in a loss of column capacity because some exchange sites are no longer available for the sample ions. Polyvalent anions or cations may be concentrating on the column.

Impurities in the chemicals or deionized water used to prepare eluent can contaminate the column. When preparing eluent, be sure to use reagent-grade chemicals or chemicals of the purity recommended in the column manual. Use only ASTM Type II (18.0 M Ω /cm resistance or 1 μ S) deionized water.

The symptoms below apply to loss of resolution for early-eluting peaks only.

- **Incorrect eluent concentration**

Remake the eluent. If gradient elution is being used, verify that the pump Method in PeakNet-PA is correct.

If an eluent generator is installed, verify that the correct concentration and flow rate combination is programmed in the PeakNet-PA Method. Refer to the eluent generator manual for details.

- **Column overloading**

Strongly-retained ions may be rinsing more weakly-retained ions off the concentrator column. Concentrate 5 mL less of sample; if linearity improves, continue decreasing the sample amount until linearity fails to improve. Increase the sample amount in increments of 5 mL, if desired, as long as linearity remains acceptable.

10.23 Small Peaks Detected

- **Air leak from loading pump**

Listen for air leaks from the air cylinder of the loading pump. If there is a leak, the two O-ring/seal sets (P/N 049317) on the air piston must be replaced. Contact Dionex for assistance.

10.24 Peak Height Greater Than Expected

- **Air leak from dilution pump**

Listen for air leaks from the air cylinder of the loading pump. If there is a leak, the two O-ring/seal sets (P/N 049317) on the air piston must be replaced. Contact Dionex for assistance.

10.25 Poor Peak Area (or Height) Precision

Poor peak area or height precision is indicated by the inability to reproduce results from injection to injection.

- **Leaks in sample or standard flow path**

Tighten any leaking fittings. If a fitting continues to leak, replace it. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

- **Leaking load/inject (LI) valve**

Make sure all liquid line connections are tight. If this does not stop the leak, replace the fitting. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions. If the leak persists, contact Dionex for assistance.

- **Insufficient or inconsistent sample or standard flow to loading pump**

For calibration standards prepared in the dilution vessel, maintain a pressure of at least 170 kPa (25 psi) on the dilution vessel when delivering standard to the loading pump.

For sample, maintain a pressure of at least 100 kPa (15 psi) on the sample inlet line to the analyzer.

- **Leaking loading pump**

A defective piston seal allows leaks between the pump head and end plate and/or from the piston rinse ports on the pump head. Check the pump head for leaks and replace the seal if necessary.

- **Sample carryover**

Edit the CC80 Method to increase the time allowed for rinsing the sample line from the analyzer to the SS80 (or the sample source for a single-sample analyzer). In general, rinse the lines with at least five tubing volumes. For ultra-trace analyses or certain samples, larger rinse volumes may be required. The sample tubing is 1.5 mm (0.060-in) ID with a volume of about 0.55 mL/ft.

To test for sample carryover, run the highest calibration standard or most concentrated sample, followed by a blank. (For ultra-low samples, use the highest calibration standard **only**.) The carryover should be less than 1%.

The following items pertain to calibration standards.

- **Incorrect valve positions in flow path after dilution vessel**

Verify that the state of the three-way valve(s) restricts flow out of the dilution vessel while it is being filled.

- **Dilution vessel does not empty between runs**

Sample may be carrying over from injection to injection. Edit the CC80 Method to increase the time allowed for purging the dilution vessel.

- **Preparing and analyzing a more concentrated standard prior to a dilute standard**

Always prepare and analyze calibration standards in the order of increasing concentration. When a Method that involves dilution is terminated while the dilution vessel is being filled or emptied, diluted sample or standard may remain in the vessel. To purge the vessel, press the **Local/Remote** button on the CC80 front panel to select Local mode, and then press the **Drain** button. The drain operation is programmed to take 3 minutes.

If the vessel empties in less than 3 minutes, you may manually terminate the drain operation. Before doing so, check the dilution vessel waste line: If there is no liquid in the line, the vessel is empty. To terminate the drain operation, press the **Drain** button again, or press the **Gas Valve** button (this vents the dilution vessel).

- **Loading pump not operating properly**

Refer to Sections 10.10 through 10.13.

- **Poor chromatography**

Check the column manual for troubleshooting advice.

10.26 Nonreproducible Peak Area and/or Retention Time

- **Insufficient time between sample injections**

Wait until the previous sample has been completely eluted before making another injection.

- **Insufficient rinse between samples**

Edit the CC80 Method to increase the time allowed for rinsing the sample line from the analyzer to the SS80 (or the sample source for a single-sample analyzer). In general, rinse the lines with at least five tubing volumes. For ultra-trace analyses or certain samples, larger rinse volumes may be required. The sample tubing is 1.5 mm (0.060 in) ID, with a volume of about 0.55 mL/ft.

To test for sample carryover, run the highest calibration standard or most concentrated sample, followed by a blank. (For ultra-low samples, use the highest calibration standard **only**.) The carryover should be less than 1%.

- **Sample concentration too high**

Install a smaller volume sample loop. For available loop sizes, refer to the *Dionex Product Selection Guide* or contact Dionex.

- **Liquid leaks**

Refer to Section 10.2 for instructions on how to locate and eliminate leaks.

10.27 Poor Retention Time Precision

- **Leaking piston seal**

Change the piston seal as instructed in the analytical pump manual.

- **Insufficient equilibration time**

The system is equilibrated when the detector background returns to the value for the initial eluent composition. Although this usually takes 15 minutes, equilibration times vary widely, depending on the Method in use.

After a gradient elution, the system is typically equilibrated in 5 to 10 minutes. However, after a substantial eluent change (for example, after column cleaning), 30 to 45 minutes is normal.

Begin increasing the equilibration time (in increments of several minutes) until consecutive injections of a standard give reproducible retention times.

- **Malfunctioning proportioning valve (gradient pump only)**

Run an isocratic Method, using a container with the correct eluent composition. If the retention time is recovered, the proportioning valve is defective. Replace the valve as instructed in the pump manual.

- **Analytical pump requires validation**

Follow the instructions in the PeakNet-PA manual or online Help to validate the pump performance. Pump validation is included in the PeakNet-PA qualification programs. If this does not isolate the cause of the problem, contact Dionex for assistance.

10.28 Abnormal Shift in Retention Time

- **Incorrect eluent composition or concentration**

1. Prepare fresh eluent. An eluent that is excessively strong causes peaks to elute more quickly (and vice versa).

When preparing eluents, use reagent-grade chemicals or chemicals of the purity recommended in the column manual. Use only ASTM Type II (18.0 M Ω /cm resistance or 1 μ S) deionized water.

2. The gradient pump proportioning valve may be malfunctioning (applicable when the gradient pump is proportioning the eluent by combining eluents from more than one container).

Check the proportioning valve operation as follows:

Run an isocratic Method, using a container with the correct eluent composition. If the retention time is recovered, the

proportioning valve is defective. Replace the valve as instructed in the pump manual.

3. If an eluent generator is installed, verify that the correct concentration and flow rate combination is programmed in the PeakNet-PA Method. Refer to the eluent generator manual for details.

- **Incorrect flow rate through system**

1. Verify that the selected flow rate is the one specified by the analytical protocol.
2. Verify that the analytical pump is delivering the correct flow rate by measuring the eluent flow rate after the column, using a stopwatch and graduated cylinder. If the pump needs to be recalibrated, go to the **FLOW CALIBRATION** screen (from the **DIAGNOSTIC MENU**, select **CALIBRATION MENU**, and then select the **FLOW CALIBRATION** option) and follow the on-screen instructions.
3. Locate and eliminate any liquid leaks in the chromatography flow path.
4. Validate the analytical pump performance. For instructions, refer to the PeakNet-PA manual or online Help. Pump validation is included in the PeakNet-PA qualification programs.

- **System not equilibrated after eluent change**

The system is equilibrated when the detector background returns to the value for the initial eluent composition. Although this usually takes 15 minutes, equilibration times vary widely, depending on the Method in use.

After a gradient elution, the system is typically equilibrated in 5 to 10 minutes. However, after a substantial eluent change (for example, after column cleaning), 30 to 45 minutes is normal.

- **Contaminated or expended trap column**

A trap column removes trace impurities from eluents (mobile phases) and/or carbonate from hydroxide eluents. Over time, the column will become expended. Follow the instructions in the column manual to clean or regenerate the column.

Contaminated guard or analytical column

Column contamination can also lead to a loss of column capacity, which will result in poor resolution and shortened retention times.

Clean the column as instructed in the column manual. If this does not eliminate the problem, replace the column.

System temperature differs from that normally used for analysis (applicable when using a column heater)

Verify that the system temperature is correct.

10.29 Poor Linear Curve

Poor linearity is indicated by a poor correlation coefficient from the linear regression data.

- **Dilution vessel does not empty between runs**

1. Sample may be carrying over from injection to injection. Edit the Method to allow more time for the dilution vessel to drain.
2. When a Method that involves dilution is terminated while the dilution vessel is being filled or emptied, diluted sample or standard may remain in the vessel. To purge the vessel, press the **Local/Remote** button on the CC80 front panel to select Local mode, and then press the **Drain** button. The drain operation is programmed to take 3 minutes.

If the vessel empties in less than 3 minutes, you may manually terminate the drain operation. Before doing so, check the dilution vessel waste line: If there is no liquid in the line, the vessel is empty. To terminate the drain operation, press the **Drain** button

again or press the **Gas Valve** button (this vents the dilution vessel).

- **Loading pump not operating properly**

Refer to Sections 10.10 through 10.13.

- **Dilution pump not operating properly**

Refer to Sections 10.14 and 10.15.

- **Too much sample loaded onto concentrator**

Strongly-retained ions may be rinsing more weakly-retained ions off the concentrator column. Concentrate 5 mL less of sample; if linearity improves, continue decreasing the sample amount until linearity fails to improve. Increase the sample amount in increments of 5 mL, if desired, as long as linearity remains acceptable.

- **Chromatographic system overloaded**

1. Decrease the concentration of analytes in the standard.
2. Check the column manual for troubleshooting advice.

10.30 Baseline Drift

- **Incorrect regenerant flow rate**

If the baseline drifts steadily upward or downward, adjust the regenerant flow rate to level out the baseline. Decreasing the flow rate usually raises the baseline, while increasing the flow rate usually lowers the baseline.

- **Eluents or reagents improperly made**

Remake the eluent and reagent. When preparing eluents, use reagent-grade chemicals or chemicals of the purity recommended in the column manual. Use only ASTM Type II (18.0 M Ω /cm resistance or 1 μ S) deionized water. Make sure the deionized water used to prepare reagents is 18 M Ω /cm resistance.

- **Inappropriate SRS operating conditions**

Refer to the SRS manual for the correct power setting and other operating conditions.

10.31 Baseline Noise—Conductivity Detection System

- **Trapped air in detector cell**

Refer to the detector manual for corrective action.

- **Detector requires validation**

Contact Dionex for assistance.

10.32 Baseline Noise—Absorbance Detection System

- **Trapped air in detector cell**

Refer to the detector manual for corrective action.

- **Detector requires validation**

Contact Dionex for assistance.

10.33 High Background—Conductivity Detection System

A system with a high background generally has excessive noise, also, with a resulting decrease in sensitivity.

- **Incorrect eluent**

Remake the eluent. Verify that the selected flow rate is the one specified by the analytical protocol.

Check the column manual for typical background values.

If an eluent generator is installed, verify that the correct concentration and flow rate combination is programmed in the PeakNet-PA Method. Refer to the eluent generator manual for details.

- **Contaminated eluent**

Remake the eluent and reagent. When preparing eluents, use reagent-grade chemicals or chemicals of the purity recommended

in the column manual. Use only ASTM Type II (18.0 M Ω /cm resistance or 1 μ S) deionized water. Make sure the deionized water used to prepare reagents is 18 M Ω /cm resistance.

- **SRS operating incorrectly**

Refer to the suppressor manual for corrective action.

- **Contaminated or expended trap column**

A trap column removes trace impurities from eluents (mobile phases) and/or carbonate from hydroxide eluents. Over time, the column will become expended. Follow the instructions in the column manual to clean or regenerate the column.

- **Contaminated analytical column or hardware**

To determine whether the column is causing the high background, remove the column from the system.

If the background then returns to normal, either replace the column or clean it as instructed in the column manual.

If the background remains high, the hardware may be contaminated. To check for this, use deionized water as eluent. The background should be less than 2 μ S. Isolate the contaminated component by removing components from the system, one by one. Replace or clean the contaminated component.

10.34 High Background—Absorbance Detection System

A system with a high background generally has excessive noise, also, with a resulting decrease in sensitivity.

- **Detector requires validation**

Contact Dionex for assistance.

10.35 Loss of Sensitivity

- **Liquid leaks**

Tighten any leaking fittings. If a fitting continues to leak, replace it. If the connection is made with 10-32 ferrule fittings, see *Installation of Dionex Ferrule Fittings* for instructions.

- **Load/inject (LI) valve not operating correctly**

Refer to the manual included with the Rheodyne rebuild kit for troubleshooting information.

- **Sample loop not filled**

When using an injection loop, the sample flow should be sufficient to flush the loop with several loop volumes.

- **SRS needs cleaning**

Clean the suppressor as instructed in the SRS manual.

- **Contaminated concentrator**

The concentrator column may have been contaminated by impurities in the sample streams, causing a loss of capacity. When this occurs, sample ions are not concentrated as effectively. This results in less sample being injected and an apparent loss of sensitivity.

To check the concentrator for contamination:

- a. Prepare a standard consisting of two components: a peak that elutes close to the void volume and a peak that is strongly retained. For anions, use 1 ppm F⁻ and 10 ppm SO₄⁻². For cations, use 1 ppm Li⁺ and 5 ppm K⁺.
- b. Disconnect the concentrator from the load/inject (LI) valve and replace it with a 50 µL sample loop.
- c. Install the concentrator in place of the analytical column.

- d. Load the injection loop manually, using a syringe filled with the standard, and inject it. (Note: If preferred, create a Method to automate the analysis and data reduction of this capacity test.)
- e. Calculate the column capacity (see Dionex Technical Note 2R). This provides a relative measure of how the current capacity of the concentrator compares to the capacity when the concentrator is new. If the capacity has decreased by more than 30%, clean or replace the concentrator.
- f. If the concentrator column is not the cause of poor sensitivity, clean the suppressor as instructed in the SRS manual. If the problem persists, contact Dionex for assistance.

DX-800 Process Analyzer

11 • Service

This chapter describes service and repair procedures that the user may perform. Any procedures not included here, including electronics-related service procedures, must be performed by Dionex personnel. For assistance, contact Dionex Technical Support. In the U.S., call 1-800-346-6390. Outside the U.S., call the nearest Dionex office.

Before replacing any part, refer to the troubleshooting information in Chapter 10 to isolate the cause of the problem. This chapter describes routine service procedures that the user can perform. All other procedures must be performed by Dionex personnel.



The CPU card contains a lithium battery. If it is necessary to replace the CPU card, dispose of the used battery according to the manufacturer's instructions.

IMPORTANT

Substituting non-Dionex parts may impair DX-800 performance, thereby voiding the product warranty. Refer to the warranty statement in the Dionex Terms and Conditions for more information.

11.1 Eliminating a Fluid System Restriction

A restriction in the fluid system (crimped tubing, etc.) can cause excessive system backpressure. This, in turn, may cause leaks or irreparable damage to system components.

1. Begin pumping eluent (mobile phase) through the system (including the columns) at the flow rate normally used.
2. Refer to the appropriate fluid schematic (P/N 051833 for the SP81, P/N 051860 for the SP82, P/N 051861 for the SP83, P/N 051862 for the SP84). Work backward through the system, beginning at the cell exit. One at a time, loosen each fitting and check the pressure. The connection at which the pressure drops indicates the point of restriction.
3. Remove the restriction, either by flushing or by replacing the section of tubing.

11.2 Cleaning/Replacing Precision Displacement Pump Check Valves

The loading pump and dilution pump are precision displacement pumps. Depending on the DX-800 configuration, one or both of these pumps is installed on the SP80 Sample Preparation panel. A dirty or worn check valve will cause an erratic flow rate and prevent the pump from delivering the expected volume of liquid.

Removing the Inlet Check Valve

1. Open the front door of the AE80 enclosure. Press the **POWER** button on the CC80 front panel to turn off the power.
2. Disconnect the Teflon fitting from the inlet check valve (see Figure 11-1).
3. Use a 1/2-inch wrench to loosen the check valve housing. Remove the housing and carefully remove the check valve cartridge from the housing.




Figure 11-1. Precision Displacement Pump Check Valves

Removing the Outlet Check Valve

1. Open the front door of the AE80 enclosure. Press the **POWER** button on the CC80 front panel to turn off the power.



To prevent the pump from inadvertently starting, unplug it from the distribution board and turn off the air supply to the pump.



Pour empêcher la pompe de démarrer par mégarde, débranchez-la du tableau de distribution et coupez-en l'alimentation en air.

DX-800 Process Analyzer



Ziehen Sie den Stecker auf der Verteilerkarte und schalten Sie die Luftzufuhr zur Pumpe ab, um zu verhindern, daß die Pumpe versehentlich startet.

2. Disconnect the fitting from the outlet check valve.
3. Use a 1/2-inch wrench to loosen the check valve housing. Remove the housing and carefully remove the check valve cartridge from the housing.

Cleaning the Check Valves

1. Place the check valve housings and cartridges in a beaker with methanol. Sonicate or agitate for several minutes.
2. Rinse each check valve housing and cartridge thoroughly with filtered, deionized water.

Replacing the Inlet Check Valves

1. The inlet check valve assembly housing has a 1/4-28 port. Replace the cartridge in the inlet check valve housing so that the double-hole end of the cartridge is visible. Liquid flows through the check valve in the large single hole and out the small double holes.
2. Reinstall the check valve. Tighten only enough to seat (25 in-lb torque). If the cartridge leaks again, it is cracked and should be replaced.



Overtightening may damage the pump head and the check valve housing and crush the check valve seats.

Replacing the Outlet Check Valves

1. The outlet check valve assembly housing has a 10-32 port. Replace the cartridge in the outlet check valve housing so that the single-hole end of the cartridge is visible. Liquid flows through the check valve in the large single hole and out the small double holes.

2. Reinstall the check valve. Tighten only enough to seat (25 in-lb torque). If the cartridge leaks again, it is cracked and should be replaced.



Overtightening may damage the pump head and the check valve housing and crush the check valve seats.

Completing the Procedure

1. Reconnect the liquid lines.
2. Press the **POWER** button on the CC80 to restore power. Close the front door of the AE80 enclosure.
3. Prime the system. If the system will not prime and all other possible causes of the problem have been eliminated, replace the check valve cartridge (P/N 047755).

11.3 Changing a CC80 Fuse

1. Press the **Emergency Off** button on the front door of the AE80 enclosure to turn off the main power.



HIGH VOLTAGE—Disconnect the main power cord from its source and also from the **POWER IN** connector on top of the DX-800 enclosure.



HAUTE TENSION—Débranchez le cordon d'alimentation électrique principale de sa source et aussi du connecteur **POWER IN** (Entrée du courant) sur le dessus de l'enceinte du DX-800.



HOCHSPANNUNG—Ziehen Sie das Netzkabel aus der Steckdose und der **POWER IN**-Buchse oben auf dem Gehäuse des DX-800.

DX-800 Process Analyzer

2. Open the enclosure door. Loosen the two retaining screws in the module chassis. Grasp the chassis by the sides and pull it forward just until the first stop on the sliders is engaged.



Do not pull the chassis beyond the slider stops. The chassis can become disengaged from the enclosure.




Ne tirez pas le châssis au-delà des butées du mécanisme coulissant. Le châssis peut se détacher de l'enceinte.



Ziehen Sie die Montageplatte nicht über die Arretierungen hinaus. Sie könnte sich sonst vom Gehäuse lösen.

3. Use a Phillips screwdriver to loosen, but not remove, the screw in each corner of the top cover of the CC80.
4. The fuse holders are near the front left corner of the CC80 (see Figure 11-2). One at a time, remove each of the five fuse holders from the DIN rail and check the fuse with an ohmmeter. If the fuse is not open, reinstall the fuse holder. If the fuse is open, replace it with a new fuse.




*Figure 11-2. Changing a CC80 Fuse
CC80 Top Cover Removed—Top View*

To replace a fuse:

- a. Hold the fuse holder with the cover (which bears the words “Open Here”) facing upward. Insert the end of a small screwdriver into the recessed lock of the cover and pop it open.
 - b. Remove the old fuse (see Figure 11-3) and insert a new 3.15 amp fast-blow IEC127 fuse (P/N 954745). Close the cover.
 - c. Hold the fuse holder with the cover facing to the right and push it into the DIN rail. The fuse holder is keyed to fit only in its proper orientation.
5. Replace the top cover of the CC80. Push the module chassis back into the enclosure and tighten the retaining screws.

DX-800 Process Analyzer

6. Close the enclosure door. Reconnect the main power cord and push the **Power Reset** button to restore power.



*Figure 11-3. Changing a CC80 Fuse
CC80 Top Cover Removed--Cutout View*

12 • TTL and Relay Control

The CC80 provides eight TTL inputs and two relay outputs. The TTL and relay connectors are on the CC80 rear panel. An internal cable connects the CC80 connectors to the **TTL/Relay** connector on top of the AE80 enclosure. See Section 12.3 for connection instructions.

- The CC80 TTL inputs allow external devices to activate various analyzer functions (see Section 12.1).
- The CC80 relay outputs allow the CC80 to control functions in external devices (see Section 12.2).

12.1 TTL Input Control

The eight TTL inputs allow an external device to trigger one or more of the following analyzer actions:

- Display an alarm on the CC80 front panel
- Turn on a relay
- Put the channel in standby
- Shut down the channel
- Bypass an injection from a specified sample.

For example, you can connect a flow sensor on a sample pipe to one of the TTL inputs. If the sample stops flowing, the sensor signals the TTL and the analyzer bypasses that sample source.


12.1.1 Configuring TTL Input Actions

NOTE

Refer to the PeakNet-PA user's guide or online Help for details about configuring analyzers.

PeakNet-PA is used to assign actions to CC80 TTL inputs. To do this:

1. Open the PeakNet-PA Analyze program and select **Analyzer Configuration** from the **Edit** menu.
2. Select the analyzer and then click the **Edit** button.
3. Select the tab for the channel whose TTL inputs you are configuring.



4. On the channel tab page, enter the desired action(s) for each connected TTL input. **Note:** Cells that are gray cannot be edited.
 - **CC80 Alarm:** Enter a **1**, **2**, **3**, or **4** to turn on the corresponding CC80 front panel alarm light.

- **Relay Out:** Enter a **1** or **2** to open the corresponding relay output.
- **Standby:** Select **Yes** to place the channel in standby.
- **Shutdown:** Select **Yes** to shut down the channel.
- **Bypass:** Select a sample valve position from the drop-down list to skip injections for that sample.

NOTE

If the relay outputs are configured to change states with alarm events, do not configure them in CC80 Methods to control devices. If they are configured with alarms, always include a sample preparation step at the beginning of the CC80 Method that closes the relays (see Section 12.2). Include the step in all CC80 Method files (*.met) or sample preparation files (*.spr) used with the channel.

12.1.2 TTL Input Signal Mode

On CC80s with Moduleware Version 5.26 and later installed, the TTL inputs are, by default, active on a positive edge signal (see Figure 12-1). If the input goes from low (0 V, closed) to high (+5 V, open), the configured action occurs. This follows the logic that a closed circuit with continuity is normal, and an open circuit is an exception (alarm).

NOTE

To check the CC80 Moduleware version, open the PeakNet-PA Configuration Editor and click the Moduleware button on the toolbar.

For the CC80 alarm, relay control, and sample bypass actions, the action remains on as long as the TTL input is high. If the input returns to low, the action is turned off. For example, if an external device triggers a CC80 front panel alarm, the alarm stays on as long as the input is at +5 V (open). When the input returns to 0 V (closed), the alarm goes off. For the channel standby and shutdown actions, returning the TTL to low has no effect.

DX-800 Process Analyzer

In the Run program, on the channel's system window, the arrow icon for each CC80 TTL input indicates the TTL state:

 indicates the TTL input is closed.

 indicates the TTL input is open.

The table below summarizes the TTL input mode functions:

TTL Inputs	Input Signal Mode	Voltage	PeakNet-PA Display	PeakNet-PA Programmed Action
1-8	Positive Edge	0 → +5		Turned on
		+5 → 0		Turned off* or no effect**

*The CC80 alarm, relay, and valve bypass actions are turned off.

**The standby and shutdown actions are not affected.




Figure 12-1. Positive Edge TTL Input Signal Mode

NOTE

A DIP switch can be used to change the TTL input signal mode. For details, see Section 12.4.

12.2 Relay Output Control

Relay outputs 1 and 2 can be programmed to switch any low-voltage control. Switched current must be less than 200 mA and 60 V peak blocking.




Relay loads in excess of 200 mA or with included power supplies over 60 V may damage the relay drivers on the CPU module.


PeakNet-PA is used to control the relay output states. You can (1) change the relay's state immediately (Direct Control), (2) program the relay's state change into a CC80 Method (Method Control), or (3) configure the relay's state change as a function of an alarm condition or TTL input (see Section 12.1.1).

To control relay output states using Direct Control:

1. On the InTouch Analyzer Status screen, click the picture of the channel whose relays you are setting.




2. On the Run program system window, select the **Direct Control** toolbar button.



Direct Control button

DX-800 Process Analyzer

3. The Direct Control dialog box appears. Select the **CC80** tab.



4. Click the **Relay1** and/or **Relay2** check boxes to close or open the relays. With the default settings for CC80s with Moduleware Version 5.26 and later installed, checking the box closes the relay (current flows) and clearing the box opens the relay (no current flows).
5. Click **OK** or **Apply**.


To control relay output states using Method Control:

1. Open the PeakNet-PA Method program and open the Method that will control the relay outputs.
2. In the Method window, double-click the CC80 module. The CC80 Preparation Steps dialog box appears.
3. To add the relay output control step to the end of the Method, select the blank step at the end of the list. Otherwise, select the step below where you want to add the new relay control step.

4. Select **Relay Output** from the **Function** drop-down list. Click the **Relay1** and/or **Relay2** check boxes to open or close the relays. With the default settings for CC80s with Moduleware Version 5.26 and later installed, checking the box closes the relay (current flows) and clearing the box opens the relay (no current flows). In the list of preparation steps, a closed relay (checked box) is indicated as “On.”
5. Click the **Insert** button.

NOTE

If the relay outputs are configured to change states with alarm events (see Section 12.1.1), do not configure them in CC80 Methods to control devices. If they are configured with alarms, always include a sample preparation step at the beginning of the CC80 Method that closes the relays. Include the step in all CC80 Method files (*.met) or sample preparation files (*.spr) used with the relay’s channel.



12.3 TTL and Relay Connections

Two TTL input and relay output connectors (P1 and P2) are on the CC80 rear panel. An internal cable connects the CC80 connectors to the **TTL/Relay** connector on top of the AE80 enclosure (see Figures 12-2 and 12-3).




Figure 12-2. Top View of AE80 Enclosure Electrical I/O Panel with BNC Connectors




Figure 12-3. Top View of AE80 Enclosure Electrical I/O Panel with RJ-45 Connector

To connect a TTL input or Relay output to an external device:

1. Connect the DX-800 TTL/Relay cable (P/N 052899) (see Figure 12-4) to the **TTL/Relay** connector on the AE80 top panel.




Figure 12-4. TTL/Relay Cable and Connector Pinout

2. Connect the required wires (refer to the table below) to the TTL or relay connector on the external device. Check the polarity of the connection. Connect signal wires to signal (+) pins and ground wires to ground (-) pins.

TTL/Relay		Wire Color	TTL/Relay		Wire Color
TTL Input 1	+	Orange	TTL Input 6	+	Purple
	-	Red/Black		-	Black
TTL Input 2	+	Red	TTL Input 7	+	Green
	-	Red/Black		-	Black
TTL Input 3	+	Blue	TTL Input 8	+	Tan
	-	Red/Black		-	Black
TTL Input 4	+	Brown	Relay Output 1	+	White
	-	Red/Black		-	Gray
TTL Input 5	+	Yellow	Relay Output 2	+	Pink
	-	Red/Black		-	White/Red

12.4 Alternate TTL Input Signal and Relay Output Modes

The functionality of the CC80 TTL inputs and relay outputs changed with the release of PeakNet-PA Version 5.21 and CC80 Moduleware Version 5.26. The new default modes are appropriate for nearly all installations. However, exceptions do exist and the information in this section is provided for reference if needed.

DIP switch #2 on the rear of the CC80 front panel controls the TTL input signal and relay output modes.

TTL Input Signal Modes

DIP switch #2, position 2 (SW2.2), controls the signal mode for the TTL inputs. When SW2.2 is ON (the default), all eight TTL inputs are active on a positive edge signal, as described in Section 12.1.2.

When SW2.2 is OFF, TTL inputs 1 – 4 are active on a negative edge and TTL inputs 5 – 8 are active on a negative pulse.

NOTE

In many Dionex modules, such as a DX-600 pump, the negative edge and negative pulse modes are referred to as *normal* edge and *normal* pulse. The positive edge mode is referred to as *inverted* edge. To comply with industrial installation requirements, the CC80 is configured to function differently.

Relay Output Modes

DIP switch #2, position 3 (SW2.3), controls the mode for the relay outputs. When SW2.3 is OFF (the default), the relay contact closures are normally open. Activating a relay closes it (current flows).

When SW2.3 is ON, the relay contact closures are normally closed. Activating a relay opens it (no current flows).

NOTE

In the PeakNet-PA Run program, the graphical display for the relays is the same for both switch settings. A closed relay icon means the relay is off. An open relay icon means the relay is on.

A • Specifications

NOTE

For specifications for the analytical pump, detector, eluent generator, and SRS, refer to the relevant user manual.

For installation requirements, refer to *Installation Requirements and Customer Responsibilities* (Document No. 031176).

A.1 Electrical

Main Power Requirements	100 to 120 Vac, 15 A, 50/60 Hz maximum; 220 to 240 Vac, 7.5 A, 50/60 Hz maximum The power supply is main voltage auto-sensing and requires no adjustment.
Air Conditioner or Blower	115 Vac, 12 A, 60 Hz maximum
Fuse Requirements	Five 3.15 amp fast-blow IEC127 fuses (P/N 954745)

A.2 Environmental

Ambient Operating Temperature	(with ventilation) 4 to 40 °C (40 to 105 °F)
	(with air conditioner) 10 to 55 °C (50 to 130 °F)
Operating Humidity	5 to 95% relative humidity, noncondensing

DX-800 Process Analyzer

A.3 Physical

Dimensions	(with ventilation) 93 x 58 x 58 cm (37 x 23 x 23 in)
	(with purge-and-pressurization unit) 118 x 72 x 58 cm (46 x 28 x 23 in)
	(with air conditioner) 93 x 81 x 58 cm (37 x 32 x 23 in)
	1 m (40 in) of clearance in front of the enclosure is required
Weight	(with ventilation) 124 kg (275 lb)
	(with purge-and-pressurization unit) 132 kg (290 lb)
	(with air conditioner) 168 kg (370 lb)
Decibel Level	75 db (at “A WEIGHING” setting)

A.4 Pumps

A.4.1 Dilution Pump

Type	Single-piston
Operating Pressure	3.5 to 21 MPa (500 to 3000 psi)
Flow Rate	15 mL/min maximum

A.4.2 Loading Pump

Type	Single-piston
Operating Pressure	3.5 to 21 MPa (500 to 3000 psi)
Flow Rate	3 mL/min maximum

A.4.3 PC80 Post-Column Reagent Pump (Optional)

Type	Single-piston
Operating Pressure	3.5 to 13 MPa (500 to 1900 psi)
Flow Rate	0.2 to 1 mL/min maximum

A.5 Valves

A.5.1 Check Standard (CS) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.2 Column Switching Valve (Optional)

Type	Electrically-actuated, 10-port Rheodyne 9650E valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.3 Diluent (DI) Valve

Type	2-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

DX-800 Process Analyzer

A.5.4 Diluent Select (DS) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.5 Dilution Vessel (DV) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.6 Gas Valve

Type	3-way, electrically-actuated gas/liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.7 Load/Inject (LI) Valve

Type	Electrically-actuated, 6-port Rheodyne 9750E06 valve
Operating Pressure	30 MPa (4000 psi), maximum

A.5.8 Metering (ME) Valve

Type	Electrically-actuated, 10-port Rheodyne 9650E liquid valve
Operating Pressure	30 MPa (4000 psi), maximum

A.5.9 Sample Select (SM) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.10 Sample/Standard (SS) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.5.11 Standard (ST) Valve

Type	3-way, electrically-actuated liquid solenoid valve
Operating Pressure	0.7 MPa (100 psi), maximum

A.6 Dilution Vessel

Material	PFA (perfluoroalkoxy) Teflon
Capacity	250 mL
Pressure Relief Valve	Opens at 0.34 MPa (50 psi)

A.7 SS80 Sample Selector (Optional)

Valve(s)	Electrically-actuated, multiport PEEK valves from Rheodyne rated to 5.5 MPa (800 psi); 1, 2, and 3 valves required for control of 7, 14, and 21 sample sources, respectively
Dimensions	50 x 26 x 25 cm (20 x 11 x 10 in)
Weight	14 kg (30 lb)

A.8 Column Heater (Optional)

Power Rating	30 W
Materials	Anodized aluminum
Operating Temperature	From ambient + 5 °C up to 80 °C ± 1 °C
Dimensions	Accommodates one 6 or 8 mm OD x 100, 150, or 250 mm ID column

A

- Absorbance detector, 1-13
 - background noise, 10-35
 - cell, 1-10
 - cell contains trapped air, 10-34
 - front panel display, 1-13
 - lamps, 1-13
 - manual, 1-13
 - performance validation, 10-35
- AC outlets, CC80 rear panel, 3-5
- AE80 Enclosure, 1-9, 2-1
 - See also* Enclosure
- Air conditioner, 2-7
- Air consumption, excessive, 10-6
- Air filters, cleaning, 2-7
- Alarm conditions, 3-1, 3-4
 - warning tone, 3-4
- Analytical column
 - See* Columns
- Analytical pump, 1-1, 1-12
 - connection to CC80 rear panel, 3-5
- Analyzer
 - definition, 1-1
 - leaks, 3-4
 - PeakNet-PA control, 1-16
- Analyzer Leak LED, 2-7
- Auto offset, 1-13

B

- Background, high, 4-11
 - absorbance detector, 10-35
 - conductivity detection system, 10-34
- Backpressure, system
 - excessive backpressure, 10-2, 10-7
 - recommended value, 10-7
- Baseline drift, 4-11
- Baseline noise, 10-34
- BIOS version number, 3-10
- Blower, 2-7
- Bottles, 7-1

C

- Calibration
 - precision displacement pump(s), 9-4
- Calibration standards, 1-9, 4-1, 4-3, 4-5, 4-9
 - order of analysis, 10-28
 - pressure requirements, 10-27
- Calibration Wizard, 4-10
- Calibration, standard loop, 9-5 - 9-6
- Cartridge, eluent generator, 1-10, 5-4, 6-2
- CC80 Component Controller, 1-2, 1-9, 3-1
 - diagnostic tests, 3-10
 - DIP switches, 3-6
 - front panel description, 1-9
 - functions, 1-9, 3-1
 - operating modes, 3-9
- Cell, absorbance, 1-10
 - trapped air, 10-34
- Cell, conductivity, 1-10, 5-3
 - thermal control, 1-13
 - trapped air, 10-34
- Cell, electrochemical, 5-3
- Channel
 - definition, 1-1
 - diluent supply, 10-16
 - does not run, 10-8
 - features, 1-9
 - initial startup procedure, 9-3
 - PeakNet-PA control, 1-16
 - waste lines, 10-3
- Check standard (CS) valve
 - specifications, A-3
- Chromatography components, 5-1
- Circuit breaker, CC80, 3-5, 10-10
- Column capacity
 - how to calculate, 10-37
 - loss of, 10-25, 10-32, 10-36
- Column heater, 1-10, 3-1
 - advantages of, 1-15
 - columns accommodated, 1-15
 - control of, 1-9
 - description, 5-4
 - does not heat, 10-20

DX-800 Process Analyzer

- specifications, A-6
- temperature controls, 1-9
- temperature range, 1-15
- Column switching valve, 1-10, 1-15
 - description, 5-4
 - specifications, A-3
- Columns, 1-10
 - See also* Concentrator column
 - contamination of, 4-11, 10-25, 10-32, 10-35
 - degradation of, 10-24
 - excessive flow rate through, 10-7
 - mounting, 5-3
 - overloading, 10-26
- Compressed gas cylinders, 1-6
- Concentration configuration, 1-1, 4-2
 - CC80 flow chart, 4-2
 - components installed, 4-3
- Concentration with reagent addition configuration, 1-1, 4-8
 - CC80 flow chart, 4-8
 - components installed, 4-9
- Concentrator column, 4-10, 5-3, 10-33
 - capacity, 10-37
 - contamination, 10-36
 - in trace ion analysis, 4-2
- Conditional responses, 1-1
- Conductivity detection, 1-10, 1-14
 - baseline noise, 10-34
 - high background, 10-34
- Conductivity detector
 - cell, 1-10, 5-3
 - front panel display, 1-13
- Configurations, 4-1
- Configurations, DX-800, 1-1
- Connections
 - DX-LAN, 9-2
 - electrical, 9-1
 - sample inlet, 9-2
 - water sources, 9-2
- Consumables, 5-3
- Contacting Dionex, 3-10
- Controller board
 - dilution pump, 10-17
 - precision displacement pumps, 10-12
- CPU
 - battery disposal, 11-1
 - diagnostic test, 3-10

D

- DC amperometry, 1-14
- Degas tubing assembly, 6-2
- Deionized water
 - for cleaning, 1-6, 11-4
 - for eluents, 10-25, 10-30, 10-35
 - troubleshooting high background, 10-35
- Detection
 - absorbance, 1-13
 - amperometry, 1-14
 - conductivity, 1-13
 - electrochemical, 1-14
 - ultraviolet, 1-13
 - visible, 1-13
- Detector cell
 - See also* Cell
 - installation site, 5-3
- Detectors, 1-1
 - absorbance, 1-13
 - conductivity, 1-13
 - connection to CC80 rear panel, 3-5
 - descriptions of, 1-13
 - electrochemical, 1-14
- Diagnostics, CC80, 3-10
- Digital autoranging, 1-13
- Diluent
 - no diluent supplied to channel, 10-16
 - pressure requirements, 10-17
- Diluent (DI) valve
 - specifications, A-3
 - verifying correct operation, 10-16
- Diluent select (DS) valve
 - specifications, A-4
- Dilution configuration, 1-1, 4-4
 - CC80 flow chart, 4-4
 - components installed, 4-5
- Dilution pump, 4-10
 - air leak, 10-26
 - calibration procedure, 4-10
 - cleaning and replacing check valves, 11-2 - 11-5
 - controller board, 4-10, 10-17
 - does not pump, 10-16
 - electrical connections, 10-13
 - gas pressure requirements, 10-16
 - leaks, 3-4, 10-3

- operation, 4-10
 - piston seal, 10-3
 - priming, 10-12 - 10-13
 - specifications, A-2
 - Dilution vessel
 - air pressure requirements, 10-13
 - description, 4-11
 - does not empty, 10-18, 10-32
 - draining, 10-4
 - fails to empty between runs, 10-28
 - gas leak, 10-15
 - gas pressure requirements, 10-4
 - inadequately pressurized, 10-18
 - liquid leaks, 10-4
 - pressure relief valve, 4-11
 - pressure requirements, 4-11
 - specifications, A-5
 - waste line, 10-4, 10-8, 10-19, 10-28, 10-32
 - Dilution vessel (DV) valve
 - malfunction, 10-20
 - specifications, A-4
 - Dilution with reagent addition configuration, 1-1, 4-6
 - CC80 flow chart, 4-6
 - components installed, 4-7
 - Dionex Technical Support, 3-10, 10-1
 - DIP switches, CC80, 3-6, 3-8
 - resetting, 3-6
 - sample select valve settings, 10-14
 - Display Refresh button, 2-3
 - Distribution board, 10-22
 - LC80, 5-1, 10-12
 - SP80, 10-12 - 10-13
 - Drip tray, analyzer
 - emptying, 2-6
 - leak sensor, 2-7
 - DS3 Detection Stabilizer, 5-3
 - Dual-column system, 1-15
 - See also* Column switching valve
 - DX-500 modules
 - display brightness, 2-2
 - DX-800 Process Analyzer, 1-9
 - configurations, 1-1
 - description, 1-1
 - electrical connections, 2-5 - 2-6
 - initial power-up, 2-3
 - options, 1-14 - 1-15
 - physical specifications, A-2
 - DX-LAN interface, 1-16, 3-9
 - electrical cable connection, 3-5, 9-2
 - eluent generator connections, 6-4
 - PeakNet-PA workstation connections, 2-6
- E**
- EEMAC, 2-1
 - EG40
 - See* Eluent generator
 - Electrical connections, 2-5 - 2-6, 9-1
 - for LC80 components, 5-1
 - precision displacement pumps, 10-13
 - Electrical precautions, 1-7
 - Electrical specifications, A-1
 - Electrochemical detector, 1-13 - 1-14
 - cell, 5-3
 - detection modes, 1-14
 - front panel display, 1-13
 - Electronics
 - EG40-PA, 6-4
 - Eluent generator, 1-10
 - cartridge, 5-4
 - main components, 6-2
 - overview, 6-1
 - power supply, 6-4
 - power switch, 6-4
 - TTL connectors, 6-4
 - Eluents, 7-2
 - See also* NOWPAK containers
 - eluent generator, 1-10
 - problem with excessively strong eluent, 10-25
 - purity requirements, 10-25, 10-30, 10-33 - 10-34
 - EluGen cartridge, 1-10
 - cartridge lifetime, 6-4
 - on LC80 panel, 5-4, 6-2
 - replacing, 6-4
 - Emergency Off switch, 2-2
 - Emergency shutdown
 - restoring power, 2-3
 - Enclosure, 1-9
 - clearance required, 2-1
 - dimensions, 2-1
 - drip tray, 2-6 - 2-7

DX-800 Process Analyzer

- liquid and gas connections, 2-6
- NEMA standards, 2-1
- types, 2-1
- Environmental specifications, A-1
- Equilibration time, 10-31

F

- Ferrule fittings, 10-3
 - manual, 1-5
 - overtightened, 10-17
- Fittings
 - repairing leaking fittings, 10-2, 10-6, 10-12
- Flow rate, erratic, 11-2
- Flushing the flow path, 9-3
- Fuses
 - replacement instructions, 11-5

G

- Gas connection
 - LM80 Liquids Manager, 7-2
- Gas consumption, excessive, 10-6
- Gas valve
 - specifications, A-4
- Gradient pump, 1-12
 - calibration, 10-7
 - front panel display, 1-12
 - performance validation, 10-25
 - test of proportioning valve operation, 10-25
- Guard columns, 5-3
 - See also* Columns

H

- Heater, column, 5-4
 - See also* Column heater
- Helium leaks, 10-6, 10-18
- Hydraulic specifications, A-2

I

- Installation, 9-1
- Integrated amperometry, 1-14
- InTouch, 1-1 - 1-2
- Ion chromatography modules
 - detectors, 1-13
 - pumps, 1-12
- Isocratic pump, 1-12
 - calibration, 10-7
 - front panel display, 1-12
 - manual, 1-12

L

- LC80 Liquid Chromatography panel, 1-10
 - components, 1-10, 1-15
 - distribution board connections, 10-20, 10-22
- Leak detection
 - air and gas, 10-6
- Leak sensor
 - calibration instructions, 10-12
 - SS80, 8-2
- Leaks
 - air and gas, 10-6
 - gas leak from dilution vessel, 10-15
 - helium, 10-18
 - liquid, 2-6, 10-2 - 10-5
- Linearity
 - improving, 10-26, 10-32 - 10-33
 - symptoms of poor linearity, 10-32
- LM80 Liquids Manager, 1-11, 7-1
 - control panel, 7-1 - 7-2
 - drip tray for NOWPAK containers, 1-11
 - pressure relief valve, 7-2
- Load/inject (LI) valve, 1-10, 5-2
 - description, 5-3
 - diagnostic test, 3-10
 - specifications, A-4
 - verify correct operation, 10-22
- Loading pump, 4-10
 - air leak, 10-26
 - calibration procedure, 4-10
 - cleaning and replacing check valves, 11-2 - 11-5
 - controller board, 4-10, 10-14

- electrical connections, 10-13
- leaks, 3-4
- no sample delivered, 10-13
- operation, 4-10
- priming, 10-12 - 10-13
- pump head leaks, 10-27
- specifications, A-3

Local mode, 3-9

Log, maintenance, 10-1 - 10-2

Loop

- See* Sample loop

M

Maintenance

- log, 10-1 - 10-2
- operational, 9-13
- routine, 9-13

Manuals

- Dionex products, 1-5
- NOWPAK containers, 1-5

Metering (ME) valve, 4-10

- diagnostic test, 3-10
- leaks, 10-17
- specifications, A-4

Methods

- downloading, 3-9
- dual-column, 5-4
- ended before dilution vessel drains, 3-1
- when aborted, 3-9

Mobile phase, 7-2

- See also* Eluents

Moduleware, CC80

- downloading, 3-10
- version number, 3-10

N

NEMA enclosures, 2-1

NOWPAK, 1-2

NOWPAK containers, 1-11, 7-2

- See also* Eluents
- documentation, 1-5
- gas pressure, 7-2 - 7-3
- maximum number, 1-11

O

Operation

- routine, 9-10

P

PC80 Post-Column Pump, 1-15, 3-1, 4-11

- control of, 1-9
- specifications, A-3
- versions available, 4-11

PeakNet-PA software, 1-1, 1-16

- alarm conditions, 3-1, 3-4
- Analyze program, 3-9
- Configuration Editor, 4-10
- features, 1-1
- initial setup, 9-8
- manual, 1-5, 1-16
- methods, 3-9
- online Help, 1-16
- startup, 9-3
- user interface, 1-1

PeakNet-PA workstation, 1-16, 2-6

Peaks

- nonreproducible peak area, 10-29
- spurious, 10-23 - 10-24
- undetected by PeakNet-PA, 10-21

Pneumatic requirements, 7-3

Post-column applications, 1-15

Post-column reagents, 1-15, 4-11

Power cord

- DX-500 modules, 10-9
- post-column pump, 10-9
- system, 10-9

Power supply

- eluent generator, 6-4

Power switch

- eluent generator, 6-4
- emergency off switch, 2-2
- power reset, 2-3

Precision displacement pumps, 4-10, 10-7

- See also* Dilution pump

Precision, poor, 10-26, 10-29

Pressure relief valve

- dilution vessel, 4-11
- LM80, 7-2

DX-800 Process Analyzer

SP80, 10-19
Problems, help for, 10-1
Proportioning valve (gradient pump), 10-30
 verification of correct operation, 10-25
Pulse damper, 4-11
Pump name, 1-1
Pumps
 See also Gradient pump, Isocratic pump
 dilution pump, 4-10
 loading pump, 4-10
 PC80 Post-Column Pump, 1-15, 3-1, 4-11
 precision displacement, 4-10
Purge-and-pressurization unit
 installation site, 2-8
Purging the dilution vessel, 10-4

R

Reaction coil, 4-11
Reagent Pump (RP-1), 4-11
Reagents, 1-11
 adding to samples, 4-6, 4-8
 post-column, 1-15, 4-11
 preparation of, 10-33, 10-35
 reservoirs, 7-1
 slow delivery, 10-6
Remote mode, 3-9
 selecting, 3-9
Reservoirs, 1-11
 gas, 7-3
 slow liquid delivery, 10-6
Reset button, 2-2 - 2-3, 11-7
Resolution, poor, 10-32
Restriction in hydraulic system, 11-2
Retention time
 decrease in, 10-24, 10-32
 shift in, 10-30
Routine operation, 9-10

S

Safety messages, 1-3
Safety practices, 1-6
Sample
 See also SS80 Sample Selector

 concentration too high, 10-29
 excess sample loaded, 10-33
 multiple, 8-1
 no sample injected, 10-22
 not delivered to loading pump, 10-13
 test for carryover, 10-23, 10-27
Sample (SM) valve
 specifications, A-4
Sample inlet line connections, 9-2
Sample loop, 4-10, 5-3
 not filled, 10-36
 sizes available, 10-29
Sample preparation, 4-1
Sample select valve
 malfunction, 10-11, 10-14
 switching positions, 10-11
Sample streams, 1-14
Sample/standard (SS) valve
 specifications, A-5
Sampling leak, 10-5
Sampling system, 1-1
Schedules
 downloading, 3-9
Sensitivity, loss of, 10-34, 10-36
Sensor, leak, 2-7
Separator column
 See Columns
Service procedures
 eliminating fluid system restrictions, 11-2
 for precision displacement pump check
 valves, 11-2 - 11-5
 fuse replacement, 11-5
Shutdown
 long-term, 9-12
 short-term, 9-11
Snoop, 10-6, 10-15
Solvent, 7-2
SP80 Sample Preparation panel, 1-9, 4-1
 access to rear of panel, 4-1
 components on front of panel, 1-9
 components on rear of panel, 1-10, 1-15
 concentration configuration, 4-1 - 4-2
 concentration w/ reagent configuration, 4-8
 concentration w/reagent configuration, 4-1
 dilution configuration, 4-1, 4-4
 dilution with reagent configuration, 4-1, 4-6
 pressure relief valve adjustment, 10-19

Specifications
 column heater, A-6
 dilution vessel, A-5
 electrical, A-1
 environmental, A-1
 hydraulic, A-2
 physical, A-2
 SS80 Sample Selector, A-5
 valves, A-3

SRS (Self-Regenerating Suppressor), 5-3
 See also Suppressor
 cleaning, 10-37
 troubleshooting, 10-8, 10-35
 waste line, 10-3

SS80 Sample Selector, 1-14, 8-1
 CC80 controls, 3-3
 connection to analyzer, 2-6
 electronics, 8-2
 installation site, 8-2
 leak sensor, 8-2
 liquid leaks, 10-5
 valves, 3-1, 8-2

Standard (ST) valve
 specifications, A-5

Standard loop
 calibration, 9-5 - 9-6
 volume, 9-5

Standards, 1-11
 reservoirs, 7-1
 slow delivery, 10-6

Startup
 initial, 9-3
 routine, 9-9

Suppressor, 1-10, 5-3
 regenerant shutoff valve, 5-3

Switches, CC80, 3-6, 3-8

System backpressure
 excessive, 10-2, 10-7, 11-2
 recommended value, 10-7

System configuration, 1-1

T

Temperature compensation, 1-13
Trace ion analysis, 4-2

Trademarks, 1-2
Transition metals applications, 4-8, 7-3
Trap column, 10-32
TTL control
 cables, 2-6
 eluent generator, 6-4
Tubing
 compromised by chemical fumes, 10-6
 damaged, 10-5
 replacing, 10-3
 Teflon, 10-3, 10-17

U

Ultra-trace analyses
 rinse between samples, 10-23

V

Valves
 See also Valve names
 column switching, 5-4
 concentration configuration, 4-3
 concentration with reagent configuration, 4-9
 dilution configuration, 4-5
 dilution with reagent configuration, 4-7
 precision displacement pump check valves,
 11-2 - 11-5
 sample select, 8-1 - 8-2
Voltammetry, 1-14

W

Waste disposal, 1-6
Waste lines
 cell, 10-3
 channel, 10-3
 problems, 10-3
Water connections, 9-2
Windows NT, 1-1 - 1-2
Wonderware, 1-1
 See also InTouch

DX-800 Process Analyzer
