



The Eight-Position Terminal Strip on the Rear of the LSC-2

### PIN OUTS

#### LSC-2: SOCKET

3	GRN	GROUND
2	WHT	READY TO PURGE OUTPUT
1	RED	READY TO PURGE OUTPUT

#### LSC-2: TERMINAL STRIP

8	ORG/BRN	ADVANCE TO PURGE INPUT
---	---------	------------------------

#### LSC-2: T-2 CONNECTOR

1	BLK	BEGIN DESORB OUTPUT
2	BLU	BEGIN DESORB OUTPUT

#### AQUATEK 50

5	GRN	GROUND
5	WHT	GROUND
1	RED	ADVANCE TO PREPURGE INPUT

#### AQUATEK 50

7	ORG/BRN	PURGE PERMISSION OUTPUT
---	---------	-------------------------

#### AQUATEK 50

6	BLK	ADVANCE TO DESORB INPUT
5	BLU	GROUND

**INTERFACE CABLE FROM LSC-2 TO AQUATEK 50**  
**TEKMAR PART NUMBER 14-3783-000**

# Installation Instructions

---

## Interface Part Number 14-3783-000

Interfacing a Tekmar™ LSC-2 to an Aquatek 50

---

This interface allows automatic operation of the Aquatek 50 with the LSC-2. The two units "communicate" with each other in the following manner:

1. The LSC-2 sends a Purge Ready signal to the Aquatek 50, causing it to step to Prepurge.
2. The Aquatek 50 sends a signal to the LSC-2, causing it to step to Purge.
3. The LSC-2 sends a begin Desorb signal to the Aquatek 50, causing it to step to Desorb.

	<b>WARNING</b>	
<b>Electrical shock hazard inside. Unplug power cord before removing panels.</b>		

	<b>CAUTION</b>
<b>Make sure that the LSC-2 auto drain function is turned on before running the Aquatek 50. Failure to do so could damage the equipment! The auto drain switch is located on the back of the LSC-2 in the lower right-hand corner (rear view).</b>	

	<b>CAUTION</b>
<b>Circuit board components can be damaged by static discharge. Avoid touching the components unless otherwise noted.</b>	

### ***Connection to the Aquatek 50***

1. Plug the nine-pin connector into the receptacle labeled "concentrator I/O" on the rear of the Aquatek 50.

### ***Connection to the LSC-2***

1. Locate the eight-position terminal strip on the rear of the LSC-2.
2. Attach the orange or brown wire to position eight on the eight-position terminal strip.
3. Remove the back panel on the LSC-2.
4. Locate the two-position terminal strip. On newer units, it is labeled "T-2". It sits on a horizontal plane directly above the terminal strip where you just attached a wire. See Figure 9 in the back of your manual.
5. Connect the two wires to the two-position terminal strip. It does not matter to which screw each wire attaches.
6. Replace the back panel.
7. Locate the socket labeled "computer interface" on the rear of the LSC-2. It lies just below the eight-position terminal strip.
8. Insert the eight-pin connector into the socket, matching pin number one on the connector to pin number one on the socket. Pin orientation is correct if all three of the connector wires are on the left side of the socket (looking at the rear of the LSC-2).



7143 East Kemper Road, Cincinnati, Ohio 45242-9576  
 (800) 543-4461 • Outside the U.S. (513) 247-7000 • Service (800) 874-2004  
 Telefax (513) 247-7050