

Table 1. Typical Fittings for Columns and Inlet/Detector Liners, Adaptors, and Inserts.

Type	Description	Typical Use	HP Part No.
1/4- inch swage stainless steel, pkg. of 20 each	nut front ferrule back ferrule	1/4- inch packed metal columns	5080- 8753
1/8- inch swage stainless steel, pkg. of 20 each	nut front ferrule back ferrule	1/8- inch packed metal columns	5080- 8751
1/4- inch swage brass, pkg. of 20 each	nut front ferrule back ferrule	1/4- inch packed metal columns	5080- 8752
1/8- inch swage brass, pkg. of 20 each	nut front ferrule back ferrule	1/8- inch packed metal columns	5080- 8750
Vespel pkg. of 10	1/4- inch ferrule	inlet/detector liners, 1/4- inch glass packed columns	5080- 8774
Vespel pkg. of 10	1/8- inch ferrule	metal columns	0100- 1107
Graphite, pkg. of 10	1.0- mm ferrule	capillary columns	5080- 8773
Graphite, pkg. of 10	0.5- mm ferrule	capillary columns	5080- 8853
Graphite	6.35- mm O- ring	inlet/detector liners, 1/4- inch glass packed columns split capillary inlet insert	0905- 0767
Graphite	6.52- mm O- ring	splitless capillary inlet insert	0905- 1004
Silicone	6.0- mm O- ring	inlet/detector liners 1/4- inch glass packed columns split/splitless capillary inserts	0905- 0322
Silicone	1.0- mm O- ring	capillary columns	0905- 0759

NOTE: Dimensions given are ids of O- rings/Ferrules.