

## Vacuum System

**This table lists the replacement parts available for the vacuum system. If there is a part you need that is not listed in this chapter, contact your Hewlett-Packard Customer Engineer or your local Hewlett-Packard office.**

Description	HP Part Number
Blank flange, KF25	0100-1396
Calibration valve and vial	05971-60550
Clamp	
KF10/16, for foreline hose fitting	0100-1397
KF25 hinged, for GC/MSD interface	0100-1398
KF25 stamped, for calibration valve	0100-0549
KF40, for top of foreline trap	0100-1426
KF50 for diffusion pump	0100-1395
Diffusion pump, complete	
110/120 V ac	05971-60530
220/240 V ac	05971-60531
Fan	
110/120Vac	05971-60540
220/240 V ac	05971-60541
Foreline hose	05971-60119
Foreline pump, Edwards E2M2	
all 60 Hz and 100 V ac, 50 Hz MSDs	3162-0149
all other 50 Hz MSDs	05989-80024
Foreline trap	
body	05971-20562
body screen	05971-00560
cap	05971-20561
cap screen	05971-00561

continued

Description	HP Part Number
Foreline vacuum gauge	
fitting	05971-20530
gauge tube	0960-0540
O-rings and seals	
KF10 O-ring assembly	0100-1183
KF25 O-ring assembly	3162-0110
KF40 O-ring assembly	0100-1425
NW25 seal, vented	0100-1551
NW50 seal	05971-20177
top plate	05971-40100
Shipping clamp	
front	05971-00210
rear	05971-00211
Thermal switch, for diffusion pump	
too cold (bottom)	3103-0145
too hot (top)	3103-0146
Top plate (includes feed-throughs)	05971-60120
Vacuum manifold	05971-20550

---