

Analytical Instrument Management Methods: Cannabis Potency

The cannabis potency method supplied with the turnkey Agilent HPLC system from AIM provides fast, reliable, and accurate results; analysis after analysis.

Get industry-leading performance in cannabis potency testing with a refurbished Agilent HPLC from AIM at an amazing price. Analyze for THC, THCA, CBD, CBN, CBC, and more with a single 8-minute analysis (Fig. 1)

Enjoy the flexibility and peace of mind offered by the great linear response (0.5-500ppm) offered by proven hardware and the optimized AIM cannabis potency method (Fig. 2)

Don't see the analytes you're testing for? Contact us for custom method development and custom reporting options!




Fig 1. 10ppm Standard Sample









Figure 2. Linear response factor of selected analytes.



8392 S. Continental Divide Rd
Suite 102
Littleton, CO 80127
303.730.4136
Sales@aimanalytical.com